# Capo Natural

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product. As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

**Schütz Dental GmbH**  
 Dieselstr. 5-6 · 61191 Rosbach/Germany  
 Tel.: + 49 (0) 6003 814-0 · Fax: + 49 (0) 6003 814-906  
 info@schuetz-dental.de · www.schuetz-dental.com

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### Instructions for use

Capo Natural is a light curing composite containing an ultrafine, radiopaque glass filler and is indicated for placing fillings using advanced techniques. It can be polished to a high lustre. Due to the ultrafine particle filler, extremely homogeneous restorations can be placed which are easily polished to a high lustre. The chamelion effect matches the shade of the filling optimally to the tooth structure. The guidelines of EN ISO 4049 have been compiled with. Capo Natural is available in syringes and compules. The compules are for single use. Please do not reuse them, as this makes it impossible to rule out contamination and germ formation.

**Composition:**
Glass powder, diurethane dimethacrylate, silicon dioxide, Bis-GMA, tetramethylene dimethacrylate

**Total filler:**
75% by weight (53% by volume) anorganic filler (0.005-3.0 µm)

#### Indications

- Anterior and posterior restorations in Black's class I, II, III, IV, and V cavities.
- Inlays, onlays and laminate veneers
- Extended fissure sealing in molars and premolars
- Stump build-up
- Splinting loose teeth
- Adjusting the contours and shades to improve aesthetics

#### Application - Pretreatment

Before commencing the treatment, clean the tooth with non-fluoride polishing paste. Use a Vita® shade guide to select the shade while the tooth is still moist.

- Cavity preparation**
Minimally-invasive preparation of the cavity as generally required for adhesive techniques. All enamel margins in the anterior region must be bevelled. Do not bevel the margins in the posterior region and avoid slice preparations. Spray the cavity with water to clean it, remove all debris and dry it. The cavity must be isolated. It is advisable to place a rubber dam.
- Pulp protection / Cavity liner**
If an enamel-dentin adhesive is used, no cavity liner is required. In very deep cavities those areas in close proximity to the pulp must be coated with a calcium hydroxide material.
- Approximalk contact areas**
When filling cavities with approximal sections, place a transparent matrix and fix it in place.
- Adhesive system**
Etch (e.g. Capo Etch) and bond (e.g. Capo Bond) according to manufacturer's instructions.
- Application of the composite from the syringe**
Extrude the required amount of composite from the screw syringe, place it in the cavity with a standard metal instrument and contour it. The layers must not be thicker than 2 mm.
- Application of the composite from compules**
Place the compules in the dispenser. Remove the sealing cap. Position the compules in such a way that the opening is at a suitable angle for application within the cavity. Insert the material into the cavity while slowly and evenly applying pressure. Do not use excessive force! The layers must not be thicker than 2 mm. Once finished, pull back the punch in order to remove the compules from the dispenser. The compules can then be removed.

Note: for hygiene reasons, the compules are only intended for single use.
- Curing**

The curing time for all shades is 40 seconds per layer with a conventional halogen curing lamp or an LED curing lamp. With a plasma curing system, the curing time is 2 x 3 seconds. Hold the waveguide as close to the surface of the filling as possible. Fillings with more than one surface must be cured from the direction of each surface separately. Due to the influence of the oxygen in the ambient air, a thin, unpolymerized film remains on the surface, the so-called dispersion layer. This dispersion layer forms the chemical bond between the layers and must not be touched or contaminated with moisture.

#### 7. Trimming

Capo Natural can be trimmed and polished immediately after curing using finishing diamonds, flexible disks, silicone polishers and polishing brushes. Check the occlusion and articulation and spot grind to eliminate high spots or undesirable paths of articulation from the surface of the filling. Due to the influence of the oxygen in the ambient air, a thin, unpolymerized film remains on the surface, the so-called dispersion layer. This dispersion layer forms the chemical bond between the layers and must not be touched or contaminated with moisture.

#### Inlays, onlays and laminate veneers – Cavity preparation

The cavity should be prepared as minimally invasively as possible with only slightly diverging sides. To prevent the material fracturing, the layer must have a minimum thickness of 1.5 mm in the lateral and vertical aspects. All internal edges and angles must be bevelled. Avoid slice preparations. Prepare a flat cervical shoulder-do not bevel it. Any unavoidable undercuts must be blocked out with glass ionomer cement. Use slightly tapering diamonds with rounded tips for the preparation. Coat those areas of dentin in close proximity to the pulp with a thin layer of calcium hydroxide material. Cavity liners containing eugenol are contraindicated.

#### Impression and temporary restoration

Once the impression has been taken, a composite temporary restoration is fabricated (e.g. with Temdent, Schütz Dental). This may only be cemented with a non-eugenol cement.

#### Fabricating an inlay, onlay or laminate veneer

Cast the impression with hard stone plaster. Allow the model to set and pull off the impression. Block out the undercuts and apply an oil-free separating agent to the model. Build up the inlay on the model layer-by-layer. Build up the approximal and deep occlusal sections first. Each individual layer may be thicker than 2 mm and is cured separately with a commercially available light curing lamp (e.g. HiLite Power, Heraeus Kulzer, intermediate polymerization 90 seconds / final polymerization 180 seconds). The finished inlay is then released from the die. Trim and polish to a high lustre. Clean the inlay thoroughly with soap and water, rinse with air / water spray and dry.

#### Placing the inlay, onlay or laminate veneer

Remove the temporary restoration and clean the cavity. Place a rubber dam before cleaning and drying the prepared surfaces of the tooth. Exert gentle pressure on the inlay to check for fitting accuracy. Do not use force. If necessary, trim the fitting surfaces to improve the fit. The occlusion and articulation may not be checked when trying to fit the inlay as this could cause fractures.

Etching (e.g. with Capo Etch) and bonding (e.g. with Capo Bond) according to the manufacturer's instructions.

#### Fixing the restoration

The restoration is fixed with a commercially available dual-curing fixing composite (e.g. Alphalink Cem, Schütz Dental). Please adhere to the manufacturer's instructions.

#### Special notes

- The working time under a surgical lamp is approximately 2 minutes.
- In case of time-consuming restorations, the surgical lamp should be either temporarily moved away from the working area or the material should be covered by an opaque foil in order to prevent the composite from curing too early.
- Use a light-curing unit with an emission spectrum of 350- 500 nm for the polymerization of this material. As the required physical properties can only be achieved when the lamp works correctly, its luminous intensity must be checked regularly as described by the manufacturer.

Light intensity for curing	≥ 650 mW/cm²
Wavelength for curing	350 - 500 nm
Curing time	40 sec.

#### Hazard and Precautionary statements

#### Contains tetramethylene dimethacrylate

**Warning:** May cause an allergic skin reaction. Avoid breathing vapours / spray. Wear protective gloves / protective clothing / eye protection / face protection. If skin irritation or rash occurs: Get medical advice / attention.

#### Contraindications / interactions

If a patient has know allergies against or hypersensitivities to a component of this product, it may not be used or only under strict medical supervision by the doctor/dentist. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth before using the product.

#### Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us. To prevent possible reactions of the pulp in cavities where the dentine is exposed, the pulp must be protected adequately (e.g. calcium hydroxide preparation).

#### Interactions with other substances

As phenolic substances (such as eugenol) inhibit polymerization, do not use cavity liners (such as zincoxide eugenol cements) containing such substances.

#### Storage

Store at 10-25 °C (50 -77 °F). Avoid direct sunlight. Close the screw syringes tightly immediately after use. The material should be at room temperature before use. Retract the plunger of the syringe slightly to prevent the apertures becoming blocked. Do not use after expiry date (or label on syringe). For use by dentists only. Keep out of reach of children. This product was developed specifically for the described range of applications. It must be used as described in the instructions. The manufacturer is not liable for damage caused by handling or processing the material incorrectly.

\* Vita is a registered trademark of Vita® Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen.

Problem	Cause	Remedy
Composite does not cure properly	Light output of the light-curing lamp is inadequate <p>Emitted wavelength range of the lightcuring lamp is inadequate</p>	Check the light output and change the light source if required <p>Consult the manufacturer of the lightcuring lamp. <p>Recommended wavelength range: 350- 500 nm</p></p>
Composite in the syringe is sticky and soft, colourless liquid separates in the syringe	Material has been stored for a longer period at > 25 °C <p>Material has been kept in a syringe warmer for too much time</p>	Adhere to storage temperature. Store at 10- 25 °C (60- 77 °F). <p>Never keep a syringe in a syringe warmer for more than one hour per application</p>
Composite appears too hard and firm in the syringe	Material stored at temperatures < 10 °C (50 °F) for a longer period of time	Allow the composite to heat to room temperature before use; use a syringe warmer if necessary
Inlay / onlay is not properly retained when fitted	Syringe not properly sealed, composite partially cured	Always seal the syringe properly with the cap after taking out composite
Composite does not cure completely (dark or opaque shades)	Composite layers applied too thickly for each curing cycle	Adhere to a max. thickness of 2.0 mm per layer
Restoration appears too yellow compared with the shade guide	Inadequate curing of the composite layer	Repeat the exposure cycle several times; min. 40 seconds

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\*Vita es una marca registrada de Vita Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen