




Tizian Blank NEM

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Instruction for use

Intended use

Tizian Blank NEM is a beryllium- and nickel-free non-precious bonding alloy for the production of milled crowns and bridges.

General guidelines for handling

This instruction for use includes important processing steps and recommendations for Tizian Blank NEM.

Indication

- Crowns and wide-span bridges up to 16 units (max. 4 pontics in the anterior region and max. 3 pontics in the posterior region),
- Implant retained superstructures
- Implant bars.

Milling

Tizian Blank NEM is construed for CNC milling machines. Please follow the instructions and parameters of the respective manufacturer of CAM Software and the CNC milling machine.

CAD

The design should be done with appropriate CAD software. Please consider an anatomically reduced framework design for the veneering with ceramic. The wall thickness should not be less than 0.3 mm. Choose a sufficient connector dimension (6 - 9 mm²). Sharp edges and undercuts should be avoided.

Cutting out the frameworks from the Blank

Remove the milled frameworks with suitable cutting tools and smoothing the supports.

Soldering / Laser welding

Tizian Blank NEM can be soldered with all suitable solder. Tizian Blank NEM parts should not be soldered with gold or palladium solder. Tizian Blank NEM is also ideally suitable for laser-welding.

Preparation before ceramic veneering

The frameworks can be elaborated with standard carbide cutters, look for smooth transitions and avoid overlapping material. Please use the same cutter for one alloy to avoid contamination. The minimum thickness of the prepared coping should not be less than 0.3 mm. It's recommended to sandblast the frames with minimum 110 µm of Aluminium oxide with 3 - 4 bar and clean with steam cleaner. Oxide firing is not mandatory but can be done as an option for 5 minutes at 980 °C (cleaning firing). The frame needs to be sandblasted again to remove the present oxide layer thoroughly. In the end the cleaning by steam cleaner is mandatory. If you use a ceramic bonder please consider the instruction for use of the manufacturer.

Handling conditions / Safety

Metal dust is harmful to health. Use when grinding and sandblasting dust extraction and respirator with filter FFP3 – DIN EN 149.

Contraindications and side effects

If the instructions are observed during the production processes, incompatibilities with non-precious dental alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist regarding the contra-indications and side effects.

Single-use

Used milling disc should be not processed for further use as a medical device.

Disposal Instructions

Consult the material safety data sheets or national regulations for disposal.

Storage conditions

Temperature, humidity or light has no effect on the product properties.

Quantity

Please consider the label on the package.

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.

Description

NPM – Cobalt-based dental alloy for milling, type 4

Dimension

Ø 98,3 mm: 8 mm, 10 mm, 12 mm, 13,5 mm, 15 mm, 16 mm, 18 mm, 20 mm

Chemical Composition (typical values)

Co %	Cr %	W %	Si %	Mn %	Fe %	Sonstige %
61,65	27,75	8,45	1,61	0,25	0,2	< 0,1

Determination, extent and validity of the chemical composition according to DIN EN 10 204 – 3.1

Typical technical Data

Yield strength 0,2%	440 MPa
Elongation	14%
Tensile strength	819 MPa
E-module	240 GPa
Density	8,55 g/m ²
Corrosion resistance	< 200 µg/cm ²
Hardness*	275 HV 10/30
CTE (25 - 500 °C)	14,1 x 10 ⁻⁶ K ⁻¹
Max. firing temp.	~ 980 °C

*with appropriate sample preparation

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Schütz Dental GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany