dialog Bonding Fluid

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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
Instructions for use

Application as repair bonding for composite veneers

Roughen the veneer to be repaired up to 2 mm above the correction margin respectively the repair margin. Evenly apply the bonding fluid with a long-haired brush. Intermediately polymerize in a light polymerization unit to create a dispersion layer:

Spektra 2000 / 2000 fast	90 sec.
Spektra LED	30 sec.
Uni-XS / HiLite Power	90 sec.

If the intermediately polymerized surface appears white, it was polymerized excessively and must be renewed. Repeat the last step with reduced polymerization time.

Next, start layering the composite material and polymerize.

Please comply with the manufacturer's instructions.

Application as bonding agent between acrylic framework and composite veneers

The light-curing Bonding Fluid is suited to use on PMMA-based acrylics and on PEEK / PEKK.

Prepare the surfaces of the framework according to the manufacturer's instructions.

Evenly apply the bonding fluid with a long-haired brush. Intermediately polymerize in a light polymerization unit to create a dispersion layer:

Spektra 2000 / 2000 fast	90 sec.
Spektra LED	30 sec.
Uni-XS / HiLite Power	90 sec.

If the intermediately polymerized surface appears white, it was polymerized excessively and must be renewed. Repeat the last step with reduced polymerization time.

Next, start layering the composite material and polymerize.

dialog Bonding Fluid is not suited to the conditioning of frameworks made of metal, ceramics or zirconium dioxide.

Composition

Methyl methacrylate, urethane dimethacrylate, initiators, stabilizers.

Storage

Store at temperatures between 10 °C and 25 °C (50 °F and 77 °F) to ensure product stability.

Shelf-life

The maximum shelf-life is printed on the product label. Do not use after the shelf-life has expired.

Guarantee

Our technical instructions, regardless of whether they are provided verbally, in writing or during practical demonstrations, are based on our own experience and should only be considered as guidelines. As our products are subject to continued development, we reserve the right to modify them.

Side effects

Side effects of this medical device are to be expected only in extremely rare cases if the product is handled and used properly. However, immune reactions (e.g. allergies) or local discomfort cannot be excluded completely. If you become aware of any side effects, please notify us, even in case of doubt.

Contra-indications/interactions

If a patient has hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and cross reactions of the product with other materials already in the patient's mouth before using the product.

Note

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

Troubleshooting

Trouble	Cause	Remedy
Insufficient bond	Material applied unevenly, patchily	Ensure to apply an even, uniform film
Material evaporated, dried up	Bottle was left open	Close bottle tightly after each use

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