# Alphalink

---

CE 0297

**SCHÜTZ DENTAL**  
Micerium Group

Schütz Dental GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany  
Telefon +49 (0) 6003 814-0 · Telefax +49 (0) 6003 814-906  
info@schuetz-dental.de · www.schuetz-dental.de

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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## Directions for use

### General product description

Two-component composite for bonding prosthetic retentive elements or telescopic crowns or conical crowns and fabricated retentive elements to the framework casting. The material can be applied in the laboratory on the model or directly in the patient's mouth. With Alphalink, stress-free connections of metal parts can be created. A high proportion of inorganic fillers guarantees a high final hardness producing the stability required for laboratory preparation.

### Composition

Filler, Bis-GMA, tetramethylene dimethacrylate, dibenzoyl peroxide

### Indications

- Connecting prosthetic retaining elements
- Fixing telescopic or conical crowns to frameworks.
- Fixing pre-fabricated attachments to frameworks.

### Contraindications

- Do not use in case of a known allergy against one of the components.

### Application / processing of Alphalink

**Processing Alphalink to connect prosthetic retaining elements or telescopic crowns respectively conical crowns, and prefabricated retaining elements with the framework.**

### Preparatory steps

#### Construction of the bonding areas

Design retentive and stable bonding surfaces. When placing the secondary component onto the primary component, make sure that there is an even, thin gap between the two components. This gap ensures a passive fit and it is filled with Alphalink. When using prefabricated retaining elements, create undercuts in the joint area. Observe the manufacturer's instructions.

#### Preparing the bonding areas

After ensuring that both the primary and secondary components fit passively, create undercuts in them and sandblast with aluminium oxide (max. 125 µm).

After steam cleaning, condition the metal surfaces with Sebond Smart metal bonding system.

### Dosage

Alphalink is available in the user-friendly dual-compartment cartridge as well as in two individual syringes.

When using the **automix syringe**, please proceed as follows:

Place a disposable mixing tip onto the dual-compartment cartridge. Discard the first bit of mixed material discharged from the tip.

- Press Alphalink from the automix tip and apply the desired amount directly onto the restoration. Spread evenly.
- As the material will cure inside the used mixing tip, leave the tip on the cartridge as a closure until the next use.
- By using the automix cartridge with the mixing tip, Alphalink is always mixed in the ideal ratio of 1 : 1.

When using **individual syringes**, please proceed as follows:

Two equal lengths of components A and B are dispensed from the two syringes by pressing the plunger and are mixed homogeneously.

Working time after starting to mix:  
Approx. 3 min.

Setting commences after approx. 4.5 min.  
Setting completed after approx. 8 min.  
These times apply to an ambient temperature of 22 °C (71 °F).

When working inside the patient's mouth:

Working time after starting to mix: 45 sec.  
Setting commences after approx. 1 min.  
Setting completed after approx. 1.5 min.  
Higher temperatures accelerate the setting reaction and lower temperatures retard it.

#### Application of Alphalink at the laboratory

Alphalink applied to matrix and patrix. The metal components are repositioned on the model stress-free and fixed into place. Ensure that the areas to be bonded are completely filled with Alphalink and that there are not any air bubbles or contamination. Do not move the object from the time setting starts until the time when the polymerization is completed.

#### Application directly in the patient's mouth

Apply Alphalink with a spatula to the surfaces to be joined. The metal components are repositioned on the model stress-free and fixed into place. Do not move the object from the time setting starts until the time when the polymerization is completed.

## Troubleshooting

Working time too short	<ul style="list-style-type: none"> <li>- High ambient temperatures (e. g. summer heat or central heating) decrease the working time.</li> </ul>
Working time too long	<ul style="list-style-type: none"> <li>- If the material is used directly after taking it from the refrigerator, curing will be delayed.</li> </ul>
Material does not cure sufficiently	<ul style="list-style-type: none"> <li>- Substances containing eugenol or oil of wintergreen prevent curing.</li> <li>- Material components are insufficiently mixed; only use original mixing tip.</li> </ul>
Material is cured	<ul style="list-style-type: none"> <li>- After use, the cap was replaced instead of leaving the mixing tip on the cartridge (contamination of component A + B).</li> </ul>
Construction does not fit	<ul style="list-style-type: none"> <li>- Construction was moved before the material had cured.</li> <li>- Material is already too firm: Please join fewer parts at a time.</li> </ul>

### General Notes

- Alphalink is dimensionally stable up to 120 °C (248 °F).
- Alphalink is easy to grind.

### Storage

Refrigerate at 3 - 9 °C (37 - 48 °F).

### Shelf-life:

The maximum shelf-life is printed on the label of each packaging. Do not use after the expiration date.

### Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (e. g. allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

### Contra-indications / interactions

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

### Note

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets. Avoid contact of mucous membranes and skin with uncured material.

### Disposal

In accordance with local and/or statutory regulations.

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