Capo Natural

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Instructions for use

Capo Natural is a light curing, Capo Natural containing an ultrafine, radiopaque glass filler and is indicated for placing fillings using adhesive techniques. It can be polished to a high lustre. Due to the ultrafine particle filler, extremely homogeneous restorations can be placed which are easily polished to a high lustre. The chameleon effect matches the shade of the filling perfectly to the tooth structure. The guidelines of EN ISO 4049 have been compiled with. The Capo Natural is available in syringes and compules. The compules are for single use. Please do not reuse them, as this makes it impossible to rule out contamination and germ formation.

Composition:

Monomer matrix

Diurethane dimethacrylate, tetramethylene dimethacrylate, isopropylidene bis[2(3)-hydroxy-3(2)-(4-phenoxy)propyl] bismethacrylate

Total filler:

75% by weight (52% by volume) anorganic filler (0.04–3.0 µm)

Indications

- Direct anterior and posterior restorations in Black’s class I, II, III, IV, and V cavities.
- Indirect restorations such as inlays, onlays and laminate veneers
- Extended fissure sealing in molars and premolars
- Spinting opaque posts
- Endodontic mobile teeth
- Adjusting the contours and shades to improve aesthetics

Contraindications/interactions

If a patient has known allergies against or hypersensitivities to a component of this product, it may not be used or only under strict medical supervision by the doctor/dentist. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient’s mouth before using the product. Unpolymerized composite may cause skin allergies. The user must take adequate precautions. In case of irritation or allergy due to one of the constituents listed under “Composition”, do not use this material.

Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us. To prevent possible reactions of the pulp in cavities where the dentine is exposed, the pulp must be protected adequately (e.g. calcium hydroxide preparation).

Interactions with other substances

As phenolic substances (such as eugenol) inhibit polymerization, do not use cavity liners (such as zincoxide eugenol cements) containing such substances.

Application – Pretreatment

Before commencing the treatment, clean the tooth with non-fluoride polishing paste. Use a Vita® shade guide to select the shade while the tooth is still moist.

1. **Cavity preparation**

Minimal-invasive preparation of the cavity as generally required for adhesive techniques. All enamel margins in the anterior region must be beveled. Do not bevel the margins in the posterior region and avoid slice preparations. Spray the cavity with water to clean it, remove all debris and dry it. The cavity must be isolated. It is advisable to place a rubber dam.

2. **Pulp protection/Cavity liner**

If an enamel-dentine adhesive is used, no cavity liner is required. In very deep cavities those areas in close proximity to the pulp must be coated with a calcium hydroxide material.

3. **Approximal contact areas**

When filling cavities with approximal sections, place a transparent matrix and fix it in place.

4. **Etching gel**

Apply Etching first to the enamel areas of the cavity and let it take effect for 30 seconds. Then fill the whole cavity (dentine) with etching gel and let it take effect for a further 15 seconds. The etching time in the dentine should not exceed 20 seconds. Then rinse the cavity thoroughly with a water spray and dry it with oil-free compressed air. Do not overdry the dentine. Dried, etched enamel surfaces have a chalky-white appearance and must not be contaminated before the bondor is applied. If the surface becomes contaminated with saliva, rinse and dry again and re-etch if necessary.

5. **Applying the bonding agent**

Refer to the Capo Bond instructions for details.

6. **Placing the composite in the cavity**

Extrude the required amount of composite from the screw syringe, place it in the cavity with a standard metal instrument and contour it. The layers may not be thicker than 2 mm. Due to the effect of the oxygen in the air, a thin smear layer of unpolymerized material remains on the surface of each layer. This bonds the layers chemically and may not be touched or contaminated with moisture.

6a. **Using Capo Natural Composite Compules**

Place the compules in the dispenser. Remove the sealing cap. Position the compules in such a way that the opening is at a suitable angle for application within the cavity. Insert the material into the cavity while slowly and evenly applying pressure. Do not use excessive force! Once finished, pull back the punch in order to remove the compules from the dispenser. The compules can then be removed. Note: for hygiene reasons, the compules are only intended for single use.

7. **Curing**

All shades must be light cured for 40 seconds per layer using a commercially available light curing unit. The light guide must be held as close as possible to the surface of the filling. Multiple surface fillings should be light cured from each side.

8. **Trimming**

Capo Natural can be trimmed and polished immediately after curing using finishing diamonds, flexible disks, silicone polishers and polishing brushes. Check the occlusion and articulation and spot grind to eliminate high spots or undesirable paths of articulation from the surface of the filling.

Indirect method – Cavity preparation

The cavity should be prepared as minimally invasively as possible with only slightly diverging sides. To prevent the material fracturing, the layers must have a minimum thickness of 1.5 mm in the lateral and vertical aspects. All internal edges and angles must be rounded. Avoid slice preparations. Prepare a flat cervical shoulder – do not bevel it. Any un-avoidable undercut must be blocked out with glass ionomer cement. Use slightly tapering diamonds with rounded tips for the preparation. Coat those areas of dentine in close proximity to the pulp with a thin layer of calcium hydroxide material. Cavity liners containing eugenol are contraindicated.

Impression and temporary restoration

Once the impression has been taken (e.g. with Trisident, Schütz Dental), a composite temporary restoration is fabricated (e.g. with Temdent, Schütz Dental). This may only be cemented with a non-eugenol cement.

Fabricating an inlay

Cast the impression with hard stone plaster in the laboratory. Allow the model to set and pull off the impression. Block out the undercuts and apply an oil-free separating agent to the model. Build up the inlay on the model layer-by-layer. Build up the approximal and deep occlusal sections first. Each individual layer may not be thicker than 2 mm and is cured separately with a commercially available light curing lamp. The finished inlay is then released from the die and cured fully. Trim the occlusal surface with fissure burs and polish to a high lustre with silicone polishers and diamond paste. Clean the inlay thoroughly with soap and water, rinse with air/water spray and dry.

Placing the inlay, onlay or laminate veneer

Remove the temporary restoration and clean the cavity. Place a rubber dam before cleaning and drying the prepared surfaces of the tooth. Exert gentle pressure on the inlay to check for fitting accuracy. Do not use force. If necessary, trim the fitting surfaces to improve the fit. The occlusion and articulation may not be checked when trying to fit the inlay as this could cause fractures. Etch the tooth with Etching, spray it thoroughly with water to rinse it and blow the etched surfaces dry with oil-free compressed air. After etching and drying, the enamel surfaces appear chalky-white and may not be contaminated prior to applying Capo Bond. Should they become contaminated with saliva, rinse and dry them again – re-etch if necessary. Brush a thin coat of Capo Bond bonding agent to the etched enamel and cavity sides before curing it for 40 seconds with a light curing unit. A dispersion layer remains after curing which may not be removed as it bonds chemically to the filling material. Mix a commercially available dual-curing luting composite and use a spatula to apply it to the adhesive surfaces of the inlay. Carefully press the inlay into position. Use the spatula to remove the bulk of the excess material. Excess material on the approximal aspects is removed with a

probe and dental floss. To make certain that the inlay does not slip out of its proper position, press it gently with a ball-end plunger until it has been fully light-cured. Use a commercially available light curing unit to cure the composite for 40 seconds from each side, beginning with the approximal aspects. Remove the excess material with fine-grit diamonds and diamond finishing strips. Check the occlusion and adjust if necessary. Finish and polish with a finishing and polishing set. Caution: In case a layer of the inlay is thicker than 2 mm, a dual-curing composite should be used.

Special notes

- The working time under a surgical lamp is approximately 2 minutes.
- In case of time-consuming restorations, the surgical lamp should be either temporarily moved away from the working area or the material should be covered by an opaque foil in order to prevent the composite from curing too early.
- Use a light-curing unit with an emission spectrum of 350–500 nm for the polymerization this material. As the required physical properties can only be achieved when the lamp works correctly, its luminous intensity must be checked regularly as described by the manufacturer.

Hazard statement

May cause an allergic skin reaction.

Contains 1,4-Butanediol dimethacrylate

Caution: May cause an allergic skin reaction. Avoid inhaling steam/aerosols. Wear protective gloves/protective clothing/face protection. In the event of skin irritation or a rash, consult with a doctor or seek medical assistance.

Precautionary statements

Wear protective gloves/protective clothing/eye protection/face protection. If skin irritation or rash occurs: Get medical advice/attention.

Storage

Store at 10- 25 °C. Avoid direct sunlight. Close the screw syringes tightly immediately after use. The material should be at room temperature before use. Retract the plunger of the syringe slightly to prevent the apertures becoming blocked. Do not use after expiry date (refer to label on syringe). For use by dentists only. Keep out of reach of children. This product was developed specifically for the described range of applications. It must be used as described in the instructions. The manufacturer is not liable for damage caused by handling or processing the material incorrectly.

Troubleshooting

Problem	Cause	Remedy
Composite does not cure properly	Light output of the light-curing lamp is inadequate <p>Emitted wavelength range of the lightcuring lamp is inadequate</p>	Check the light output and change the light source if required <p>HConsult the manufacturer of the lightcuring lamp. Recommended wavelength range: 350- 500 nm</p>
Composite in the syringe is sticky and soft, colourless liquid separates in the syringe	Material has been stored for a longer period at > 25 °C <p>Material has been kept in a syringe warmer for too much time</p>	Adhere to storage temperature, store in a refrigerator <p>Never keep a syringe longer in a syringe warmer for more than one hour per application</p>
Composite appears too hard and firm in the syringe	The material has not been heated to room temperature after taking it out of the refrigerator	Allow the composite to heat to room temperature before use; use a syringe warmer if necessary
Inlay/onlay is not properly retained when fitted	Syringe not properly sealed, composite partially cured <p>Restoration is too opaque to be cemented using only light-curing composite</p>	Always seal the syringe properly with the cap after taking out composite <p>Use dual-curing luting composite</p>
Composite does not cure completely (dark or opaque shades)	Composite layers applied too thickly for each curing cycle	Adhere to a max. thickness of 2.0 mm per layer
Restoration appears too yellow compared with the shade guide	Inadequate curing of the composite layer	Repeat the exposure cycle several times; min. 40 seconds

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