Tizian Blanks PMMA

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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Instructions for using Tizian Blanks PMMA

Indications:

Tizian Blanks PMMA consist of millable PMMA (acrylic). They are available in different tooth shades as well as in uncolored, transparent.

The tooth-colored Tizian Blanks PMMA are suited to the production of crowns and bridges as long-term temporaries.

The transparent Tizian Blanks PMMA will burn out residue-free. They are suited to the production of lost-material casting molds as well as to the production of functional splints and drilling templates.

All Tizian Blanks PMMA can be used with or without veneers. In addition, it is possible to check fit and function of the restorations on a plaster model.

The constructions are modelled with Tizian CAD/CAM as usual and are then milled from Tizian Blanks PMMA.

Framework:

- Max. 2 pontics and max. 1 end pontic is permitted for the posterior teeth.
- The minimum wall thickness (without veneer) is: 1.0 mm (occlusal), 0.6 mm (cervical)
- Connector cross sections
incisors $>10 \text{ mm}^2$
posterior teeth $> 15 \text{ mm}^2$

Milling

- Always mill the material in an ambient temperature of $20 \text{ }^\circ\text{C} - 25 \text{ }^\circ\text{C}$ ($68 \text{ }^\circ\text{F} - 77 \text{ }^\circ\text{F}$).
- Please adhere to the instructions of your milling system.
- Tizian Blanks PMMA are suited to dry processing only.

- All **tools** used must be suited to processing PMMA. We recommend sharp-edged radius mills.
- The milling parameters strongly depend on the milling system, on the CAM software and on the **milling geometries**. If you are uncertain, please use your parameters for PMMA.

Separating the finished work

The finished work should be carefully separated from the blanks using a mill to stab out the material (stitching mill).

Finishing / polishing:

- Finish with clean, cross-toothed hard-metal mills.
- In order to avoid thermal damage ensure low speed, a minimal contact pressure force and a suitable amount of cooling.
- To polish, use conventionally available, high-gloss composite polishing pastes (e.g. dialog Vario Polish by Schütz Dental), goat hair brushes and cotton buffs.

Caution:

- Do not overheat the acrylic when using mills, steam cleaners etc.
- When grinding, the thickness of the framework must not fall below 0.4 mm.

Cleaning the frameworks

- Alcohol in high concentrations generally causes tiny cracks (so-called crazing). Therefore, do not use alcohol-based cleaning agents or disinfectants.
- When steam-cleaning, do not overheat the surfaces.
- Alternative cleaning methods are mechanical cleaning methods.

Veneer

The veneer is made with veneer composites, e.g. dialog Vario (Schütz Dental). A bonding agent must be used.

Fixation**when used as a temporary:**

The bridge can be fixed with all temporary cements.

Wearing time for temporaries made of this material:

The wearing time is limited to 12 months maximum.

Storage:

Store in a dry place.

Disposal:

The material can be disposed of with household waste.

Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-indications / interactions:

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

Note:

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

Guarantee:

Our processing instructions, regardless of whether they are issued verbally, in writing or in the form of a practical demonstration, are based upon our own experience and can therefore be regarded solely as guidelines.

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Schütz Dental GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany