
IMPLA

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Implantat-System

Artikel-Nr.: 608025

SCHÜTZ DENTAL
Micerium Group

CE 0297

Schütz Dental GmbH • Dieselstr. 5-6 • 61191 Rosbach / Germany • Tel.: +49 (0) 6003 814-0 • Fax: +49 (0) 6003 814-906
www.schuetz-dental.de • info@schuetz-dental.de

Mandler 09 / 2018 – 200

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


List of contents

Page

IMPLA Implant Systems Dual Surface and Micro Retention	3–6
IMPLA – Interim implants	7–8
IMPLA – Implant Systems Cylindrical	9–12
IMPLA – Mini-Implants conetop/balltop	13–14
General notes on implantation	15–16
General notes on IMPLA Abutments	17
Ball Anchors	18–21
Locator® Implant Abutment System	22–26
HSL-Abutments	27–28
Multi-Unit Abutments	29–30
IMPLA PreFace Abutments	31
IMPLA surgery kit incl. surgery tray	
IMPLA surgery kit incl. surgery tray and IMPLA 3D Navi surgery tray	32–33
Disinfection, cleaning and sterilization of surgical and prosthetic tools and instruments of the IMPLA system	34–37
User and safety instructions	38
Manufacturer’s information on rotating surgical instruments	
Standardized manual reprocessing of rotating instruments	39–40
Information on storage, disinfection, cleaning procedures and sterilization of rotating instruments	41
Annex: Forms and information sheets for patients	42–53
„Medical examination for a dental implantation“	43–44
„Informed Consent“	45–47
„Information Sheet“	48
„Surgery Log File“	49–50
„Implantological aftercare in terms of a recall“	51
„Instructions“	52–53

IMPLA Implant Systems Dual Surface and Micro Retention

IMPLA implants are made of pure titanium (ISO 5832-2). The IMPLA system comprises surgical, prosthetic and laboratory components and instruments.

Important notes, warning

IMPLA implants are part of a complete system and may only be used with the original components and instruments of the system. The use of components which are not part of this system voids the guarantee and rules out parts being replaced. The IMPLA system may only be used by dentists, doctors and dental technicians familiar with dental implantology including diagnosis, preoperative planning, surgical techniques and prosthodontics.

Consultation is done in written form or by electronic media as well as with demonstrations during special workshops. All information is based on the state of the art and latest scientific findings known at the time of launching the product. However, we urge the users of our products to check the suitability of the product in each case. This is especially important for techniques not expressly recommended by the manufacturer. In case of doubt, we encourage our customers to confer with Schütz Dental Group/Section IMPLA.

Reuse and sterilization

All abutments and components are for single use only and must not be reprocessed. Any reuse involves the risk of cross-contamination! All instruments, abutments and components are packaged non-sterile and must be sterilized (EN 13060) prior to placement (see page 26–35). Using unsterile components increases the infection risk for the patient.

Indications

IMPLA implants Dual Surface and Micro Retention are indicated for:

- restoring individual tooth gaps
- bridgework spanning interdental and free-end gaps (only 1 cantilever unit)
- restoration of edentulous mandibles and maxillae

Further indications:

Implants with a diameter of 3.3 mm

These implants are an alternative in case of limited jaw ridge width. They are, because of their lower mechanical strength when compared with the standard diameter (4.2 mm), intended for the following indications:

- a) restoring individual anterior teeth of the mandible
- b) restoring individual posterior teeth of the maxilla
- c) four interforaminal implants to restore the edentulous mandible with a bar construction
- d) **Do not** use with ball anchors, magnets or telescopes!

Implants with a diameter of 4.2 mm

The IMPLA standard implants with a diameter of 4.2 mm are used for all indications of oral enossal implants for mandibles and maxillae. A prerequisite is an adequate minimal amount of bone base with a width of about 5.5 mm.

Implants with a diameter of 5.3 mm

These implants are indicated for the use in the posterior tooth areas of mandible and maxilla (especially for the restoration of molars). The jaw ridge width should not be less than 6.5 mm.

The implants are primarily indicated for delayed implantations.

Contraindications

When selecting patients, the general contraindications for dental surgery must be taken into account. These include: Vascular diseases, uncontrolled diabetes mellitus and other metabolic or systemic diseases affecting bone regeneration, blood coagulation disorder, anticoagulative therapy, bone metabolism diseases.

Temporary contraindications: Acute and chronic inflammations, inadequate soft tissue present, chemotherapy or radiotherapy.

Oral contraindications:

- patient not really prepared to undergo total oral rehabilitation
- inadequate oral hygiene
- untreated parafunctions, e. g. bruxism
- inadequate amount of bone available
- inadequate interocclusal distance
- inadequate occlusion and/or articulation

Important notes

When working with our implants and instruments intraorally, always use a rubber dam to protect against aspiration!

Primary stability directly after the insertion of the implant is absolutely necessary for a successful healing. During the healing phase, IMPLA implants are generally never exposed to a masticatory load.

Surgical procedure

To create optimum conditions for the implant to heal, it is essential to treat the hard and soft tissues gently. It is advisable to use the countersink and bone thread cutter in order to ensure that the implant can be placed atraumatically applying as little torque as possible (maximum 30 Ncm). As thermal trauma may impede or prevent the implant healing, heat generation should be kept to a minimum as follows:

- recommended drill speeds:
2.0 mm Ø: approx. 800 rpm, 2.5 mm Ø: approx. 700 rpm, 3.15 mm Ø: approx. 500 rpm, Ø 4.2 mm, Ø 4.8 mm, Ø 5.3 mm: approx. 400 rpm
- recommended speeds: pilot drill Ø 1.2 mm and Ø 2.0 mm - approx. 800 rpm, enlargement drill Ø 3.15 mm - approx. 500 rpm, cavity drills of all Ø - approx. 400 rpm
- intermittent drilling and cutting using sharp drills and cutters (use 10 times only)
- cool the instrument externally with sterile physiological saline (NaCl) or Ringer's solution.

Important

To determine the optimum implant length and prevent nerv injuries, check in advance that adequate bone is available (X-ray).

Placing the implant after preparing the implant bone supporting area

We advise to study our detailed surgery documentation CD thoroughly.

Removing the insertion aid

In order to remove the insertion aid (blue), pull off the plastic implant carrier and loosen the vertical screw (blue) with the hexagonal screwdriver AF 1.2 mm and dispose the vertical screw. Secure the implant against rotation with help of a counter wrench in cases of low primary stability and concomitant augmentation in the according alveole. Afterwards, the insertion aid can be lifted from the implant easily.

Closing the soft tissue

The cover/healing aid must be screwed onto the implant before the wound is closed with atraumatic suturing material.

Recovering the implant

Once the individual healing period (minimum of three months for the mandible and six months for the maxilla) has elapsed, the implant should be recovered and the cover/healing screw removed. In case of augmentation techniques, the healing period must be prolonged accordingly. The interior of the implant is then cleaned and the correct height of gingiva former for the oral mucosa is screwed in. Make sure that the gingiva former fits properly! The restoration should not be fitted to the implant until the soft tissue has healed fully and shows no irritations.

Impression taking

Open impression technique:

Once the gingiva former has been removed, the impression post should be screwed into the implant. When taking the impression, the impression material is first syringed around the impression post. Once the impression material has cured, loosen the screw, remove the impression from the mouth together with the impression post and close the implant again with the gingiva former. Carefully screw an implant analog into the impression post and cast the model – preferably including a gingival mask. Once the cast has set, remove the impression post, respectively transfer coping, and replace it with the required abutment.

Closed impression technique:

The impression is taken with a standard impression tray and with help of plastic caps which audibly click onto the implant. The plastic caps remain in the impression and are fitted with lab implants to create a plaster model.

Laboratory components

When selecting the abutments, the angulation of the implant and the thickness of the gingiva must be taken into account. Those areas of the abutment which are in contact with the implant may not be blasted with abrasives or treated in any other manner.

Documentation

Every IMPLA implant can be traced via its article and lot number. The label on the outer blister package includes all the data for recording in the patient's file. Any claims can only be asserted when you provide us with the requested information, indicating the above mentioned data.

Packaging and sterility

IMPLA implants are packed sterilely in a special clear acrylic container and are color-coded. The chosen implant is removed from the sterile container with the insertion aid only directly before insertion and is inserted into the prepared bone cavity without touching the implant. A special attachment on the bottom of the implant container holds the sealing, covering and healing screw which can be easily removed with the screwdriver 1.2 mm wide hexagon wrench key. The outer cardboard box shows the user important data such as implant length, lot and CE certification.

General information

- All implants are for single use only.
- Each is packed sterilely in a closed container, please note the expiration date of sterility.
- The implant is sterile (beta sterilization) as long as the packing remains unopened and undamaged. A subsequent sterilization by the user is not permitted by the manufacturer. A chemical sterilization of the implant (chemiclave) may not be carried out under any circumstances, as this could damage the implant surface. In case of damaged packaging, the implant is no longer sterile and may not be used anymore.
- Do not touch the implant.
- Store the sealed and undamaged implant container in a dry place at room temperature. Avoid direct light incidence.
- System errors due to mixing up instruments and implants should be avoided. Please always check the color codings and labelings.
- Thermal damages of the jaw bone by rotating instruments (such as pilot drills, thread cutters or enlargement drills) must be avoided.
- On the following pages, please find information on the preoperative diagnostics, indications and contraindications, risk factors, post-operative treatment of the implant as well as patient information sheets, the „Informed consent“ form and surgery record form „surgery log file“.

Torque of screw connections

- 10 Ncm horizontal screw connections
- 20 Ncm vertical screw connections
- 30 Ncm insertion of the implant, Locator Abutment screw connections

IMPLA – Interim implants

Note:

Before performing implantology treatment, it is essential to be familiar with pertinent professional literature. It is absolutely necessary to have experience in surgery and in the assessment of the normal and pathological course of treatment.

The „General notes on implantation“ (see page 15–16) must be considered before treatment.

Description:

The one-part IMPLA interim implant is made of grade 5 titanium in accordance with ISO 5832-3.

Warning: If the implant breaks, it is possible that fragments remaining cannot be removed without damage to the surrounding bone tissue!

Implant threads: The IMPLA interim implant has primary self-cutting threads, allowing primary stable insertion into the jaw bone while preserving bone material. The screw design and resulting compression stability together with the pilot drill ensure sufficient primary stability even in spongy bone tissue. Implant surface: for the intended use as interim implant, the implant surface has been machine polished. This is to make removal of the implant after the desired period possible.

Important: Due to a potentially high insertion torque, the IMPLA positioning key (IPK) is to be used with caution to avoid the danger of fracturing or deforming the implant.

Indications:

The IMPLA interim implant is a temporary support for drill templates or for a temporary bridge/denture during the healing phase of the permanent implants for patients who desire a stable, esthetically pleasing temporary prosthesis. The soft tissue around conventional implants is protected. Any existing augmented area is protected from strain. The mucosa must be healthy with good circulation.

Contraindications:

Absolute:

- A sufficient number of temporary immediate implants cannot be placed in the implant area to support the temporary denture.
- Patients whose systemic condition does not allow implantation.

Relative:

- Patients who are not capable of follow-up care for psychological or physical reasons and whose compliance is not anticipated
- Patients undergoing treatment for a malignant tumor
- Patients with uncontrolled diabetes
- Bruxism
- Heavy smokers
- Insufficiently available bone tissue and/or immediate vicinity to endangered structures (nerve, sinuses)

Diagnostics / patient information:

Diagnostics consist of the patient history, clinical examination, X-ray examination with OPG and CT if necessary. A medical check-up is recommended. Implantations are elective procedures for which the patient must be given comprehensive information. Treatment alternatives and consequences of the implantation must be indicated. In particular, the importance of oral hygiene and costs of further treatment must be discussed. The relevant laws indicate the kind and extent of information required.

Side effects:

- Even in a professional surgical procedure, potential horizontal and vertical bone loss must be taken into account as with all dental implants; the extent and kind of bone loss cannot be estimated in advance.
- Reversible or irreversible impairment or interaction is possible when implants are inserted in the vicinity of certain structures (e.g. adjacent teeth, maxillary sinus, or nerves).
- In rare cases, chronic pain can occur in association with dental implants.

To minimize the potential risk of infection, antibiotic prophylaxis can be recommended. In rare cases, allergies to the components of the implant may prevent successful implantation.

Procedure:

The entire course of treatment has three stages:

1. Pilot drilling;
2. Initial insertion of the implant;
3. Further insertion using the IMPLA positioning key (IPK)

For stage 1) A pilot drilling is made through the mucosa at the optimum site, guided and controlled by the axis and length of the implant. A minimum distance of about 2 mm from the final implants must be maintained. The pilot drill must be selected depending on the bone quality (D1–D2 => 2.0 mm bit in the lower jaw; D3–D4 and upper jaw => 2.0 mm bit with less depth). Cooling with sterile saline solution is required. To prevent potential thermal damage, for bone quality D1–D2, the holes should be pre-drilled with the pre-drill and the cavity widened with the 2.0 mm pilot drill if necessary.

Important: Drill alternately at a speed of 500–600 rpm.

The entire implant length must be pre-drilled. In the upper jaw, less if the bone quality is poor (30–50%). For stage 2) The packaging may not be opened until immediately before insertion of the implant. Hold the implant at the thread using the packaging as protection and open the side facing the implant head. Place the positioning key (IPK) firmly on the implant head and make sure it is gripping firmly. Do not pull the implant out of the packaging until you are near the implantation site and insert it there immediately using the IMPLA positioning key (IPK) to prevent contamination.

For stage 3) Lower the implant with the IMPLA positioning key (IPK) until the upper edge of the thread is flush with the gingiva. If necessary, the torque wrench must be used to reach the final position.

Important: Due to great torque, the IMPLA positioning key (IPK) is to be used with utmost caution to avoid fracturing the implant.

Removal of the temporary implant:

The temporary implant must be removed prior to or at the time the permanent implants are inserted. Detach the temporary implant and remove all remaining cement. Put the IMPLA positioning key (IPK) and the ratchet on the implant and inspect for correct fit. Paying attention to torque, turn the implant counterclockwise to unscrew. Remove any granulation tissue with a curette, as conditions similar to tooth extraction can be expected. **The temporary implant is usually worn for 3 to 6 months. Implant healing must be monitored (e. g. X-ray) regularly. Infected or loose implants must be removed promptly.**

Notes:

- After implantation, the type of implant used and LOT number must be recorded in the documentation (sticker).
- Implants must be stored in a closed container in a dry place. The seal may not be broken until immediately before insertion. Before insertion, any contact with foreign materials must be avoided. The endosseous part of the implant may not be touched. In the event the implant, superstructure, or an accessory is swallowed, the location of the object must be determined and the appropriate medical steps are to be taken.
- Healing at the implant site must be monitored regularly, e. g. by X-ray. Implants that are infected or not osseointegrated must be removed promptly. The dentist must determine the time for removal.

Scope of delivery:

The IMPLA interim implant is delivered in non-sterile condition.

The user must sterilize the implants in an autoclave at 134°C before use. We accept no responsibility for the correct handling and functioning of the sterilizer.

The implant is sealed in a sterilizable film bag bearing a label for entering the implant data such as length and LOT number that can be put in the patient file. This film wrapping is the actual packaging. The IMPLA interim implant is intended to be used one time only!

If the packaging is damaged or not airtight, there is a risk of contamination with bacteria. The implant may not be used.

IMPLA-Implant Systems Cylindrical

IMPLA implants are made of pure titanium (ISO 5832-2). The IMPLA system comprises surgical, prosthetic and laboratory components and instruments.

Important notes, warning

IMPLA implants are part of a complete system and may only be used with the original components and instruments of the system. The use of components which are not part of this system voids the guarantee and rules out parts being replaced. The IMPLA system may only be used by dentists, doctors and dental technicians familiar with dental implantology including diagnosis, preoperative planning, surgical techniques and prosthodontics.

Consultation is done in written form or by electronic media as well as with demonstrations during special workshops. All information is based on the state of the art and latest scientific findings known at the time of launching the product. However, we urge the users of our products to check the suitability of the product in each case. This is especially important for techniques not expressly recommended by the manufacturer. In case of doubt, we encourage our customers to confer with Schütz Dental GmbH/Section IMPLA.

Reuse and sterilization

All abutments and components are for single use only and must not be reprocessed. Any reuse involves the risk of cross-contamination! All instruments, abutments and components are packaged non-sterile and must be sterilized (EN 13060) prior to placement (see page 26–35). Using unsterile components increases the infection risk for the patient.

Indications

IMPLA implants Cylindrical Implants are indicated for:

- restoring individual tooth gaps
- bridgework spanning interdental and free-end gaps, in accordance with conventional and special principles of prosthetic treatment
- restoration of edentulous mandibles and maxillae

Further indications:

Implants with a diameter of 3.6 mm

These implants are an alternative in case of limited jaw ridge width. They are, because of their lower mechanical strength when compared with the standard diameter (4.5 mm), intended for the following indications:

- a) restoring individual anterior teeth of the mandible
- b) restoring individual lateral anterior teeth of the maxilla
- c) four interforaminal implants to restore the edentulous mandible with a bar construction

Implants with a diameter of 4.5 mm

The IMPLA standard implants with a diameter of 4.5 mm are used for all indications of oral enossal implants for mandibles and maxillae. A prerequisite is an adequate minimal amount of bone base with a width of about 5.8 mm.

Implants with a diameter of 5.5 mm

These implants are indicated for the use in the posterior tooth areas of mandible and maxilla (especially for the restoration of molars). The jaw ridge width should not be less than 6.8 mm.

The implants are primarily indicated for delayed implantations.

Contraindications

When selecting patients, the general contraindications for dental surgery must be taken into account. These include: Vascular diseases, uncontrolled diabetes mellitus and other metabolic or systemic diseases affecting bone regeneration, blood coagulation disorder, anticoagulative therapy, bone metabolism diseases.

Temporary contraindications: Acute and chronic inflammations, inadequate soft tissue present, chemotherapy or radiotherapy.

Oral contraindications:

- patient not really prepared to undergo total oral rehabilitation
- inadequate oral hygiene
- untreated parafunctions, e. g. bruxism
- inadequate amount of bone available
- inadequate interocclusal distance
- inadequate occlusion and/or articulation

Important notes

When working with our implants and instruments intraorally, always take increased care to protect against aspiration!

Primary stability directly after the insertion of the implant is absolutely necessary for a successful healing. During the healing phase, IMPLA implants are generally never exposed to a masticatory load.

Surgical procedure

To create optimum conditions for the implant to heal, it is essential to treat the hard and soft tissues gently. It is advisable to use the countersink and bone thread cutter in order to ensure that the implant can be placed atraumatically applying as little torque as possible (maximum 30 Ncm). The so-called „Dense Bone Drills“ have been designed especially for the cortical bone preparation. They are meant to reduce the insertion torque for longer implants. As thermal trauma may impede or prevent the implant healing, heat generation should be kept to a minimum as follows:

- recommended drill speeds:
2.0 mm Ø: approx. 800 rpm, 2.75 mm Ø: approx. 700 rpm, 3.0 mm Ø and 3.4 mm Ø: approx. 500 rpm, 3.7 mm Ø, 3.9 mm Ø, 4.25 mm Ø, 4.9 mm Ø and 5.25 mm Ø: approx. 400 rpm
- intermittent drilling and cutting using sharp drills and cutters (use 10 times only)
- cool the instrument externally with sterile physiological saline (NaCl) or Ringer's solution.

Important

To determine the optimum implant length and prevent nerv injuries, check in advance that adequate bone is available (X-ray).

Placing the implant after preparing the implant bone area

We advise to study our detailed surgery documentation CD/System overview thoroughly.

Removing the insertion post

In order to remove the insertion post (red, blue, brown), pull off the plastic implant carrier and loosen the vertical screw (blue) with the hexagonal screwdriver 1.2 mm. Dispose of the vertical screw. In cases with concomitant augmentation or where the primary stability is low, secure the implant against rotation with help of a counter wrench. Afterwards, the insertion aid can be lifted from the implant easily.

Closing the soft tissue

The cover/healing aid must be screwed onto the implant before the wound is closed with atraumatic suturing material.

Recovering the implant

Once the individual healing period (minimum of three months for the mandible and six months for the maxilla) has elapsed, the implant should be recovered and the cover/healing screw removed. In case of augmentation techniques, the healing period must be prolonged accordingly. The interior of the implant is then cleaned and the correct height of gingiva former for the oral mucosa is screwed in. Make sure that the gingiva former fits properly! The restoration should not be fitted to the implant until the soft tissue has healed fully and shows no irritations.

Impression taking

Once the gingiva former has been removed, the impression post should be screwed into the implant. When taking the impression, the impression material is first syringed around the impression post. Once the impression material has cured, loosen the screw, remove the impression from the mouth together with the impression post and close the implant again with the gingiva former. Carefully screw an implant analog into the impression post and cast the model – preferably including a gingival mask. Once the cast has set, remove the impression post, respectively transfer coping, and replace it with the required abutment.

Laboratory components

When selecting the abutments, the angulation of the implant and the thickness of the gingiva must be taken into account. Those areas of the abutment which are in contact with the implant may not be blasted with abrasives or treated in any other manner.

Documentation

Every IMPLA implant can be traced via its lot number. The label on the outer blister package includes all the data for recording in the patient's file. Any claims can only be asserted when you provide us with the requested information, indicating the above mentioned data.

Packaging and sterility

IMPLA implants are packed sterile in a special glass container.

The chosen implant is removed from the sterile container with the insertion aid only directly before insertion and is inserted into the prepared bone cavity without touching the implant.

A special attachment on the bottom of the implant container holds the sealing, covering and healing screw which can be easily removed with the screwdriver 1.2 mm wide hexagon wrench key.

The outer cardboard box shows the user important data such as implant length, implant diameter, article no., lot and CE certification.

General information

- All implants are for single use only.
- Each is packed sterilely in a closed container, please note the expiration date of sterility.
- The implant is sterile (beta sterilization) as long as the packing remains unopened and undamaged. A subsequent sterilization by the user is not permitted by the manufacturer. A chemical sterilization of the implant (chemiclave) may not be carried out under any circumstances, as this could damage the implant surface. In case of damaged packaging, the implant is no longer sterile and may not be used anymore.
- Do not touch the implant.
- Store the sealed and undamaged implant container in a dry place at room temperature. Avoid direct light incidence.
- System errors due to mixing up instruments and implants should be avoided. Please always check the color codings and labelings.
- Thermal damages of the jaw bone by rotating instruments (such as pilot drills and extension drills) must be avoided.
- Please find information on the preoperative diagnostics, indications and contraindications, risk factors, post-operative treatment of the implant as well as patient information sheets, the „Informed consent“ form and surgery record form „surgery log file“ in our IMPLA implants instruction manual.

Torque of screw connections

10 Ncm horizontal screw connections

20 Ncm vertical screw connections

30 Ncm insertion of the implant, Locator Abutment screw connections

IMPLA-Mini-Implants conetop / balltop

Note:

Before performing implantology treatment, it is essential to be familiar with pertinent professional literature. It is absolutely necessary to have experience in surgery and in the assessment of the normal and pathological course of treatment.

The „General notes on implantation“ (see page 15–16) must be considered before treatment.

Description:

The one-piece IMPLA Mini-Implant is made of unalloyed grade 4 titanium in accordance with ISO 5832-2.

Warning: If the implant breaks, it is possible that fragments remaining cannot be removed without damage to the surrounding bone tissue!

Implant threads: The IMPLA Mini-Implant has primary self-cutting threads, allowing primary stable insertion into the jaw bone while preserving bone material. The screw design and resulting compression stability together with the pilot drill 2.0 mm and the extension drill 2.5 mm ensure sufficient primary stability even in spongy bone tissue.

Implant surface: The implant surface was designed in accordance with the implant's application as permanent implant (see indications).

Indications:

The IMPLA Mini Implant was developed for the following indications:

Mini implant cone „conetop“:

- Fixing full dentures
- Bar-supported dentures (on a minimum of 4 implants)
- Support implant for fixing bridges

Mini implant ball-head „balltop“:

- Fixing full dentures

Contraindications:

Absolute:

- Patients whose systemic condition does not allow implantation.

Relative:

- Patients who are not capable of follow-up care for psychological or physical reasons and whose compliance is not anticipated
- Patients undergoing treatment for a malignant tumor
- Patients with uncontrolled diabetes
- Bruxism
- Heavy smokers
- Insufficiently available bone tissue and/or immediate vicinity to endangered structures (nerve, sinuses)

Diagnostics / patient information:

Diagnostics consist of the patient history, clinical examination, X-ray examination with OPG and CT if necessary. A medical check-up is recommended. Implantations are elective procedures for which the patient must be given comprehensive information. Treatment alternatives and consequences of the implantation must be indicated. In particular, the importance of oral hygiene and costs of further treatment must be discussed. The relevant laws indicate the kind and extent of information required.

Side effects:

- Even in a professional surgical procedure, potential horizontal and vertical bone loss must be taken into account as with all dental implants; the extent and kind of bone loss cannot be estimated in advance.
- Reversible or irreversible impairment or interaction is possible when implants are inserted in the vicinity of certain structures (e. g. adjacent teeth, maxillary sinus, or nerves).
- In rare cases, chronic pain can occur in association with dental implants.

To minimize the potential risk of infection, antibiotic prophylaxis can be recommended. In rare cases, allergies to the prosthetic components (immediate implantation) may prevent successful implantation.

Procedure:

The entire course of treatment has four stages:

1. Pilot drilling, 2.0 mm
2. Extension drilling, 2.5 mm
3. Initial implant insertion
4. Further insertion all the way to the implant's final position with help of the IMPLA screwdriver.

Point 1) A pilot drilling is made through the mucosa at the optimum site, guided and controlled by the axis and length of the implant. Subsequently, the cavity is further prepared with the extension drill 2.5 mm. A minimum distance of about 2 mm from implants already in place must be maintained.

The treatment must be adapted to the bone situation:

(D1–D2 => 2.0 mm => 2.5 mm drill in the lower jaw

**D3–D4 and upper jaw => 2.0 mm => 2.5 mm
drill with decreased drilling depth)**

Cooling with sterile saline solution is required.

Important: Drill alternately at a speed of 500–600 rpm.

The entire implant length must be predrilled.

Point 2) The preparation depth with the extension drill 2.5 mm should be adapted to the bone quality (see point 1).

Point 3) The packaging may not be opened until immediately before insertion of the implant. Remove the glass tube with the implant inside and remove the implant with the insertion aid. With help of the insertion aid, insert the implant into the cavity. Once an increasing torque is reached, the plastic insertion aid is removed.

Point 4) Use the IMPLA screwdriver to screw in the implant until the top edge of the structured implant area is level with the bone. If necessary, use the torque wrench to reach the final position.

Important: Because of the thread configuration of the implant, a torque of min. 25 to max. 50 Ncm should be reached, depending on the indication.

Notes:

- After implantation, the type of implant used and LOT number must be recorded in the documentation (sticker).
- Implants must be stored in a closed container in a dry place. The seal may not be broken until immediately before insertion. Before insertion, any contact with foreign materials must be avoided. The endosseous part of the implant may not be touched. In the event the implant, superstructure, or an accessory is swallowed, the location of the object must be determined and the appropriate medical steps are to be taken.
- Healing at the implant site must be monitored regularly, e. g. by x-ray. Implants that are infected or not osseointegrated must be removed promptly. The dentist must determine the time for removal.

Scope of delivery:

The IMPLA Mini-Implant is delivered sterile.

The implant is sealed in a film bag bearing a label for entering the implant data such as length and LOT number that can be put in the patient file. This film wrapping is the actual packaging.

The IMPLA Mini-Implant is intended to be used one time only!

If the packaging is damaged or not airtight, there is a risk of contamination with bacteria. The implant may not be used in this case.

General notes on implantation

Attention

All rotating instruments are subject to wear. Therefore, we cannot warrant the instruments' (pilot drills, extension drills, countersinks, cortical drills and thread cutters) functionality for more than 10 times of use, depending on the bone density. If the instruments are used more often, the manufacturer cannot guarantee a proper cutting efficiency and minimal temperature development any longer.

- Insufficient primary stability can be caused by poor preparation of the implant site or by insufficient bone quality. In this case, the loose implant should be removed and be replaced by an implant with a larger diameter. Implants should always be screwed into the bone all the way, even with their thread-free neck, when a covered healing process is considered. Above all, a primary stability of the implant must be achieved.
- Implants should only be implanted when all prerequisites for a complication-free treatment and a longterm success are given. If there is any concern regarding the success of the implantation and/or the stability of the desired implant-based denture, another solution without implantation should be favored.
- All rotating instruments can, in principle, be used several times. It is the user's/implantologist's personal responsibility to decide for each instrument whether and how often the instrument should be re-used or if it is to be discarded, considering the particular case at hand and the possible wear of the instrument.

Care of instruments and parts

- Used instruments/parts must always be disinfected before cleaning.
- Damaged instruments/parts must be disinfected, cleaned, and disposed of.

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Please consider our information on each material on page 28–31.

Disinfection

- Instruments/parts must be submerged in an adequate disinfection solution.
- Disinfection may be carried out manually or automatically. Please observe the manufacturer's instructions concerning dosing/concentration, residence time, temperature.
- In order to avoid damages, the instruments/parts may not have any contact during disinfection.
- It is highly recommended to wear protective gloves when working with contaminated instruments/parts.

Cleaning

- Only use adequate cleaning agents.
- Always observe the manufacturer's instructions concerning concentration of the agent, temperature, etc.
- It is highly recommended to wear protective gloves (DIN EN 455) when cleaning contaminated instruments/ parts manually.
- In case of automatic cleaning and disinfection, the instruments/parts may not have any contact in order to avoid damages.
- Only instruments/parts made of the same material should be cleaned together.

- In case of heavy contamination of the instruments/parts, ultrasound cleaning techniques are recommended.
- After cleaning, rinse the instruments/parts thoroughly with water and dry them immediately.
- Damaged instruments/parts must be removed and disposed of separately.

Sterilization

- All surgical instruments and parts must be sterilized after cleaning.

Please consider our information on each material on page 28–31.

Storage

- The sterilized instruments/parts should be stored in a dry and dust-free atmosphere at room temperature.
- After the expiration of the maximum storage period, the instruments/parts must be sterilized again.

Dose 25 kGy/min

General notes on IMPLA Abutments

Important notes / warnings:

IMPLA Abutments are a part of the IMPLA system and must only be used with the proper original IMPLA components. For the insertion of the abutments, only use the IMPLA instruments and tools. Otherwise, a damage of the original components cannot be ruled out.

Delivery:

All IMPLA Abutments are supplied in unsterile packaging. Therefore, the components must be sterilized accordingly before use.

Material:

IMPLA Abutments are manufactured from titanium grade 5 (Ti6Al4V ASTM 136).

Side-effects:

Allergies and sensitivities towards the alloy cannot be ruled out completely and may occur in very rare cases.

Interaction:

In case of occlusal contact, the presence of different alloy types in the patient's mouth may lead to galvanic reactions.

Safety advice:

Metal dust and ceramic dust are hazardous to the health. When working on these materials, always wear a protective mask and use an extraction system with a fine particle air filter commonly used in dental practices.


Indications:

Please carefully read the instructions and the list of indications of the selected abutment type.


Contra-indications:

When selecting the patients, it is extremely important to ensure that they are able to maintain the necessary oral hygiene. In addition, the general rules and limitations in surgery and prosthetics for the selection and treatment of patients apply.


Furthermore, please adhere to the instructions and according description of indications as well as to the separately listed contra-indications of each abutment type.


BALL ANCHORS

light, ca. 600 g

standard, ca. 800 g

heavy, ca. 1200 g

The new snap button system features several improvements towards the present anchor systems. The retention mechanism is based upon the newly conceived design of frico-retention. As this mechanism is obtained with help of a thermoplastic acrylic, the matrix will always keep its perfectly round shape. This system does completely without activation or deactivation with spring or multi-disk. If the inner matrix is to be replaced by a completely new one during a recall, this can be done within only a few seconds with help of a specially constructed torque socket wrench. Three different, color-coded types of inner matrix with different withdrawal forces are available. In addition, a high-precious metal inner matrix is now available for extreme loads.

Experience has shown that extreme rotative forces may occur in a two anchor system without laminar support. This fact indicates metal matrices rather than acrylic matrices. We have met these circumstances with the new high-precious metal matrix.

In constructions with laminar support, however, inner matrices with acrylic rings can be used without a problem.

If an inner matrix made of acrylic in an existing denture seems too weak, it can be replaced with a high-precious metal matrix without any problem and without replacing the retention cap! This is the great advantage of this system.


Withdrawal forces (mean reference value):

- yellow = 600 g
- green = 800 g
- red = 1200 g
- gold = 800 g

The easy and quick replacement of the inner matrix saves you valuable time. The dentist obtains the safety of reproducible withdrawal forces without any manual manipulation.


INSTRUCTIONS

Initial situation with two anchors

The following work steps for the production of a ball anchor (atrix) supported denture are illustrated with help of a case with two anchors (on implants or root caps) on teeth 21 and 27 (US dental numbering system).

Use of the mounting matrix

The inner mounting matrix (2) is screwed into the required retention cap (3) with help of the torque socket wrench (1) and set onto the patrix. (Inner mounting matrices with blue acrylic rings must not be used as permanent inner matrices for dentures!)

Placing the matrix

In a paralleling device, both matrices are placed into the same denture insertion position with help of a parallelism equalizing tool. Any undercuts are blocked out with wax.


Blocking out

Both mounting matrices are blocked out with wax.

Mounting the spacers

For stability reasons, we recommend a cast framework support. If this technique is employed, spacers are set onto the matrices with help of the respective paralleling device.

Covering the alveolar ridge

The complete alveolar ridge basis is covered with a wax layer of 0.5 mm thickness.

Please note: You may skip this workstep, if the denture does not have a framework support.

Modellation, framework support

Framework support modelled in wax.

Finished framework support

Cast framework support as seen from occlusal direction.

Attention: The retention caps must not be connected to the framework support (laser fused or welded), but may only be cemented or fixed with attachment glue.

Setting up teeth, investing and polymerizing


Finished tooth set-up as seen from labial direction.

Afterwards, the denture is invested and polymerized.

Fixing the inner matrix

After finishing and polishing the denture, the inner mounting matrix can be removed and replaced by the required, permanent inner matrix.


Finished denture

Finished denture as seen from labial direction.


Recall - replacing an inner matrix

We recommend a replacement of the inner matrix during each recall. The inner matrix can be screwed out by three revolutions with the torque socket wrench and is then replaced with a new one.

BALL ANCHOR WITH REPLACEABLE INNER MATRIX


Indication partial and hybrid dentures implants with a ball coupling of \varnothing 2.25 mm

Contraindication unilateral cantilever prostheses without transversal blocking


Patrices two different patrices available gingiva height 2 mm and 4 mm


Inner matrices four inner matrices with different withdrawal forces available
 inner matrix yellow: light, ca. 600 g
 inner matrix green: standard, ca. 800 g
 inner matrix red: strong, ca. 1200 g
 inner matrix gold: ca. 800 g

Retention cap retention cap

Torque socket wrench socket wrench for matrices with fixed torque of 4 Ncm = ca. 400 g

Materials titanium grade 5

titanium alloy

composition: Ti 90.0 %,
 Al 5.5 – 6.5 %,
 V 3.5 – 4.5 %

melting interval 1604°C – 1660°C
 hardness, cured 350 – 385 HV5

LOCATOR® IMPLANT ABUTMENT SYSTEM

IMPORTANT: This document contains the most current instructions for use. Please read and retain.

CLASSIFICATION: Universal hinge, resilient attachment for endosseous implants.

ATTACHMENT COMPONENT ORDER NUMBERS

Product numbers of Locator Implant Abutments vary with implant type, implant diameter, and tissue cuff height. It is necessary to have this information available to place your order.

PARTS IDENTIFICATION

Implant Abutment
 3.3/3.0 mm **No. 636060**
 3.3/5.0 mm **No. 636061**
 4.2/3.0 mm **No. 636062**
 4.2/5.0 mm **No. 636063**
 5.3/3.0 mm **No. 636064**
 5.3/5.0 mm **No. 636065**
 Implant Abutment Titanium alloy with TiN coating

Locator matrix set
 3 pcs.
No. 636070
 Titanium Cap with black Low Density Polyethylene Male, Block Out Spacer, and Nylon Repl. Males (clear, pink, blue)

Locator plug connector clear, 4 pcs.
No. 636071
 Nylon (clear)
 pull-off force: 2.270 g

Locator plug connector, pink slight retention, 4 pcs.
No. 636072
 Nylon (pink)
 pull-off force: 1.360 g

Locator plug connector, blue extra-light retention, 4 pcs.
No. 636076
 Nylon (blue)
 pull-off force: 680 g

Locator plug connector, green extended application range, 4 pcs.
No. 636073
 Nylon (green)
 pull-off force: 1.820 g


Locator plug connector, red extended application range, slight retention, 4 pcs.
No. 636074
 Nylon (red)
 pull-off force: 450 g

Locator impression post
No. 636067
 Aluminum housing with black Low Density Polyethylene Male

Locator lab implant (4 mm)
No. 636068
 Aluminium

Locator-Instrument, threepart (Male Removal Tool, Male Seating Tool and Abutment Driver)
No. 636066

Locator Adapter with WS connection
No. 636075

INDICATIONS

The LOCATOR® Implant Attachment System is designed for use with overdentures or partial dentures retained in whole or in part, by endosseous implants in the mandible or maxilla.

CONTRAINDICATIONS

Not appropriate where a totally rigid connection is required. Use on a single implant with divergence of greater than 20 degrees is not recommended.

CAUTION

Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed dentist.

SINGLE-USE DEVICES

Locator Males: The inadvertent re-use of Locator nylon males could cause loss of retention for the overdenture due to wear from previous use or damage during removal with the Locator Core Tool.

Locator Abutments: The inadvertent re-use of Locator abutments could contain patient contamination build-up and subsequent wear of the retention bands. This would result in the device to perform with improper fit and function which would result in loss of retention for the prosthesis.

STERILIZATION

All components and instruments are supplied **NON-STERILE**.

Titanium abutments may be sterilized by autoclave or dry heat sterilization using the following parameters:

1. Autoclave sterilize using 121°C (250°F), (15 – 20 psi at sea level), for 20 minutes minimum.
2. Dry heat sterilize using 170°C (338°F) for 2 hours minimum.

Locator Core-Tools (**the disassembled state only**) may be sterilized by Autoclave or Dry Heat sterilization using the following parameters:

1. Autoclave sterilize using 121°C (250°F), 15 – 20 psi (at sea level), for 40 minutes minimum.
2. Dry heat sterilize using 170°C (338°F) for 2 hours minimum.

FEATURES

1. **LOWEST VERTICAL HEIGHT:** The total height of the Locator Attachment (abutment plus male) is only 3.17 mm on an externally hexed implant.
2. **LOCATING DESIGN:** Self-locating design allows a patient to easily seat their overdenture without the need for accurate alignment of the attachment components.
3. **RETENTION INSIDE AND OUT:** The patented Dual Retention innovation provides the Locator Attachment with greater retention surface area than ever before available with other attachments. A combination of inside and outside retention ensures the longest lasting performance.
4. **ROTATIONAL PIVOTING ACTION:** The design of the pivoting Locator Male allows a resilient connection for the prosthesis without any resulting loss of retention. The retentive nylon male remains completely in contact with the abutment socket while its titanium denture cap has a full range of rotational movement over the male.
5. **USE WITH NON-PARALLEL IMPLANTS:** The Locator Replacement Males 636072, 636076 can be used to restore an implant with up to 10 degrees of divergence (20 degrees between implants). The Locator Extended Range Replacement Males 636073, 636074 can accommodate a divergent implant between 10 and 20 degrees (40 degrees between implants).

plug connector, red
No. 636074

plug connector, blue
No. 636076

plug connector, pink
No. 636072

plug connector, green
No. 636073

PLACEMENT OF THE LOCATOR IMPLANT ABUTMENT

1. To select the proper Locator Implant Abutment, determine the type of implant and the diameter being used. Then measure the tissue thickness from the apical rim of the implant body to the crest of the gingiva at the highest side of the implant site. Choose the corresponding abutment tissue cuff height that exactly equals the tissue measurement, or is the next closest higher size available. The exact tissue cuff height of Locator abutment will position the proper 1.5 mm of working attachment above the surrounding gingival level (**which should not be submerged below the tissue**).
2. After the secondary gingival healing period is complete, remove the healing cuff according to instructions provided by the manufacturer of the implant system being used.
3. It is imperative that all bone and soft tissue be removed from the superior aspect of the implant body to guarantee complete seating of the Locator Implant Abutment.
4. A special gold plated Abutment Driver (contained in Locator Core Tool, 636066) is designed to engage the inside diameter of the Locator Abutment and thread it into the implant. (Fig. 1)
5. Final torque tightening of the Locator Abutment to prevent screw loosening is achieved using the 30 N cm Torque Wrench Kit.

NOTE: Use the IMPLA Torque wrench with either of these options to achieve the minimum seating force of 30 N cm that will help prevent screw loosening of the locator Implant Abutment.

LOCATOR MALE PLACEMENT BY THE DENTIST

1. Insertion of the proper Locator Implant Abutment at tissue level must be completed (see Section A-1) before beginning the procedure for placement of the Locator Male.
2. Place a White Block-Out Spacer (636070 Package) over the head of each Locator Abutment. (Fig. 2) The spacer is used to block out the area immediately surrounding the abutment. The space created will allow the full resilient function of the pivoting metal denture cap over the Locator Male.

NOTE: If the White Block-Out Spacer does not completely fill the space between the tissue and the metal denture cap, it is necessary to block out any remaining undercuts to prevent the added acrylic resin from locking the denture onto the abutment. This can be accomplished by stacking more Block-Out Spacers.

3. Insert a Locator Cap with Black Processing Male (636070 Package) into each Locator Implant Abutment, leaving the White Block-Out spacer beneath it (Fig. 3). The Black Processing Male will maintain the overdenture in the upper limit of its vertical resiliency during the processing procedure.
4. Prepare a recess in the denture to accommodate the protruding Locator Male. There must be no contact between the denture and the titanium cap. If the denture rests on the metal cap, excess pressure on the implant will result.
5. Use a light curing device to light cure bond the Locator Denture Cap Male into the denture (Fig. 4 and 5), or mix a permanent self-curing acrylic and place a small amount in the recess of the denture and around the metal cap of the Processing Cap Male.
6. Insert the denture into position in the oral cavity. Guide the patient into occlusion, maintaining a proper relationship with the opposing arch. **Maintain the denture in a passive condition, without compression of the soft tissue, while the acrylic sets. Excessive occlusal pressure during the setting time may cause tissue recoil against the denture base and could contribute to dislodging and wear of the nylon males.**
7. After the acrylic resin has cured, remove the denture and discard the white spacer. Use a bur to remove excess acrylic, and polish the denture base before changing to the final male.


Fig. 1

Abutment Driver/
Locator Core Tool
No. 636066


Fig. 2

White Block-out
Spacer
No. 636070


Fig. 3

Processing Cap
Nr. 636070


Fig. 4



Fig. 5


Acrylic Resin

8. Use the Locator Male Removal Tool (attached to the Locator Core Tool, Zest order 636066) to remove the Black Processing Male from the metal denture cap. The sharp circular edge on the end of the removal tool should be wedged tightly down into the very bottom of the plastic male so that it will catch the inside of the black plastic insert and pull it at an angle out of the metal housing (Fig. 6). To discard the nylon male from the tip on the Core Tool, point the tool down and away from you and tighten the Male Removal Tool clockwise back onto the Core Tool. This will activate the removal pin and dislodge the nylon male from the tip end of the Male Removal Tool.

Male Removal Tool/
Locator Core Tool
No. 636066

9. The Locator Male Seating Tool (contained in Locator Core Tool, 636066) is used to firmly push a Locator Replacement Male into the metal Denture Cap (Fig. 7). The replacement male must be placed securely in level with the rim of the cap (Fig. 8).

Male Seating Tool/
Locator Core Tool
No. 636066


NOTE: The Replacement Male will not stay on the tool when it is turned upside down due to the varying sizes of males available. It is best to hold the denture with the base side down and snap the male into the metal denture cap.

10. Instruct the patient in the path of insertion. Have the patient insert and remove the appliance several times.

LOCATOR MALE PLACEMENT BY THE LABORATORY


In the Operatory:

1. Insertion of the proper Locator Implant Abutment at tissue level must be completed (see Section A-1) before beginning the following impression procedure.
2. Place a Locator Impression Coping with Black Processing Male (636067) onto each Locator Abutment (Fig. 9).
3. Take an impression using a firm body impression material, exercising caution not to compress the soft tissue. The Locator Impression Coping is designed with minimum retention to be picked up with the impression material.
4. Snap a Locator Female Analog (636068 – 4 mm) onto each Impression Coping in the impression. The analog female must not fall off when turned upside-down with vibration (Fig. 10).


Impression Coping
No. 8505

NOTE: An alternative reline impression technique using the patient's prosthesis is possible with use of the Locator Black Processing Cap Male (636070 Set). When the impression is withdrawn, the Processing Cap Male will remain on the abutment. Remove the Processing Cap Male from each abutment and snap it onto a Locator Female Analog. Reposition this assembly back into the impression making sure it is fully seated.


Impression Coping
No. 8505

5. Pour the master cast. Upon separation, the Locator Female Analog is a part of the master cast replicating the position of the Locator Implant Abutment in the oral cavity.
6. Before waxing and processing the appliance, place a Locator Cap with Black Processing Male into each Female Analog in the master cast (Fig. 11). Make sure the male is fully seated.
7. Set the teeth and wax the appliance. Proceed with the processing technique of your choice through the boil-out step.
8. After the boil-out, remove the Processing Cap Male. Place a White Block-Out Spacer over the head of each Female Analog. The spacer is used to block out the immediate area surrounding the Locator Implant Abutment. The space created will allow the full resilient function of the pivoting metal denture cap over the Locator Male.

Female Analog
Impression
No.

9. Re-insert the Locator Black Processing Cap Male into each Female Analog, leaving the White Block-Out Spacer beneath it. The Black Processing Male will maintain the overdenture in the upper limit of its vertical resiliency during the processing procedure.

Denture Base
Processing Model
No.

10. Complete the processing and discard the white spacer. Avoid damage to the final male by polishing the denture base before changing to the final male.
11. Use the Locator Male 636066 Removal Tool attached to the Locator Core Tool to remove the Black Processing Male from the metal denture cap. The sharp circular edge on the end of the removal tool should be wedged tightly down into the very bottom of the plastic insert so that it will catch the inside of the black plastic insert and pull it at an angle out of the metal housing (Fig. 12).
12. The Locator Male Seating Tool is used to firmly push a Locator Replacement Male into the empty metal denture cap (Fig. 13). The replacement male must seat securely into place, level with the rim of the cap (Fig. 14).

NOTE: The Replacement Male will not stay on the tool when it is turned upside down due to the varying sizes of males available. It is best to hold the denture with the base side down and snap the male into the metal denture cap.

HOW TO CHANGE THE LOCATOR MALE

1. The Locator Core Tool (636066), which contains a Locator Male Removal Tool and Locator Male Seating Tool, is used to remove the nylon male from the metal denture cap and replace it with another Locator Replacement Male.
2. Use the Male Removal Tool attached to the Locator Core Tool to remove the nylon male from the metal denture cap. The sharp circular edge on the end of the removal tool should be wedged tightly down into the very bottom of the plastic male so that it will catch the inside of the plastic insert and pull it at an angle out of the metal housing.
3. The Male Seating Tool is used to firmly push a Locator Replacement Male into the empty metal denture cap. The replacement male must sit securely in place, level with the rim of the cap. Use of multiple Locator attachments (3 or more) in the same dental arch may require use of the 3.0 lbs. (light retention) pink colored Replacement Male No. 636072 or 1.5 lbs. (extra light retention) blue colored Replacement Male No. 636076, for easier removal of the prosthesis by the patient.

NOTE: The Replacement Male will not stay on the tool when it is turned upside down due to the varying sizes of males available. It is best to hold the denture with the base side down and snap the male into the metal denture cap.

PATIENT CARE

Good oral hygiene is vital to implant success. The Locator Implant Abutments must be thoroughly cleaned each day to prevent wear of the abutments due to buildup of abrasive plaque in the socket of the abutment. The use of a soft nylon bristle or end-tufted toothbrush, and superfloss to polish the abutments should be taught. A non-abrasive gel toothpaste, and an irrigation system is recommended to keep the socket of the Locator Abutment clean.

Patients should maintain a three to four month recall for cleaning and implant evaluation. The inside socket of the Locator Abutment and the sulcus area around the implant abutment are the primary areas of concern. Use plastic instruments for scaling the abutments. Do not use metal instruments which may create scratches on the abutment surface. Examine patients for signs of inflammation around the implant abutments, and for implant mobility. Use a 30Ncm torque wrench to make sure the Locator Implant Abutment is tight before dismissal.

We continually strive to improve our products, and therefore, reserve the right to improve, modify or discontinue products and components at any time without notice or incurring obligation. Purchaser assumes all risks and liability resulting from the use of our products, whether used separately or in combination with other products not of our manufacture.


Fig. 12

Male Removal Tool/
Locator Core Tool
No. 636066


Fig. 13

Male Seating Tool/
Locator Core Tool
No. 636066

Fig. 14

HSL-Abutments

Description:

The IMPLA HSL-Abutment is used for the production of individualized, occlusally screw-fixed abutments in the cast-on technique. The abutment is composed of a base and a plastic tube (screw canal). The base is manufactured of a high-fusing cast-on alloy and the plastic tube burns out completely residue-free. During the cast-on procedure, the base is covered by the melt of a high-precious casting or cast-on alloy. The metallic link is formed by diffusion of the components of both alloys.

Important notes / warnings:

IMPLA HSL-Abutments are a part of the IMPLA system and must only be used with the proper original IMPLA components (IMPLA implants). For the insertion of the abutments, use only the IMPLA instruments and tools. Otherwise, a damage of the original components cannot be ruled out.

Indications:

Single tooth gaps, partially or fully edentulous jaws, cement-fixed crowns or bridges, screw-fixed crowns.

Contra-indications:

Do not fix inappropriate abutments on IMPLA implants. Ensure an accurate selection of the abutment version (diameter of the connector and implant geometry). Manipulating the connector geometry will lead to an imprecise fit which will render the abutment useless.

Delivery:

All HSL-Abutments are supplied in unsterile packaging. Therefore, the components must be sterilized accordingly before use.

Procedure:

Modelling	The plastic tube is modelled as customary in dental technology. To do so, it can be shortened individually, ensuring the preservation of the connector geometry towards the HSL-base and keeping the connector geometry of the abutment absolutely clean to avoid cast-on in this area.
Investing	Ensure that the investment remains bubble-free. Only use phosphate-bonded investment materials (e. g. MANI Schütz products) when casting high-precious metals.
Preheating	In order to allow for a complete cast-on, please adhere to the instructions of the preheating furnace as well as to the preheating temperatures and holding times prescribed by the manufacturer of the alloy.
Casting	Please adhere to the instructions of the casting machine. To prevent the IMPLA HSL-Abutment from melting, please ensure that the casting temperature of the ceramic alloy you use is significantly lower than 1400°C (solidus temperature).
Cooling	We recommend letting the casting cylinder cool to room temperature to avoid tension in the cast object.
Divesting	To avoid a deformation of the object, especially of the connector geometry towards the implant, only blast with beads at low pressure. Afterwards, the object can be stripped in an appropriate stripping solution.
Finishing	We recommend using ceramic grinders and hard-metal mills. Adhere to industrial safety regulations (protective goggles, extraction system, respiratory protection)!
Veneering	Please note the coefficient of thermal expansion of each veneering material to prevent cracking.

Composition of the IMPLA HSL-Base:

Gold	60 %
Platinum	24.9 %
Palladium	15 %
Iridium	< 1%

Technical data:

Density	18.1 g/cm ³
Melting interval	1350 – 1460°C
Vickers hardness HV5/30	150 – 205 – 230

Side effects:

Allergies and sensitivities towards the alloy cannot be ruled out completely and may occur in very rare cases.


Interaction:

In case of occlusal contact, the presence of different alloy types in the patient's mouth may lead to galvanic reactions.

Multi-Unit Abutments

Description:

The IMPLA Multi-Unit Abutment Concept facilitates improved usage of the available bone structure. By means of an angulated insertion of the terminal implants, the available bone is used in distal direction which improves the stability. With the according angulation, the Multi-Unit Abutments ensure a uniform insertion direction of the prosthetics to facilitate a tension-free screw fixation. For this purpose, the IMPLA Multi-Unit Abutment Concept contains straight and angulated components.

Components Multi-Unit System Cone Connection:

Abutments with the diameters 3.3 mm, 4.2 mm, 5.3 mm, each with gingiva heights of 1 mm and 3 mm, and angulations of 0°, 20° and 30° (not 5.3 mm).

Accessories IMPLA Multi-Unit System:

Impression system with impression aid and fixing screw as well as the according lab implant.

Gingiva sleeve for conditioning of the mucosa as well as the appropriate metal or plastic sleeve for the prosthetics with suitable screw for secondary screw fixation.

Important notes / warnings:

IMPLA Multi-Unit Abutments are a part of the IMPLA system and must only be used with the proper original IMPLA components (IMPLA implants). For the insertion of the abutments, use only the IMPLA instruments and tools. Otherwise, a damage of the original components cannot be ruled out.

Indications:

The IMPLA Multi-Unit Abutments are exclusively meant for a use with occlusal screw-fixed prosthetics (e. g. bridges, hybrid prostheses and cast bars).

A minimal implant length of 11.5 mm in the distal position is indicated for the use of IMPLA Multi-Unit Abutments.

The minimal number of implants for a complete reconstruction is four in the mandible and six in the maxilla.

A sufficient bone quality is a basic requirement for implantation, as it ensures the primary stability. The proper torque is between 30 and 40 Ncm.

The distal implants must be inserted with an angulation that corresponds to the angulation of the abutments inserted later, i.e. **0°, 20° or 30°**.

Presurgical planning as well as a comprehensive medical history and diagnosis (including radiological examination) are absolutely necessary.

In principle, the general rules and instructions for surgical and prosthetical procedures in dental implantology apply.

Contra-indications:

When selecting the patients, it is extremely important to ensure that they are able to maintain the necessary oral hygiene. In addition, the general rules and limitations in surgery and prosthetics for the selection and treatment of patients apply.

Delivery:

All components of the IMPLA Multi-Unit Abutment Concept are supplied in unsterile packaging. Therefore, the abutments must be sterilized accordingly before use.

Procedure:

Based upon presurgical planning and following the drill protocol for IMPLA implants, a pilot drilling for a mesial implant without angulation is performed first.

Subsequently, all other implants can be aligned to this first drilling with help of a drilling template. The drilling template improves the spacial orientation and presets the angulation of the distal implants. The drill is guided through a sleeve in the drilling template into the bone in the exact angle intended.

After the cavities have been prepared in accordance with the IMPLA drilling protocol, the implants can be inserted.

The prosthetic restoration can be carried out either in two stages or by immediate load.

In case of a two-stage restoration, we recommend to cover the implant with a surgical cover screw (covered healing) or with a gingiva former for open healing.

A primary stability of 30–40 Ncm is absolutely necessary for the immediate load approach.

Impressions are taken with help of the according impression components of the IMPLA Multi-Unit Abutment Concept.

Next, straight abutments are fixed to the mesial implants and angulated abutments are fixed to the distal implants in accordance with the previously selected implant angulation (0° , 20° , 30°). The abutments are fixed with the appropriate vertical screw with a torque of 20 Ncm.

Finally, the prosthetic construction is set upon the abutments tension-free and is fixed with a prosthetic screw with a torque of 20 Ncm.

Material:

The IMPLA Multi-Unit Abutments are made of titanium grade 5 (Ti6Al4V ASTM 136).

Side effects:

Allergies and sensitivities towards the alloy cannot be ruled out completely and may occur in very rare cases.

Interaction:

In case of occlusal contact, the presence of different alloy types in the patient's mouth may lead to galvanic reactions.

IMPLA PreFace Abutments

Description:

The IMPLA PreFace Abutment is finished with an appropriate milling system with help of the original PreFace Abutment holder. The PreFace Abutment holder serves to fix the PreFace Abutment inside the milling system in the appropriate working position. The PreFace Abutments are retained in the PreFace Abutment holder on the face. Ensure that they are exactly in end position. The PreFace Abutment holders are available for many different milling machines.

Important notes / warnings:

IMPLA PreFace Abutments are a part of the IMPLA system and must only be used with the proper original IMPLA components (IMPLA implants). For the work with the PreFace Abutments, only use the appropriate PreFace Abutment holder.

Indications:

The IMPLA PreFace Abutment serves for the production of individual titanium abutments on IMPLA implants in CAD/CAM milling machines. The individual abutments can be combined with crowns or bridges to reconstruct function and aesthetics.

Contraindications:

The IMPLA PreFace Abutments must only be combined with IMPLA implants. Manipulating the connector geometry will lead to an imprecise fit which will render the abutment useless. The recommended torque for the fixation of IMPLA PreFace Abutments is 20 Ncm.

Delivery:

The IMPLA PreFace Abutments are supplied in unsterile packaging. Therefore, the components must be sterilized accordingly before use.

Procedure:

The shape of the PreFace Abutment can be adjusted individually to the present anatomic conditions. The production of the PreFace Abutment is done digitally with an appropriate CAD software. The basic rules must be adhered to.

Design the emergence profile gently. Any further expansion of the gingiva than the expansion obtained with the gingiva former should always be discussed with the attending dentist beforehand. The gingiva can be expanded gradually by individually designed temporary abutments.

In general: The circular notch should be positioned slightly towards the subgingival direction in the vestibular area. This way, the cement between coping and abutment is easy to remove and the transition is not visible from vestibular direction. The minimal wall thickness of the abutments should not fall below 0.5 mm (screw hole to the outer surface of the abutment). The diameter of the emergence profile should not exceed 13 mm, while the height of the emergence profile should not exceed 6 mm. The total length of the PreFace Abutment should not exceed 15 mm, the maximum angle of the PreFace Abutment should not exceed 30°. The taper of the PreFace Abutment should not exceed 6°. The PreFace Abutment should be rounded occlusally and should taper off. Avoid sharp edges.

Finishing:

In case that further work steps are necessary after working on the PreFace Abutment in the CAD/CAM milling machine, please note the following:

To mill PreFace Abutments, use appropriate, flawless mills with only little pressure. The minimal wall thickness must not fall below 0.5 mm. Avoid ridges and edges. We recommend to fix the PreFace Abutments in the lab analog to protect the connection geometry.

Scan body:

A scan abutment is used to determine the exact position of the implant during the scanning procedure. The scan abutment's diameter must correspond to the diameter of the selected implant.

* PreFace is a registered trademark of Medentika GmbH

IMPLA surgery kit incl. surgery tray

Purpose of the surgery tray:

Storage and fixation of the surgery and prosthetic instruments which are required for dental implantation.

The instruments are to be stored on the surgery tray in a sterilization box. They may also be sterilized in this sterilization box (with any sterilization procedure conventionally used in a dental practice).

Furthermore, the instruments are arranged in the box as chronologically needed by the oral surgeon.

Setup of the surgery tray:

The surgical instruments are arranged in their designated places with their functional side up. The preset color code of the implants makes it possible to arrange and to rearrange the instruments chronologically.

The goal is to arrange the instruments in accordance with the conventional course of the surgery, i. e. to arrange the instruments conveniently for the oral surgeon from left to right.

IMPLA surgery kit incl. surgery tray

IMPLA 3D Navi surgery tray

Delivery systems (cases and trays) are intended for the protection, organization and the delivery to the surgical field of surgical instruments and/or other medical devices. Delivery systems (cases and trays) are not designed to maintain sterility by themselves. They are designed to facilitate the sterilization process when used in conjunction with a wrapping material or a specified filtered sterilization container system. Wrapping materials and sterilization containers are designed to allow air removal, steam penetration / evacuation (drying) and maintain the sterility of the internal components.

Recommendations for Care, Cleaning and Sterilization

Both physical and chemical (detergent) processes may be necessary to clean soiled instruments. Chemical (detergent) cleaners alone cannot remove all soil and residues; therefore, a careful manual cleaning of each item with a soft sponge or cloth is essential for maximum decontamination. For difficult access areas, a clean soft bristled brush is recommended. Once the items have been cleaned, they should be rinsed thoroughly with clean water to remove any detergent or chemical residue before sterilization. We recommend the use of a mild enzymatic detergent with a near neutral pH. Do not use solvents, abrasive cleaners, metal brushes or abrasive pads. Cases and trays may be placed in mechanical cleaning equipment.

Always inspect for cleanliness or damage before use. Make sure all latches and handles are secure and in working order. Do not overload cases and trays. Balance contents uniformly within the container and arrange to allow steam to come in contact with all objects in the container. We also recommend that the delivery systems (cases and trays) be processed according to the sterilization wrap or filtered container manufacturer's instructions prior to sterilization to maintain sterility of internal components/ items and for proper aseptic presentation to the surgical field.

Our sterilization box (cases and trays) is suited to steam sterilization of conventional surgical instruments, such as reamers, handpieces, pliers, blades and drills, including short-lumen instruments.

Complex instruments may require disassembly and extended sterilization times. Always follow instrument manufacturer instructions if their sterilization or drying recommendations exceed these guidelines.

1. Prevacuum Sterilizer

Wrapped cases, trays and instruments, or cases, trays and instruments should be exposed to 132°C to 135°C for at least four minutes. Dry for 20 to 40 minutes.

2. Gravity Displacement Sterilizer

Wrapped cases, trays and instruments should be exposed to 132°C to 135°C (270°F to 275°F) for at least 30 minutes, or 121°C to 123°C (250°F to 254°F) for at least 55 minutes. Dry for 20 to 50 minutes.

Variables that may affect drying times include: loading density of the case/tray, instrument configuration, total contents of the sterilizer, steam quality, equipment maintenance and others.

Contraindications

Delivery systems (cases and trays) have **NOT** been validated for use with flexible endoscopes, or devices with lumens or working channels longer than 4 inches (3 mm internal diameter). Always refer to instrument manufacturer instructions.

Stacking of delivery systems (cases and trays) and overloading of the units will adversely effect sterilization and drying effectiveness.

DO NOT STACK cases and trays in autoclave chamber.

Notice:

1. **DO NOT** load cases into sterilizer on their sides or upside down with lid side on the shelf or cart. Always load cases on cart or shelf, with the lid facing upward. This will allow for proper drying. Cases and trays are designed to drain in this position.
2. After the autoclave door is opened, all cases must be allowed to cool thoroughly. Place cases on a rack or shelf with linen cover until cooling is complete. The potential for condensation may increase if the case is not allowed to cool properly.
3. If condensation is observed, check to ensure steps 1 and 2 were followed. In addition, verify that the steam, which is being used for sterilization processing, has a quality of more than 97%. Furthermore, confirm that the sterilizers have been inspected for routine maintenance in accordance with manufacturer recommendations.
4. Silicone mat accessories should be placed in alignment with perforations in cases or trays so that sterilant pathways are not obstructed.
5. Only use internal baskets, trays, mats, or accessories that have been designed and tested for use in ETO processes. Trays and mats are made from metal, Radel®-R or silicone and meet these requirements.

The following table provides guidance on maximum loads. Do not overload delivery systems (cases and trays), and always follow standards for maximum loading. Follow instrument manufacturers' instructions for use.

Disinfection, cleaning and sterilization of surgical and prosthetic tools and instruments of the IMPLA system

All surgical and prosthetic instruments/tools of the IMPLA system are composed of one of the following materials:

Material groups:

1. surgical steel / titanium

x-ray spheres, drills, drill extensions, countersinks, bone thread cutters, ratchet, screw drivers, IPK, depth stop, gingiva cutter, contra angle, flat wrench, sleeve holder, fittings of IMPLA 3D sterilization box, carrier and guide sleeves, depth gauge, tweezers, sealing and healing cap, gingiva former, abutments, impression aid, horizontal rosette, vertical screw, bar screw, bone screws

2. aluminium

IMPLA 3D sterilization box

	Surgical steel / titanium	Aluminum
Please observe!	Neutralize and dry the instruments thoroughly in order to avoid precipitation on the instruments to be sterilized. It is recommended to autoclave unmixed materials only.	Some components of the IMPLA system consist of anodized aluminum. The anodized surface increases the abrasion resistance as well as the wear resistance. However, the anodized surface is only resistant in the pH neutral range and is affected by e. g. alkaline cleaning agents and acid neutralization agent, the colors fade or become spotted.
Disinfection	Always wear gloves when cleaning the components. Avoid drying of sticking residues of blood and/or tissue. Always use an appropriate disinfection bath and dose according to manufacturer instruction.	Always wear gloves when cleaning the components. Avoid drying of sticking residues of blood and/or tissue. Always use an appropriate disinfection bath and dose according to manufacturer instruction.
Cleaning	Stainless steels are to be cleaned uncritically by means of all cleaning agents suitable for the cleaning of surgical instruments and are resistant to any kind of cleaning process (e. g. ultrasonic bath). Please avoid the use of metal brushes and abrasives. After cleaning, please clean with pH neutral cleaning agent and then flush away cleaning agent and particle residues with demineralized water and dry immediately.	The cleaning facilities must be suitable for the respective cleaning. Do not use dishwashers or ultrasonic devices for the cleaning of aluminum parts. Only use neutral cleaning agents. Please observe manufacturer instructions, do not use abrasives or metal brushes. Containers must be opened for cleaning. After cleaning, please clean with pH neutral cleaning agent and then flush away of cleaning agent and particle residues with demineralized water and dry immediately.
Sterilization	Surgical steels are resistant to any kind of sterilization process. Make sure that all deposits have been removed by cleaning. The parts may not have any contact during sterilization. The temperature should not exceed 176°C (248°F). Please observe the instructions of the autoclave manufacturer.	Make sure that all deposits have been removed by the cleaning. The parts may not have any contact during sterilization. Please observe the instructions of the autoclave manufacturer. The following sterilization methods may be used: Sterilization in the autoclave at 134°C (273°F), in the chemoclave at 132°C (270°F) and sterilization with dry heat at 180°C (356°F) – please observe the respective cycle times as stated in the manufacturer instruction.

Instructions for the practice:

1. Do not let surgical residues (blood, tissue, etc.) dry on the instruments/parts. Submerge the instruments/parts in disinfective agent directly after the surgery.
2. Clean the instruments/parts only with a nylon brush. (Do not use metal brushes, such as brass brushes, etc.)
3. Rinse off thoroughly with clear water.
4. Dry instruments/parts directly after cleaning and rinsing.
5. Always sterilize instruments/parts (e. g. ratchet) made of different materials separately!
6. Please note the instructions for use of each cleaning agent, as such agents are not always suitable for all materials (see material groups above).

1 Use of the torque ratchet:

The torque ratchet is build up of the ratchet head (Fig. 1-b), the ratchet handle (Fig. 1-c) and the setting screw (Fig. 1-e). Inside the ratchet head, there is a tool fitting (Fig. 1-a).


Fig. 1

The ratchet head is lockable to lock the torque limit (see 1.1 below). If the ratchet head is not locked, the torque can be adjusted infinitely between 10 and 30 Ncm with the setting screw (see 1.2).

The ratchet must be disassembled, cleaned, dried and **immediately before its use, it must be sterilized** (see 3).


Fig. 2


Fig. 3

Fig. 4


Fig. 5

1.1 Adjusting the ratchet

In the factory setting, the ratchet head is not locked. If you would like to lock it, turn the setting screw to less than 10 Ncm (Fig. 4), to relieve all pressure on the spring.

Subsequently, pull the handle from the ratchet head (Fig. 2). Then, turn the handle clockwise by 90° and let it click into place in the ratchet head. Finally, preset a torque of 20 Ncm to ensure that the ratchet remains locked.

Attention: In the fixed condition, the maximum load is 100 Ncm!

1.2 Setting the torque function

First, ensure that the ratchet is not locked. Turn the setting screw until you have reached the necessary torque. The correct torque is set once the according marking ring is just still visible (Fig. 4). The ratchet head will now bend once the preset torque has been reached (Fig. 5). Afterwards, do not keep on turning, as this could lead to exceeding the preset torque significantly.

Important notice:

Whenever the ratchet is not in use, turn back the setting screw to relieve the pressure from the spring.

1.3 Change tool fitting


(See 2.1 [1.] and 2.2 [5.]

2 Maintenance of the torque ratchet

2.1 Disassembly

1. The tool fitting (Fig. 6-a) is removed from the ratchet head (Fig. 6-b), by pulling on the two pin ends (arrows) towards the ratchet handle. Please note that the tool fitting can only be pulled from the ratchet head towards one side. **Do not use force.**
2. Completely screw off the setting screw (Fig. 6-e) from the ratchet handle (Fig. 6-c).
3. Pull off the ratchet handle (Fig. 6-c) from the ratchet head (Fig. 6-b) and remove the spring.
4. Disinfect, clean and dry all parts (see 3 – reconditioning).

Afterwards and **before the sterilisation**, the ratchet is reassembled.

2.2 Reassembly

1. Before reassembling the ratchet, all parts must be disinfected, cleaned and dried completely. Cover the areas marked with ● with oil (see 3.4).
2. Push the ratchet handle (Fig. 6-c) over the threaded bolt into the ratchet head (Fig. 6-b).
3. Push the spring (Fig. 6-d) into the ratchet handle (Fig. 6-c).
4. Screw the setting screw into the ratchet handle (Fig. 6-e) until the resistance of the spring becomes evident.
5. The tool fitting (Fig. 6-a) is set back into the ratchet head (Fig. 6-b), by pulling the two pin ends (arrows) towards the ratchet handle. Please consider that the tool fitting can only be set into the ratchet head from one side. **Do not use force**

Important information:

Always relieve pressure on the spring when ratchet is not in use. Screw off the setting screw until the pressure on the spring is relieved completely.

3 Reconditioning

3.1 Disinfection

Before first use and immediately after each use, all ratchet parts must be disinfected in a suitable disinfectant. Adhere to the instructions of the disinfectant's manufacturer.

Chemical as well as thermal disinfection techniques are suitable, if they come up to the requirements of medical instrument disinfection.

3.2 Cleaning

Immediately after disinfection, all parts of the ratchet must be cleaned with a suitable cleaning agent. Adhere to the instructions of the cleaning agent's manufacturer. Otherwise, impurities will dry off onto the instrument and make cleaning difficult.

In order to avoid corrosion and discoloration, closely adhere to the instructions on concentration, submergence time and temperature given by the manufacturer of the cleaning agent.

Remove any coating, impurification, etc. thoroughly with a soft brush (do not use metal brushes, scouring agents, wire wool or too much pressure).

Please note:

When disinfecting or cleaning the device with an ultrasound unit, please adhere to the instructions of the unit's manufacturer.

3.3 Rinsing and drying

Directly after disinfection, cleaning or ultrasound treatment, the separate parts must be rinsed carefully with clear desalinated or distilled water (do not use tap water). Afterwards, dry all parts thoroughly to avoid water stains.

If the separate parts are not rinsed, residue is not removed and/or the parts are not dried completely, this could lead to discoloration or corrosion of the device during sterilization.

3.4 Care and checking

Parts that have joints must be treated with a sterilizable oil spray (paraffin-based) after all cleaning and disinfection procedures. Such a spray will reduce metal-on-metal-friction, avoids fretting corrosion and reduces the danger of oxidation during heavy use.

3.5 Sterilization

Only sterilize carefully cleaned, dried devices without any type of corrosion.

Badly cleaned and insufficiently dried parts bear the danger of burning-in or corrosion damage.

If any corroded devices are sterilized with immaculate instruments, flash rust can cause damages on the immaculate instruments.

Please note the instructions of the manufacturer.

Important note:

Only sterilize the torque ratchet directly before use.


3.6 Storage

Important information regarding the storage of the torque ratchet:

After use of the ratchet the spring shall be released completely to guarantee proper operation and durability. The ratchet must be stored in released condition.

Cleaning and sterilization of the ratchet

Screw off turn-knob, remove spring and handle.

To clean, pull back the interlock to unlock the receptacle

Sterilization of the ratched only in disassembled state, without the spring.
Attention: Autoclave springs eparately!

User and safety instructions

All rotating instruments and parts for medical use (dental implantology, surgery) were developed for their specific application. For that reason, incorrect use may cause premature wear of the instruments and endanger the patients.

All rotating instruments for medical use should be deployed exclusively by medical personnel or other specialists who are qualified through their professional training and experience in the practical use of these instruments.

Correct application

- Ensure that only technically and hygienically well-serviced contra angles are used with the rotating instruments.
- The instruments must be inserted correctly.
- The instrument must revolve at operating speed before being placed on the bone surface (see page 4).
- Avoid wedging or levering the instrument as this may increase the risk of fracture.
- Thermal damage through the use of rotating instruments must be avoided (we recommend the slow-speed mode with sufficient cooling, please see page 4).

Pressing forces

- Excessive pressing forces must be avoided at all times. Excessive pressing forces may damage the working sections of the instruments or fracture their blades. This also results in an increased generation of heat.

Cooling

- To prevent excessive heat generation during preparation, provide for adequate air or water cooling (at least 50 ml/min).
- Inadequate cooling may lead to irreversible damages of the bone and/or tissue.

Manufacturer's information on rotating surgical instruments

Products:

The present manufacturer's information applies to all instruments that are used for surgical application. These include rotating instruments made of stainless steel.

Instruments delivered in non-sterile condition must be prepared prior to first use.

Limited number of reprocessing cycles:

The end of a product's service life depends on its degree of damage and wear. Do not exceed limitations of the permitted frequency of use in certain instruments, if these are known. The manufacturer confirms that the above detailed reprocessing methods are suitable for preparing the above named instrument group to enable their reuse (please see page 4). Frequent reprocessing does not affect the performance of these instruments.

Work station:

Hygienic precautions according to the provisions valid in your country.

Storage and transport:

Place instruments in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). It is recommended to process the instruments within one hour of use at the very latest. The instruments should be in the cleaning/disinfection tank (Fräsator) when transported to the site where the reprocessing is to take place. Endodontic instruments can be transported in a special interim support equipped with a foam insert drenched in disinfecting solution.

Standardized manual reprocessing for rotating instruments

Equipment used:

- Nylon brush
- Suitable detergent/disinfectant for rotary instruments with proven disinfecting effect
- Ultrasonic bath

Reprocessing:

1. Rinse off surface contamination under running water. Remove stubborn contamination with a nylon brush under running water, turning the instrument constantly.
2. Rinse instrument thoroughly with running water.
3. Place the instruments in a suitable sieve into the ultrasonic device filled with detergent/disinfectant.
4. Observe the instructions of the manufacturer regarding concentration and immersion time for chemical cleaning/disinfection in the ultrasonic device. Be sure to observe the full correct immersion time which does not start until the last instrument has been placed into the ultrasonic device. Attention: Do not exceed 45°C (113°F) (risk of protein coagulation)!
5. On completion of the immersion time, rinse instruments thoroughly with suitable water (preferably with demineralised water to avoid residues).
6. Dry instruments (preferably with compressed air).
7. Visual examination with a suitable magnifying device to ensure that the instrument is clean and undamaged (experience has shown that a magnification factor of 8 permits a visual examination). If there are still visible residues of contamination, repeat the cleaning and chemical disinfecting process until no visible contamination is left.

Control and functional test:

Instruments showing the following defects are to be discarded immediately:

- Blunt and chipped blades
- Deformations (e. g. bent/twisted/fractured instruments)
- Corroded surfaces

Packing:

Make sure that the packaging is suitable for the instrument and the chosen method of sterilization. Single pack: The packaging must be large enough to ensure that there is no pressure on the seal. In the set: Place instruments onto the tray provided or onto universal sterilization trays. The instruments must be protected. Use an appropriate method to pack the tray. Instruments with limited times of use are to be marked accordingly.

Sterilization:

Steam sterilization using a fractionated vacuum process at 134°C (273°F) in a device that complies with the provisions of DIN EN 13060; with validated processes.

- Fractionated pre-vacuum (type B) respectively simple pre-vacuum (type S).
- Sterilization temperature: 134°C (273°F)
- Hold time: 5 minutes minimum (full cycle)
- Drying time: 10 minutes

In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer.

Transport and storage:

The packed sterile goods must be protected from dust, moisture and recontamination during transport and storage.

Universally valid notes:

Observe the legal provisions regarding the reprocessing of medical products valid in your country. The manufacturer confirms that the above detailed reprocessing methods are suitable for preparing the above named instrument group to enable their reuse. The person carrying out the reprocessing is responsible for ensuring that the applied method is carried out with appropriate equipment, materials and trained personnel at the reprocessing site and that it actually achieves the desired result. To guarantee this, routine controls of the validated mechanical and/or manual preparation methods are necessary. Any deviation from the above detailed process (e. g. use of different chemicals) must be carefully checked by the operator to ensure effectiveness and to avoid possible adverse consequences.

Information on storage, disinfection, cleaning procedures and sterilization of rotating instruments

- Rotating instruments must be disinfected, cleaned and – whenever the instrument material permits it and it seems suitable – sterilized prior to first use on patients and immediately after use.
- Used instruments must always be disinfected before cleaning.
- Before first use, the rotating instrument should be kept in its original packaging and be stored at room temperature, protected against dust and moisture.
- All instruments and parts are packaged non-sterile and must be sterilized before use according to their application. Cleaning and sterilisation not in accordance with the instruction may lead to microbiological contamination and infection of the patient.
- All materials are sterilizable according to EN 554 (autoclaving, e.g. at 134°C (273°F) for a minimum time of 3 min. or at 121°C (250°F) for 15 min.). Please observe manufacturer's instructions with regard to the autoclave.
- Before sterilization, the instruments must be cleaned and disinfected properly.
- Disinfection may be carried out manually or by machine.
- In case of severe contamination of the instrument, ultrasound cleaning techniques are recommended.
- After cleaning rinse the instruments thoroughly with water and dry them immediately.
- Contact between the instruments during cleaning and disinfection procedures may cause damage and should be avoided.
- It is highly recommended to wear protective gloves when working with contaminated instruments.
- Thermal disinfection procedures are generally not suitable for rotating instruments.
- Avoid sterilization temperatures exceeding 180°C (356°F). Exceeding this temperature reduces the hardness of the working section and reduces the instrument's service life.
- All sterilized instrument containers should be labeled with the date of the sterilization as well as the expiry date.
- Non-rustproof instruments may corrode. This may cause discolorations and reduce the service life. Therefore, it is recommended to use anti-corrosion disinfectants and cleaning agents in these cases.
- The instructions for use, reaction time, and suitability of disinfectants and cleaning agents for certain types of instruments are provided by the manufacturer's instructions. They must be observed.
- Contact of the instruments with hydrogen peroxide (H₂O₂) must be avoided.
- The instruments should be stored in appropriate, hygienically well-kept containers. This also applies to sterilized instruments. When stored, the instruments should be protected against dust, humidity, and recontamination. Observe the maximum storage period.
- Cracked or broken cutting edges cause vibrations and high pressing forces which may lead to broken preparation edges and rough surfaces and should be discarded out immediately.
- Bent or non-concentric instruments should be sorted out immediately and discarded.

Annex:

Forms and information sheets for patients

„Medical examination for a dental implantation“	43–44
„Informed Consent“	45–47
„Information Sheet“	48
„Surgery Log File“	49–50
„Implantological aftercare in terms of a recall“	51
„Instructions“	52–53

MEDICAL EXAMINATION FOR A DENTAL IMPLANTATION

Date / Place _____

Dear Colleague,

My patient, Mr./Mrs. _____ will undergo a jawbone surgery. This surgery will include the implantation of a bio-inertial implant into the jawbone.

I would like to ask you for a detailed examination of the patient in order to execute the following diagnostic findings and exclude certain diseases.

	Yes	No
- Frequent doctor's appointments	<input type="checkbox"/>	<input type="checkbox"/>
- Heavy nicotine consumption	<input type="checkbox"/>	<input type="checkbox"/>
- Heavy medication		
Corticoid	<input type="checkbox"/>	<input type="checkbox"/>
Cytostatics	<input type="checkbox"/>	<input type="checkbox"/>
Marcumar	<input type="checkbox"/>	<input type="checkbox"/>
Others	<input type="checkbox"/>	<input type="checkbox"/>
- Neurological diseases	<input type="checkbox"/>	<input type="checkbox"/>
- Bone diseases (Morbus Paget, Osteoporosis, Recklinghausen)		
Frequent bone fractures	<input type="checkbox"/>	<input type="checkbox"/>
- Cardiovascular diseases, if yes which type:	<input type="checkbox"/>	<input type="checkbox"/>
- Cardiac arrhythmia	<input type="checkbox"/>	<input type="checkbox"/>
- Rheumatism	<input type="checkbox"/>	<input type="checkbox"/>
- Diabetes	<input type="checkbox"/>	<input type="checkbox"/>
- Thyroid diseases	<input type="checkbox"/>	<input type="checkbox"/>

- | | | |
|--|--------------------------|--------------------------|
| - Metabolic diseases | <input type="checkbox"/> | <input type="checkbox"/> |
| Allergies | <input type="checkbox"/> | <input type="checkbox"/> |
| Allergic to Penicillin (antibiotics) | <input type="checkbox"/> | <input type="checkbox"/> |
| Allergic to Procaine | <input type="checkbox"/> | <input type="checkbox"/> |
| Allergic to metals | <input type="checkbox"/> | <input type="checkbox"/> |
| - Other infectious diseases | <input type="checkbox"/> | <input type="checkbox"/> |
| - Aberrance of the normal hemogram | <input type="checkbox"/> | <input type="checkbox"/> |
| Differential hemogram | <input type="checkbox"/> | <input type="checkbox"/> |
| ESR | <input type="checkbox"/> | <input type="checkbox"/> |
| Hepatitis search results | <input type="checkbox"/> | <input type="checkbox"/> |
| - Are there general objections against an implantation surgery for this patient? | <input type="checkbox"/> | <input type="checkbox"/> |

Date / Place, Signature of the attending physician

The decision and the responsibility for this implantation surgery remain with the executing dentist / implantologist.

Please let me have the signed copy of this report back at your earliest convenience.

Many thanks and best regards

INFORMED CONSENT

1. I have been informed about the nature and the method of the implantation and I understand the surgical procedures. It has been explained to me that an implant has to be inserted into the bone, the alveolar ridge or the tooth to act as fixation.

2. All alternative therapies of the dental reconstruction have been explained to me. My dentist has carefully examined my mouth. After a discussion with the dentist, I decided on the implantation.

3. I have been comprehensively informed that possible risks or complications could occur with every surgical measurement, every local anesthesia, or in connection with additional medication.

I am informed about possible subsequent complications such as swelling, pain, infections or a general state of indisposition. Numbness of the lip, tongue, cheek, chin or teeth could also occur. The exact duration of this impairment cannot be predicted – in some cases such impairment could also become irreversible. Additional infections of the wound, the vascular system or the surrounding tissue are possible. I have been informed that bone fractures, intrusions of the maxillary sinus, delayed healing or allergic reactions due to drugs or anesthetics could occur.

4. I am aware that I have to inform the dentist immediately if there should occur any problems which go beyond the normally experienced issues.

5. I have been informed about the fact that until today there are no methods to precisely predict the healing chances of the bone and the alveolar ridge before the surgery. The healing process is individually different. However, the further medical treatment is dependent on the healing process.

6. I will receive a one-/two-phased implant. A one-phased implant means that within a short time after the surgery or immediately after the surgery, a provisional sealing will be cemented over the implant. After some 6 – 8 weeks, the additional treatment will be effected. With a two-phased implant, the healing period will last some 3 – 4 months – in rare cases even much longer than that. During this period of time, the implant will be left as it is covered or mostly covered by the mucous membrane.

7. I have been informed that there is no success guarantee for implants. In case of a failure, the implant must be removed immediately. The timing for the removal is determined by the dentist.

8. I have been informed that maximum oral hygiene must be applied especially around the implants. I will support the success of the implantation by applying this maximum oral hygiene.

9. I have been informed that excessive smoking and alcohol abuse as well as parafunctional exposures of my implants could massively influence the success of the surgery. I promise to follow the instructions of my dentist and agree to a regular check-up twice a year.

10. I have informed my dentist about my dental and medical background and also received an additional questionnaire for my doctor. I will instruct the doctor to hand the filled in questionnaire back to the dentist.

11. I do agree to x-ray pictures and photographs during the surgery as well as during subsequent examinations.

12. My dentist has explained to me that implantological treatments are private health care services which are also invoiced as such (patient payment). Cutbacks of the charges on behalf of the health insurance or the insurance are not to be borne by the dentist but will be borne by myself.

13. Cost of subsequent examinations or the exchange of implants/implant parts or changes in the overall implant construction must be covered by myself, notwithstanding the possible support payments of the health insurance.

14. I have been informed about the appraiser processes which apply for implantological services. My dentist informed me that he will answer all possible questions which might normally occur during such an appraisal process. Possibly cutbacks of the charges due to a reduced financial support by the health insurance will not be effected by me without prior consulting with the dentist.

Date, Place _____

Signature of the physician

Signature of the patient

Signature of the legal guardian

Signature of the assistant

INFORMATION SHEET

Guideline for patient counseling

Remarks:

1. What is an implant?

Material, methods,
own experience, anesthesia

2. Prosthetic possibilities

With implant
Alternative treatment – without implant

3. Indication / contraindication

See evaluation of the Medical History Sheet

4. Success – failure – risk

Cannot be guaranteed!
Paresis: paresthesia etc.

5. Oral hygiene

Conditio sine qua no
Motivation
Recall: check up

6. Billing

The implantation is not covered by the health insurance.

Date / Place _____

Stamp

Signature of the dentist

Signature of the assistant

SURGERY LOG FILE

Patient Name	:	
Surgery Date	:	
Begin of Surgery	:	
End of Surgery	:	
Pre-Medication	:	
Blood Pressure / Pulse (pre)	:	
(intra)	:	
(post)	:	
ANESTHESIA		
Local Anesthesia	:	<input type="radio"/>
Infiltration Anesthesia	:	<input type="radio"/>
Nerve Block Anesthesia	:	<input type="radio"/>
Endotracheal Anesthesia	:	<input type="radio"/>
MUCOPERIOSTEAL FLAP OPERATION (UPPER JAW)	:	<input type="radio"/>
Palatinal	:	<input type="radio"/>
Vestibular	:	<input type="radio"/>
Mid of alveolar ridge	:	<input type="radio"/>
MUCOPERIOSTEAL FLAP OPERATION (LOWER JAW)	:	<input type="radio"/>
Lingual	:	<input type="radio"/>
Vestibular	:	<input type="radio"/>
Mid of alveolar ridge	:	<input type="radio"/>
INSERTED IMPLANTS		
As planned	:	<input type="radio"/>
Differently	:	<input type="radio"/>
(Fill in the implants and augmentations)		
DEFECT FILLING		
Hydroxylapatite	:	<input type="radio"/>
Autologous bones	:	<input type="radio"/>
Homologous bones	:	<input type="radio"/>
Bone chip extraction	:	<input type="radio"/>
DEFECT COVERING		
Mucosa extraction	:	<input type="radio"/>
Mucosa membrane transplant	:	<input type="radio"/>
Membrane technique	:	<input type="radio"/>
NUMBER OF SUTURES		
Silk Suture	:	<input type="radio"/>
Catgut Suture	:	<input type="radio"/>
Tissue Glue	:	<input type="radio"/>

- UNFORESEEN OCCURANCES :
- Maxillary antrum perforation :
- Mouth nose connection :
- Bone perforation (lower jaw) :
- Nerve channel perforation :
- Bone chippings :
- Bone fracture :
- Breakage of instrument :
- Soft tissue perforation :
- Vessel rupture :
- Bleeding cannot be stopped :
- Aspiration :
- Others :

Remarks:

Physician

Assistant

IMPLANTOLOGICAL AFTERCARE IN TERMS OF A RECALL

The patient in the center of the dental treatment is still in need of treatment after a successfully executed surgical or prosthetical care.

The post-operative treatment duty does not end with the insertion of enossal implants, as it seems to be possible by the separation of the individual areas of expertise (surgical procedures, parodontological treatment and prosthetical care), but it should be executed by the physician in charge of the final treatment meaning the physician who executed the last treatment to complete the reconstruction. This treatment also includes the monitoring of the treatments and prosthetical care.

Should the aftercare mean an overstrain with regard to the capabilities of the last physician in charge, he has to ensure that the patient is being transferred to the corresponding specialist treatment.

The intervals of the regular check-ups are defined in the parodontological guidelines and will be performed accordingly. A check-up every six (6) months is generally sufficient. However, the process requires an individual consulting with the patient. For example, in case that the patient is not able to perform the required maximum care for the implants or the prosthetical construction due to physical handicaps, shorter check-up intervals could be required.

It is an inevitable standard that prior to any surgical or prosthetical treatments preventive and prophylactic measurements have to be applied. Of course any such instruction will be documented in the corresponding files. In addition to that, this education falls under the number GOZ 100 [German scale of dentist fees] which is the standard position for the health insurances to reimburse the costs according to the guidelines of the private health insurance (PKV). If the education is executed without the application of the GOZ 100, it could happen that the insurance company claims a lack of evidence and this might lead to a non-reimbursement of the costs.

During the check-up examination, it is always necessary to find out if an increase of mobility is achieved by applying a periimplantological or osteolytical or any other process.

When ordering a patient to a recall, it is of forensic importance to get the signed Informed Consent form.

INSTRUCTIONS

for the behavior after an

IMPLANTATION, TRANSPLANTATION AND REIMPLANTATION, AUGMENTATION

1. On the day of the surgery as well as the following two days, no smoking and no excessive alcohol consumption is permitted.
2. No physical strain.
3. Only when the local anesthesia is gone, eating and drinking are allowed.
4. Only have liquid and soft food – if possible, avoid dairy products as well as food which contains flour.
5. Do not chew in the area of the implant.
6. Apply maximum oral care to clean the other areas of the mouth and the jaw.
7. Only gently rinse the mouth in the area of the wound. During the healing process, no electric mouth shower or electric toothbrush should be used. After any meal shortly rinse the mouth with lukewarm water.
8. Do not touch the implant area and also avoid touching the stitches with the tongue.
9. In case of secondary bleedings use a clean paper tissue or gauze bandage and gently bite on it. In case the bleeding should not stop within 30 min., please call the dentist right away.
10. Swellings in the area of the surgery, the cheek and the chin are normal and should disappear within a few days. Cooling from the outside will help.

Behavior after the stitches have been taken out:

1. The wound should be rinsed gently and cleaned with a soft toothbrush. The better the wound has healed the stronger and more intense the toothbrush can be used to massage the surgery area.
2. In case of pain or any other unusual symptoms (like bleeding of the alveolar ridge or swellings), please see the dentist immediately.

3. If you are not sure how to clean the teeth, please ask the dentist or any of his assistants.

4. After the prosthetical care – which will start after three months at the earliest – the implant must be cleaned with extreme care.

5. A check-up of the implant is required every six (6) weeks

6. In case of any additional questions, please contact your dentist or his assistants.

Date / Place _____

Signature of the physician

Signature of the patient

Signature of the legal guardian

Signature of the assistant