Capo Bond

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Instructions Capo Bond

PART 1: Introduction

Capo Bond is a light curing, one part primer and bonding system which does not require additional mixing. Bonding is carried out in a single stage. Capo Bond was developed for firmly bonding composites and compomers to enamel and dentine as well as non-precious and precious metals. Capo Bond is also indicated for impregnating root canals prior to filling or cementing passive or threaded root posts.

Capo Bond can also be used as a primer for indirect adhesive restorations, e.g. porcelain and composite inlays, onlays, laminate veneers, crowns and bridges placed with chemically or dual-curing cement.

The high bond strength between Capo Bond and enamel or dentin is based on the same principles as apply to glass ionomer cements. The methacrylic carboxylic acid ester provides for long term high bond strength and biocompatibility. In case of irritation or methacrylate allergy, do not use this material.

Capo Bond is compatible with all commonly available brands of light curing composite. As Capo Bond is based on ethanol it is hydrophilic and it can be applied to slightly moist dentin surfaces using the "Wet Bonding" technique.

In cases where a dual-curing or self-curing priming and bonding system is recommended, Capo Bond can be dual-cured using a 1:1 mixture of Capo Bond and Capo Bond Activator. Capo Bond Activator is available separately. This special type of application is described in the Capo Bond Activator instructions.

PART 2: Using Capo Bond for filling materials

Preparation for light curing composite and compomer filling materials

1.1 Isoliation

It is advisable to use a rubber dam.

1.2 Cavity preparation

Before beginning preparation, clean the tooth with pumice and water. To ensure that as little healthy tooth structure as possible is reduced, only prepare the tooth minimally. Enamel cavity margins must be bevelled slightly (0.5–1.0 mm) to enlarge the adhesive surface in the enamel and increase the bond strength.

1.3 Pulp protection

In very deep cavities those areas in close proximity to the pulp must be coated with a thin layer of calcium hydroxide cavity liner.

1.4 Etching the enamel

Etch (e.g. with Capo Etch) according to the manufacturer's instructions. The etched enamel surface then appears chalky-white.

The etched enamel surface must not be touched or contaminated with saliva prior to applying the Capo Bond. Should it become contaminated, repeat the procedure, i.e. etch, rinse with water and dry again.

1.5 Applying Capo Bond

Shake the Capo Bond before use.

Use a brush to apply copious amounts of Capo Bond to the dentin and enamel and brush it in vigorously for 30 seconds. Then dry carefully with oil-free compressed air for approximately 15 seconds. Cure the Capo Bond coat for 20 seconds with a dental lamp before applying a second coat.

CAUTION: Do not spill Capo Bond! If Capo Bond is not used immediately, keep it in subdued light to prevent it from curing prematurely due to ambient light. Capo Bond will not cure spontaneously.

1.6 Applying the second coat of Capo Bond

Use a brush to apply a second coat of Capo Bond to the surfaces being adhered and brush it in vigorously for approximately 30 seconds. Then dry for approximately 15 seconds with oil-free compressed air and cure for 20 seconds with a dental lamp. The light curing composite filling can then be placed.

1.7 Placing, curing and finishing the filling

The Capo Hybrid instructions describe how the composite filling is placed, cured and finished.

2. Priming the root canals

2.1 Preparation

It is essential to isolate the teeth with a rubber dam. Clean the teeth with hydrogen peroxide solution (3 %) or sodium hypochlorite solution. When placing passive or threaded endodontic posts with luting composites, refer to the manufacturer's instructions. Prepare the canal to determine the length of the root post required. Adapt the post to the canal and reduce the length if necessary. Remove the post, keep it in a safe place and prepare the canal for cementing the post.

2.2 Priming the root canals

As the surfaces of the root canal cannot be reached with a brush, Capo Bond must be injected into the prepared canals with a syringe to rinse them for 60 seconds in four separate stages. Ensure that all surfaces being bonded are wetted continually. Wherever possible, use a brush to work the Capo Bond in and enhance wetting. Soak up the excess Capo Bond from the root canal using, for example, a soft paper point and dry the treated surfaces carefully for 15 seconds with oil-free compressed air. The Capo Bond is then cured for 20 seconds with a dental lamp.

2.3 Applying the luting composite

Place and cure the luting composite as described in its manufacturer's instructions.

PART 3: Special applications for Capo Bond

3. Using Capo Bond on precious/non-precious metals and porcelain for repairing crown and bridgework

3.1 Conditioning the surfaces of the restoration

Roughen the surfaces of the metal or fractured porcelain with a diamond instrument (40 mm) – if necessary, sandblast them in addition (with an intraoral sandblaster). Apply Capo Etch to the rough surfaces, rinse them with copious amounts of water and dry them with oil-free, dry compressed air.

3.2 Applying Capo Bond

Brush Capo Bond onto the metal surface carefully and dry it for 15 seconds with oil-free, dry compressed air. Cure the coat with a dental curing lamp for 20 seconds. Begin repairing the restoration.

PART 4: Storage and technical data

4. Storage

Store at 10–25 °C.

5. Technical data

Curing time (light polymerization) with a dental lamp: 20 sec. Use a light polymerization system with an emission range of 350 –500mm to polymerize the material. The required physical properties are only reached if the polymerization light functions properly. Therefore, it is necessary to check the light intensity regularly according to the manufacturer's instructions.Curing time for chemical curing ca. 3 min. (mixing ratio 1:1 of Capo Bond and Capo Bond Activator without light and under the exclusion of air).

6. Package contents:

Capo Bond set: 5 ml bottle Capo Bond, 1 dental brush holder, 50 single-use dental brushes, 2 mixing wells
Order No.: 234014

Accessories:
Capo Bond
Order No.: 234401
Capo Bond Activator
Order No.: 234402
1 Dental brush holder
Order No.: 230913
50 single-use dental brushes
Order No.: 230912
Mixing wells
Order No.: 234930

Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-indications / interactions

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and crosreactions of the product with other materials already in the patient's mouth before using the product.

Note

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

Error	Possible cause	Solution
Capo Bond does not cure	Light output of the polymerization lamp is insufficient	Check light output; replace light source, if necessary
Restoration does not stay attached	Working times, especially etching times were not kept <p>Individual layers were not polymerized separately</p> <p>Capo Bond was not brushed in intensively</p>	Adhere to the working times listed in the instructions <p>Polymerize each individual layer separately</p> <p>Adhere to the work steps listed in the application instructions</p>
Material is dried out or cured	Incorrect storage, bottle not closed tightly	Put the cap back onto the bottle after each use
Material does not create a homeogenous film	Capo Bond was applied unevenly or the layer thickness is insufficient	Apply fresh material once again

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