# Prothesen-Color-System

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## Instructions for use

### Indications

- characterization of denture base acrylics for partial and full dentures
- individualization of prefabricated acrylic teeth
- Characterization of composite veneers in the crown and bridge technique.

### Colors

The PC-System is available in the following colors:  
white, cream (off-white), ochre, apricot, orange, khaki, red, blue, brown, black

### Processing instructions

#### Characterization of denture base acrylics for partial and full dentures:

Form the denture body considering anatomical aspects and sand the edges. After polishing with pumice, polish those sections to a high shine, which are not to be painted.

Thoroughly clean the denture with an ultrasound cleaner or a steam cleaner. Please ensure that the sections to be characterized are not touched after cleaning.

The surface to be painted is moistened with Primer liquid and polymerized. In this process, a dispersion layer (sticky surface) must be created.

**Note: For other polymerization systems, please see the table at the end of these instructions.**

**Note: In case of a shiny, non-sticky surface, the selected polymerization time was too long. Surfaces that still show any pools of Primer liquid must be polymerized again.**

The necessary colors of the PC-System are put onto a mixing palette and are mixed individually, as needed.

After characterization, cure the denture for 3 min. in Spektra 2000. **An intermediate polymerization of each step is unnecessary.**

The dispersion layer may be removed with a soft tissue soaked with a bit of PCS-Cleaner.

Finally, polish to a high shine conventionally with polishing paste and linen buff.

Individualization of prefabricated acrylic teeth:  
The prefabricated teeth must be fashioned considering age, sex, and type of the patient.

After shaping the acrylic teeth, sandblast them with aluminium oxide (50 -125 µm). Clean the teeth as described in the section „Characterization of denture base acrylics for partial and full dentures“.

The acrylic teeth are moistened with the Primer liquid, which is cured for 1 min. in Spektra 2000.

The PCS-colors are applied to the dispersion layer thus produced and are then polymerized.

The removal of the dispersion layer and the high shine polish are carried out as described in the section „Characterization of denture base acrylics for partial and full dentures“.

#### Effects for crowns and bridges:

- The below steps are to be followed for liquid light-curing veneer materials but also for light-curing veneer materials with a higher viscosity:
- Once the opaquer material has been applied, the cervical, dentin and incisal materials may be characterized with the PCS colors.
  - Individualization of the veneer materials is obtained by mixing, inlaying, and painting.
  - The polymerization times correspond to those listed in the section „Characterization of denture base acrylics for partial and full dentures“ (see above).

### Guarantee

Our technical instructions, regardless of whether they are provided verbally, in writing or during practical demonstrations, are based on our own experience and should only be considered guidelines.

As our products are subject to continued development, we reserve the right to modify them.

### Intolerances

This product has been specially developed for the above mentioned indications and is only to be processed as described above. The manufacturer will not be held liable for any damages caused by improper handling or processing. The user is obligated to check before use, whether the product is adequate for use in the intended manner. A sensitization by the product may not be ruled out completely.

**The use of the product should be discontinued, if any reactions should occur or if there is a known intolerance against one or several of the components (e. g. methyl metacrylate).**

### Main components

Organic matrix  
methyl metacrylate  
Bis-GMA

### Additives

Thickener, initiators, stabilizers, pigments.

### Processing time of the colors

Up to 5 min., depending on the light conditions.

### Polymerization

The light curing unit used for polymerization must have an emission spectrum of 350 -500 nm.

**An adequate polymerization may only be obtained with properly functioning lamps. Therefore, it is necessary to check the light intensity regularly according to the instructions of the lamp manufacturer.**

**Shelf life**

The maximum shelf life is printed onto the label of each container. **Do not use after the shelf life has expired.**

**Storage**

Store in a well-ventilated area at 10 °C - 25 °C (50 °F - 71 °F). For best results, let the colors reach room temperature before use. **Shake the container well before use!**

**Keep this medical device away from children.**

Only for use in the dental technical sector.

**Immediately close the containers after use!**

Always keep containers closed to avoid unnecessary thickening of the material. The Primer liquid may be used as a thinner, if necessary.

**Avoid direct sun light.**

**Dangers**

**Avoid skin contact with the unpolymerized material and the inhalation of monomer vapors.**

**Please see the Material Safety Data Sheef for more information!**

**Packaging Size**

The PCS-set contains.

10 different colors (3.0 g each), PCS-Primer (two of 3.0 g each), PCS-Cleaner (50 ml).

**Polymerization of the PCS-colors in different light-curing units in the market**

light-curing unit	Polymerization times	
	PCS-Primer	PCS-colors
Spektra LED (Schütz Dental)	30 sec.	1 min.
Spektra 2000 / Spektra 2000 fast (Schütz Dental)	1 min.	3 min.
DC-XS / Uni-XS (Heraeus Kulzer)	180 sec.	180 sec.
Hilite power / HeraFlash (Heraeus Kulzer)	180 sec.	180 sec.
Spectramat (Ivoclar)	5 min.	5 min.
Labolight LV-II / III (GC)	5 min.	5 min.
Solidilite EX / V	5 min.	5 min.
Visio Beta Vario (Espe)	2 min. (under vacuum)	7 min. (without vacuum)

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