# FuturaJet universal

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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## Instructions for use

### Composition:

- Powder: bead polymer made of poly-methyl methacrylate (PMMA), pigments, initiators
- Liquid: methyl methacrylate, tetra-methylene dimethacrylate, initiators, stabilizers
- FuturaJet universal is free of cadmium.

### Indications:

- For the production of partial and full dentures in the injection technique or in the flask technique (UNIPRESS).
- For the production of the acrylic part of implant-supported restorations and on frameworks for dentures.
- For repairs and relinings of existing restorations.
- For the production of splints
  - occlusal splints
  - reflex / anti-grinding splints

### Contraindications / interactions with other agents:

Do not use in case of a known allergy against one of the components.

### Application / processing:

#### Preparatory steps:

Produce plaster model and wax-up as usual.

The prefabricated teeth set up in the wax are fixed with plaster or silicone (stent / flask). Subsequently, the wax is eliminated.

Roughen the acrylic teeth, provide them with undercuts and coat them with monomer or HS Cross-Liquid (Schütz Dental) to ensure that they are securely retained. The surfaces of the plaster models are insulated with the alginate separating agent FuturaSep Plus (Schütz Dental).

#### Mixing ratio:

2 parts powder by weight, 1 part liquid by weight.

45 cm<sup>3</sup> powder : 17 ml liquid (volume) or 30 g powder : 15 g liquid (weight) are required for an average partial denture.

Veins can be added, if desired. Please note the separate manufacturer's information and instructions.

### Production of partial and full dentures as well as of implant-supported restorations and framework completions in the flow technique or in the flask technique:

Mix powder and liquid (according to the mixing ratio) and fill the die as usual.

At room temperature (21 °C / 69 °F), FuturaJet universal is pourable for ca. 3 min. After 4 - 5 min., the material enters the plastic phase and after 7 - 8 min., at the very latest, the restoration must be placed in the pressure pot.

Polymerize FuturaJet universal in a pressure polymerization system for 15 min. at a pressure of 2 - 5 bar and a water temperature of 45 °C / 113 °F.

### Production of partial and full dentures as well as of implant-supported restorations and framework completions in the injection system UNIPRESS (please observe the separate instructions) (in case of other injection systems, please observe the manufacturer's instructions):

Mix powder and liquid (according to the mixing ration). After allowing a doughing time of 15 sek., stir FuturaJet universal with a spatula in a porcelain mixing bowl to produce a homogeneous mixture and fill into the UNIPRESS cylinder.

As soon as the material has reached a viscoelastic consistency, inject it into the prepared UNIPRESS flask.

These work steps must be completed after 8 min. at the very latest. Polymerization takes place under a defined pressure build-up inside the plate-and-spring system inside of UNIPRESS. The FuturaJet universal dentures can be deflasked after 20 - 30 min.

### Repairs and relining of existing appliances:

Produce a die or a stent, depending on the construction.

Clean the acrylic surfaces to be repaired, roughen them and smooth out the edges. Wet with monomer or with HS Cross-Liquid (Schütz Dental). Mix acrylic according to the mixing ratio and let the acrylic flow in.

Polymerize FuturaJet universal in a pressure polymerization system for 15 min. at a pressure of 2 - 5 bar and a water temperature of 45 °C / 113 °F.

#### **Polishing:**

Polish the material with goat hair brushes, polishing paste and soft woolen wheels. A careful surface treatment and polish are a precondition for an optimal result.

#### **Production of splints:**

##### **Preparatory steps:**

Production of a plaster model.

Blocking undercuts with suitable materials.

Carve splint in wax as desired.

Produce impression, fit it with sprue canals and remove wax.

Insulate plaster model with suitable materials (e. g. FuturaSep, Schütz Dental).

#### **Processing:**

Mix powder and liquid (according to the mixing ratio) and fill the mould as usual. Complete these worksteps by the time 7-8 min. have passed, at the very latest.

Polymerize in a pressure polymerization system for 15 min. at a pressure of 2 - 5 bar and a water temperature of 45 °C / 113 °F.

#### **Finishing the material:**

After polymerization, the splint is lifted from the model and it is finished with suitable tools.

Silicone polishers (wheel, lentil, roll) as well as hard metal cutters or diamond grinding tools are suited for finishing the material.

Prepolish with sandpaper or pumice.

#### **Polishing:**

Polish the material with goat hair brushes, polishing paste and soft woolen wheels. A careful surface treatment and polish are a precondition for an optimal result.

#### **Repairs:**

Produce a die or a stent, depending on the construction.

Clean the acrylic surfaces to be repaired, roughen them and smooth out the edges. Wet with monomer or with HS Cross-Liquid (Schütz Dental). Mix acrylic according to the mixing ratio and let the acrylic flow in.

Polymerize FuturaJet universal in a pressure polymerization system for 15 min. at a pressure of 2 - 5 bar and a water temperature of 45 °C / 113 °F.

## Troubleshooting / FAQ list

| Problem                                       | Cause                                                                                                                                                                                                                | Remedy                                                                                                                                                                                                                                           |
|-----------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| acrylic does not cure properly                | Incorrect mixing ratio                                                                                                                                                                                               | <ul style="list-style-type: none"> <li>– observe instructions for use</li> <li>– observe mixing ratio</li> </ul>                                                                                                                                 |
| acrylic shows a whitish tint along the bottom | <ul style="list-style-type: none"> <li>– incorrect or improper insulation</li> <li>– model was not soaked in water</li> <li>– first layer of insulation was already dry once the second layer was applied</li> </ul> | <ul style="list-style-type: none"> <li>– insulate properly</li> <li>– check insulation agent, if necessary</li> <li>– soak model in lukewarm water for ca. 10 min.</li> <li>– apply second layer of insulation on a moist first layer</li> </ul> |
| acrylic forms bubbles along the bottom        | <ul style="list-style-type: none"> <li>– model was not soaked in water</li> <li>– insulating agent too old or contaminated</li> </ul>                                                                                | <ul style="list-style-type: none"> <li>– soak model in lukewarm water for ca. 10 min.</li> <li>– use a new bottle of insulating agent</li> </ul>                                                                                                 |
| acrylic surface white / spongy                | too much time passed until the acrylic was injected or until the acrylic was put in the pressure pot                                                                                                                 | observe processing times                                                                                                                                                                                                                         |
| color deviations of repairs                   | incorrect mixing ratio                                                                                                                                                                                               | observe mixing ratio                                                                                                                                                                                                                             |
| acrylic fractures when lifted                 | improperly insulated<br><br>undercuts were blocked out insufficiently or not at all                                                                                                                                  | check insulation<br><br>block out undercuts sufficiently                                                                                                                                                                                         |
| acrylic / prosthesis shows whitish streaks    | <ul style="list-style-type: none"> <li>– improperly mixed</li> <li>– powder content too high</li> </ul>                                                                                                              | <ul style="list-style-type: none"> <li>– mix acrylic dough thoroughly with a spatula</li> <li>– observe mixing ratio</li> </ul>                                                                                                                  |
| splint fractures                              | design too thin / unsuitable                                                                                                                                                                                         | ensure a sufficient layer thickness (min. 2 - 3 mm)                                                                                                                                                                                              |

**Guarantee:**

Our processing instructions, regardless of whether they are issued verbally, in writing or in the form of a practical demonstration, are based upon our own experience and can therefore be regarded solely as guidelines. Our products are subject to continuous development. We therefore reserve the right to make possible changes.

**General notes:**

For deviations from 21 °C / 69 °F room temperature:

Higher temperatures reduce the above times, lower temperatures increase them.

To further reduce the risk of intolerance reactions, we recommend storing the orthodontic appliance in lukewarm water for at least 12 hours before insertion.

**Storage:**

Storage temperature at 10 - 25 °C / 50 - 77 °F. Close container carefully.

**Shelf-life:**

The maximum shelf-life is printed on the label of each packaging. Do not use after the expiration date.

**Side-effects:**

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

**Contra-indications / interactions:**

If a patient has known hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth when using the product.

**Note:**

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

**Disposal:**

Powder can be disposed of as domestic waste. Liquid must be disposed of as hazardous waste.

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Schütz Dental GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany