
**IMPLA**

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## Implantat-System

Artikel-Nr.: 608025

**SCHÜTZ**  **DENTAL**  
Micerium Group

**CE 0297**

Schütz Dental GmbH • Dieselstr. 5-6 • 61191 Rosbach / Germany  
Tel.: +49 (0) 6003 814-0 • Fax: +49 (0) 6003 814-906  
[www.schuetz-dental.de](http://www.schuetz-dental.de) • [info@schuetz-dental.de](mailto:info@schuetz-dental.de)

Mandler 04 / 2022 – 1.000

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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# List of contents

Page

IMPLA – Implant Systems Dual Surface and Micro Retention .....	3–6
IMPLA – Interim implants .....	7–8
IMPLA – Implant System Cylindrical .....	9–12
IMPLA – Mini-Implants conetop/balltop .....	13–14
General notes on implantation .....	15–16
IMPLA Abutments .....	17
Locator® Implant Abutment System .....	18–22
HSL Abutment Connectors .....	23–24
Multi-Unit Abutments .....	25–26
IMPLA PreFace Abutments .....	27
IMPLA Aesthura Abutments .....	28
Disinfection, cleaning and sterilization of surgical and prosthetic tools and instruments of the IMPLA system .....	29–32
Application and safety instructions .....	33
Manufacturer’s information on rotating surgical instruments Standardized manual reprocessing of rotating instruments .....	34–35
Information on storage, disinfection, cleaning procedures and sterilization of rotating instruments .....	36

# IMPLA – Implant Systems Dual Surface and Micro Retention

IMPLA implants are made of pure titanium (ISO 5832-2). The IMPLA system comprises surgical, prosthetic and laboratory components and instruments.

## Important notes, warning

IMPLA implants are part of a complete system and may only be used with the original components and instruments of the system. The use of components which are not part of this system voids the guarantee and rules out parts being replaced. The IMPLA system may only be used by dentists, doctors and dental technicians familiar with dental implantology including diagnosis, preoperative planning, surgical techniques and prosthodontics.

Consultation is done in written form or by electronic media as well as with demonstrations during special workshops. All information is based on the state of the art and latest scientific findings known at the time of launching the product. However, we urge the users of our products to check the suitability of the product in each case. This is especially important for techniques not expressly recommended by the manufacturer. In case of doubt, we encourage our customers to confer with Schütz Dental Group/Section IMPLA.

## Reuse and sterilization

All abutments and components are for single use only and must not be reprocessed. Any reuse involves the risk of cross-contamination! All instruments, abutments and components are packaged non-sterile and must be sterilized (EN 13060) prior to placement (see page 29–36). Using unsterile components increases the infection risk for the patient.

## Indications

IMPLA implants Dual Surface and Micro Retention are indicated for:

- restoring individual tooth gaps
- bridgework spanning interdental and free-end gaps according to general and special principles of prosthetic treatment guidelines.
- restoration of edentulous mandibles and maxillae

Further indications:

Implants with a diameter of 3.3 mm

These implants are an alternative in case of limited jaw ridge width. They are, because of their lower mechanical strength when compared with the standard diameter (4.2 mm), intended for the following indications:

- a) restoring individual anterior teeth of the mandible
- b) restoring lateral upper incisors
- c) four interforaminal implants to restore the edentulous mandible with a bar construction
- d) **Do not** use with ball anchors, magnets or telescopes

Implants with a diameter of 4.2 mm

The IMPLA standard implants with a diameter of 4.2 mm are used for all indications of oral enossal implants for mandibles and maxillae. A prerequisite is an adequate minimal amount of bone base with a width of about 5.5 mm.

Implants with a diameter of 5.3 mm

These implants are indicated for the use in the posterior tooth areas of mandible and maxilla (especially for the restoration of molars). The jaw ridge width should not be less than 6.5 mm.

The implants are primarily indicated for delayed implantations.

## **Contraindications**

When selecting patients, the general contraindications for dental surgery must be taken into account. These include: Vascular diseases, uncontrolled diabetes mellitus and other metabolic or systemic diseases affecting bone regeneration, blood coagulation disorder, anticoagulative therapy, bone metabolism diseases.

Temporary contraindications: Acute and chronic inflammations, inadequate soft tissue present, chemotherapy or radiotherapy.

Oral contraindications:

- patient not really prepared to undergo total oral rehabilitation
- inadequate oral hygiene
- untreated parafunctions, e. g. bruxism
- inadequate amount of bone available
- inadequate interocclusal distance
- poor occlusion and/or articulation

## **Diagnostics / Patient information:**

Diagnostics consist of a medical history, clinical examination, X-ray examination using OPG and, if necessary, CT. A medical check up is recommended. It is also advisable to perform biomechanical treatment planning to avoid overloading the bone, the implant or the prosthetic construction, taking into account the individual masticatory forces and occlusion. Implantations are elective procedures that require extensive patient education. Treatment alternatives and consequences of the implantations must be pointed out. In particular, the importance of oral hygiene and the costs of follow-up treatment must be pointed out. With regard to the type and scope of the information, reference is made to the relevant case law.

## **Side effects:**

- Horizontal and vertical bone loss may come about even with a medically appropriate surgical procedure. The extent and the nature of bone loss cannot be estimated in advance.
  - Reversible or irreversible impairments or interactions are possible if the implants are placed near special structures (e.g. adjacent teeth, maxillary sinus or nerves).
  - In connection with dental implants, chronic pain may also occur in rare cases.
- To minimize any risk of infection, antibiotic prophylaxis may be recommended. In rare cases, allergies to the prosthetic components (immediate implantation) can also prevent successful implantation.

## **Important notes**

When working with our implants and instruments intraorally, always use a protection against aspiration! Primary stability directly after the insertion of the implant is imperative for a successful healing. During the healing phase, IMPLA implants are generally not functionally loaded.

## **Surgical procedure**

To create optimal conditions for the implant to heal, it is essential to treat the hard and soft tissues gently. It is advisable to use the countersink and bone thread cutter in order to ensure that the implant can be placed atraumatically applying as little torque as possible (maximum 30 Ncm). As thermal trauma may impede or prevent the implant healing, heat generation should be kept to a minimum as follows:

- recommended drill speeds:  
2.0 mm Ø: approx. 800 rpm, 2.5 mm Ø: approx. 700 rpm, 3.15 mm Ø: approx. 500 rpm, Ø 4.2 mm, Ø 4.8 mm, Ø 5.3 mm: approx. 400 rpm
- recommended speeds: pilot drill Ø 1.2 mm and Ø 2.0 mm - approx. 800 rpm, enlargement drill Ø 3.15 mm - approx. 500 rpm, cavity drills of all Ø - approx. 400 rpm

- intermittent drilling and cutting using sharp drills and cutters (use 10 times only)
- cool the instrument externally with sterile physiological saline (NaCl) or Ringer's solution.

### **Important**

To determine the optimum implant length and prevent nerv injuries, check in advance that adequate space is available (X-ray).

### **Placing the implant after preparing the implant bone supporting area**

#### **Removing the insertion aid**

In order to remove the insertion aid (for implant 3.3 mm brown; 4.2 mm red; 5.3 mm blue), pull off the plastic implant carrier and loosen the vertical screw (blue) with the hexagonal screwdriver SW 1.2 mm and dispose of the vertical screw. Secure the implant against rotation with help of a counter wrench in cases of low primary stability and concomitant augmentation in the according alveole. Afterwards, the insertion aid can be lifted from the implant easily.

#### **Closing the soft tissue**

The cover / healing screw must be screwed onto the implant before the wound is closed with atraumatic suturing material.

#### **Recovering the implant**

Once the individual healing period (minimum of three months for the mandible and six months for the maxilla) has elapsed, the implant should be recovered and the cover / healing screw removed. In case of augmentation techniques, the healing period must be prolonged accordingly. The interior of the implant is then cleaned and a gingiva former suitable to the height of the mucous membrane is screwed in. Make sure that the gingiva former fits properly! The restoration should not be fitted to the implant until the soft tissue has healed fully and shows no irritations.

### **Impression taking**

#### **Open impression technique:**

After the gingiva former has been removed, a forming aid with the appropriate diameter is fixed to the implant with help of a fixation screw. When taking the impression, the impression material is first syringed around the impression post. Once the impression material has cured, loosen the screw, remove the impression from the mouth together with the impression post and close the implant again with the gingiva former. Carefully screw an implant analog into the impression post and cast the model – preferably including a gingival mask. Once the cast has set, remove the impression post, respectively transfer coping, and replace it with the required abutment.

#### **Closed technique:**

Impression posts of the appropriate diameter are screwed onto the implant using a laboratory screw. Impressions are taken using a standard impression tray with the aid of transfer caps which are clicked onto the impression posts. The transfer caps remain in the impression. To produce a plaster model, they are fitted with a lab analog after the impression posts have been reinserted.

### **Laboratory components**

When selecting the abutments, the angulation of the implant and the thickness of the gingiva must be taken into account. Those areas of the abutment which are in contact with the implant may not be blasted with abrasives or treated in any other manner.

### **Documentation**

Each IMPLA implant can be traced via its lot number which is noted on a self-adhesive label for internal user documentation. Possible claims can only be asserted by providing this data.

### **Packaging and sterility**

IMPLA implants are packed sterilely in a special clear acrylic container and are color-coded. The chosen implant is removed from the sterile container with the insertion aid only directly before insertion and is inserted into the prepared bone cavity without touching the implant. The covering and healing screw is attached to the back of the implant carrier. It can easily be removed from there with the hex screwdriver SW 1.2 mm. Important reference data on implant length, implant diameter, article number, batch, CE certification, etc. are printed on the outer packaging made of cardboard.

**General information**

- All implants are for single use only.
- Each is packed sterilely in a closed container, please note the expiration date of sterility.
- The implant is sterile (radio sterilisation) as long as the packing remains unopened and undamaged. A subsequent sterilization by the user is not permitted by the manufacturer. A chemical sterilization of the implant (chemiclave) may not be carried out under any circumstances, as this could damage the implant surface. In case of damaged packaging, the implant is no longer sterile and must not be used anymore.
- Do not touch the implant.
- Store the implant in a dry place at room temperature. Avoid direct light incidence.
- System errors due to mixing up instruments and implants should be avoided. Please always check the color codings and labelings.
- Thermal damages of the jaw bone by rotating instruments (such as pilot drills, thread cutters or enlargement drills) must be avoided.
- It is advisable that the implantation is performed only after the patient information consultation has been completed, after the patient has signed the consent form and after the patient has undergone a preliminary medical examination for implantation.

**Torque of screw connections**

15 Ncm	horizontal screw bar screw cone adapter Multi Unit secondary prosthetic screw
20 Ncm	vertical screw cone adapter Multi Unit abutment 0° Multi Unit vertical screw Aesthura fixing screw
30 Ncm	Locator abutment ball connector

# IMPLA – Interim implants

## Note:

Before performing implantology treatment, it is essential to be familiar with pertinent professional literature. It is mandatory to have experience in surgery and in the assessment of the normal and pathological course of treatment.

**The „General notes on implantation“ (see page 15–16) must be considered before treatment.**

## Description:

The one-part IMPLA interim implant is made of grade 5 titanium in accordance with ISO 5832-3.

**Warning: If the implant breaks, it is possible that fragments remaining cannot be removed without damage to the surrounding bone tissue!**

Implant threads: The IMPLA interim implant has primary self-cutting threads, allowing primarily stable insertion into the jaw bone while preserving bone material. The screw design and resulting compression stability along with the pilot drill ensure sufficient primary stability even in spongy bone tissue.

Implant surface: Due to its intended use as an interim implant (see indications), the implant surface was machine polished. This is to achieve the possibility of removing the implant after the desired retention time.

**Important: Due to a possibly occurring higher insertion torque, the IMPLA positioning key (IPK) must be used with great caution to avoid the danger of fracturing or deforming the implant.**

## Indications:

The IMPLA interim implant is a temporary support for drill templates or for a temporary bridge/denture during the healing phase of the permanent implants for patients who desire a stable, esthetically pleasing temporary prosthesis. The soft tissue around conventional implants is protected. Any existing augmented area is protected from strain. The mucosa must be healthy with good circulation.

## Contraindications:

### Absolute:

- A sufficient number of temporary immediate implants cannot be placed in the implant area to support the temporary denture.
- Patients whose systemic condition does not allow implantation.

### Relative:

- Patients who are not capable of follow-up care for psychological or physical reasons and who are expected to cause complications
- Patients undergoing treatment for a malignant tumor
- Patients with uncontrolled diabetes
- Bruxism
- Heavy smokers
- Insufficiently available bone tissue and/or immediate vicinity to endangered structures (nerve, sinuses)

## Diagnostics / patient information:

Diagnostics consist of the patient history, clinical examination, X-ray examination with OPG and CT if necessary. A medical check-up is recommended.

It is also advisable to perform biomechanical treatment planning taking into account the individual masticatory forces and occlusion in order to avoid overloading the bone, the implant or the prosthetic construction.

Implantations are elective procedures for which the patient must be given comprehensive information. Treatment alternatives and consequences of the implantation must be indicated. In particular, the importance of oral hygiene and costs of further treatment must be discussed. The relevant laws indicate the kind and extent of information required.

## Side effects:

- Even in a professional surgical procedure, potential horizontal and vertical bone loss must be taken into account as with all dental implants; the extent and kind of bone loss cannot be estimated in advance.
- Reversible or irreversible impairment or interaction is possible when implants are inserted in the vicinity of certain structures (e.g. adjacent teeth, maxillary sinus, or nerves).
- In rare cases, chronic pain can occur in association with dental implants.

To minimize the potential risk of infection, antibiotic prophylaxis can be recommended. In rare cases, allergies to the components of the implant may prevent successful implantation.

**Procedure:**

The entire course of treatment has three stages:

1. Pilot drilling;
2. Initial insertion of the implant;
3. Further insertion using the IMPLA positioning key (IPK)

For stage 1) A pilot drilling is made through the mucosa at the optimum site, guided and controlled by the axis and length of the implant. A minimum distance of about 2 mm from the final implants must be maintained. The pilot drill must be selected depending on the bone quality (D1–D2 => 2.0 mm bit in the lower jaw; D3–D4 and upper jaw => 2.0 mm bit with reduced drilling depth). Cooling with sterile saline solution is mandatory. To prevent potential thermal damage, for bone quality D1–D2, the holes should be pre-drilled with the pre-drill and the cavity widened with the 2.0 mm pilot drill if necessary.

**Important: Drill alternately at a speed of 500–600 rpm.**

The entire implant length must be pre-drilled. Less in the upper jaw, if the bone quality is poor (30–50%). For stage 2) The packaging must not be opened until immediately before insertion of the implant. Hold the implant at the thread using the packaging as protection and open the side facing the implant head. Place the positioning key (IPK) firmly on the implant head and make sure it is gripping firmly. Do not pull the implant out of the packaging until you are near the implantation site and insert it there immediately using the IMPLA positioning key (IPK) to prevent contamination.

For stage 3) Lower the implant with the IMPLA positioning key (IPK) until the upper edge of the thread is flush with the gingiva. If necessary, the torque wrench must be used to reach the final position.

**Important: Due to great torque, the IMPLA positioning key (IPK) is to be used with utmost caution to avoid fracturing the implant.**

**Removal of the temporary implant:**

The temporary implant must be removed prior to or at the time the permanent implants are inserted. Detach the temporary implant and remove all remaining cement. Put the IMPLA positioning key (IPK) and the ratchet on the implant and inspect for correct fit. Paying attention to torque, turn the implant counterclockwise to unscrew. Remove any granulation tissue with a curette, as conditions similar to tooth extraction can be expected. **The temporary implant is usually kept for 3 to 6 months. Implant healing must be monitored (e. g. X-ray) regularly. Infected or loose implants must be removed promptly.**

**Notes:**

- After implantation, the type of implant used and LOT number must be recorded in the documentation (sticker).
- Implants must be stored in a closed container in a dry place. The seal may not be broken until immediately before insertion. Before insertion, any contact with foreign substances must be avoided. The endosseous part of the implant must not be touched. In the event the implant, superstructure, or an accessory is swallowed, the location of the object must be determined and the appropriate medical steps are to be taken.
- Healing at the implant site must be monitored regularly, e. g. by X-ray. Implants that are infected or not osseointegrated must be removed promptly. The dentist must determine the time for removal.

**Scope of delivery:**

**The IMPLA interim implant is delivered in non-sterile condition.**

**The user must sterilize the implants in an autoclave at 134°C before use. We accept no responsibility for the correct handling and functioning of the sterilizer.**

**The implant is sealed in a sterilizable foil pouch equipped with a label for transferring the implant data, such as length and the lot number, to the patient's file. This foil pouch is the actual packaging. The IMPLA interim implant is intended to be used one time only!**

**If the packaging is damaged or not airtight, there is a risk of contamination with bacteria. In this case, the implant must not be used.**

# IMPLA – Implant System Cylindrical

IMPLA implants are made of medical grade 4 pure titanium (ISO 5832-2). The IMPLA system comprises surgical, prosthetic and laboratory components and instruments.

## **Important notes, warning**

IMPLA implants are part of a complete system and may only be used with the original components and instruments of the system. The use of components which are not part of this system voids the guarantee and rules out parts being replaced. The IMPLA system may only be used by dentists, doctors and dental technicians familiar with dental implantology including diagnosis, preoperative planning, surgical techniques and prosthodontics.

Consultation is done in written form or by electronic media, with demonstrations during special workshops and by the implantology field staff. All information is based on the state of the art and latest scientific findings known at the time of launching the product. However, regardless of this, the user of our products is required to carry out a personal suitability test in the specific case and based on the indication. This is especially important for techniques not expressly recommended by the manufacturer. In case of doubt, we encourage our customers to confer with Schütz Dental GmbH/Section IMPLA.

## **Reuse and sterilization**

All abutments and components are for single use only and must not be reprocessed. Any reuse involves the risk of cross-contamination! All instruments, abutments and components are packaged non-sterile and must be sterilized (autoclaved EN 13060) prior to placement (see page 29–36). Using unsterile components increases the infection risk for the patient.

## **Indications**

IMPLA implants Cylindrical Implants are indicated for:

- restoring individual tooth gaps
- bridgework spanning interdental gaps and free-end gaps, in accordance with conventional and special principles of prosthetic treatment
- restoration of edentulous mandibles and maxillae

### **Further indications:**

Implants with a diameter of 3.6 mm

These implants are an alternative in case of limited jaw ridge width. They are, because of their lower mechanical strength when compared with the standard diameter (4.5 mm), intended for the following indications:

- a) restoring individual anterior teeth of the mandible
- b) restoring individual lateral anterior teeth of the maxilla
- c) four interforaminal implants to restore the edentulous mandible with a bar construction

Implants with a diameter of 4.0 mm and 4.5 mm

The IMPLA standard implants with a diameter of 4.0 mm and 4.5 mm are used for all indications of oral enossal implants for mandibles and maxillae. A prerequisite is an adequate minimal amount of bone base with a width of about 5.3 mm and 5.8 mm.

Implants with a diameter of 5.5 mm

These implants are indicated for the use in the posterior tooth areas of mandible and maxilla (especially for the restoration of molars). The jaw ridge width should not be less than 6.8 mm.

The implants are primarily indicated for delayed implantations.

The implants with a length of 6.5 mm are generally not intended for single-tooth replacement or for accommodating a stand-alone prosthetic restoration.

### **Contraindications**

When selecting patients, the general contraindications for dental surgery must be taken into account. These include: Vascular diseases, uncontrolled diabetes mellitus and other metabolic or systemic diseases affecting bone regeneration, blood coagulation disorder, anticoagulative therapy, bone metabolism diseases.

Temporary contraindications: Acute and chronic inflammations, inadequate soft tissue present, chemotherapy or radiotherapy.

Oral contraindications:

- patient not really prepared to undergo total oral rehabilitation
- inadequate oral hygiene
- untreated parafunctions, e. g. bruxism
- inadequate amount of bone available
- inadequate interocclusal distance
- inadequate occlusion and/or articulation

### **Diagnostics / Education:**

Diagnostics consist of medical history, clinical examination, X-ray examination using OPG and CT, if necessary. It is recommended to have a medical check-up. It is also advisable to perform biomechanical treatment planning to avoid overloading the bone, implant or prosthetic construction, taking into account the individual occlusal forces and occlusion. Implants are elective procedures, about which extensive information must be provided to the patient. Treatment alternatives and consequences of the implantation must be pointed out. In particular, the importance of oral hygiene and the costs of follow-up treatment must be pointed out. With regard to the type and scope of the information, reference is made to the relevant case law.

### **Side effects:**

- Horizontal and vertical bone resorption may occur even with a medically professional surgical procedure, and the extent and nature of bone resorption cannot be estimated in advance.
  - Reversible or irreversible impairments or interactions are possible when implants are placed near special structures (e. g. adjacent teeth, maxillary sinus or nerves).
  - Chronic pain may also occur in rare cases in connection with dental implants.
- To minimize any risk of infection, antibiotic prophylaxis may be recommended. In rare cases, allergies to the prosthetic components (immediate implantation) can also prevent successful implantation.

### **Important notes**

When using implants and instruments intraorally, it is generally important to ensure that they are secured against aspiration!

Primary stability directly after the insertion of the implant is absolutely necessary for a successful healing. During the healing phase, IMPLA implants are generally never exposed to a masticatory load.

### **Surgical procedure**

To create optimum conditions for the implant to heal, it is essential to treat the hard and soft tissues gently. It is advisable to use the countersink in order to ensure that the implant can be placed atraumatically applying as little torque as possible (maximum 30 Ncm). The so-called „Dense Bone Drills“ have been designed especially for the cortical bone preparation. They are meant to reduce the insertion torque for longer implants. As thermal trauma may impede or prevent the implant healing, heat generation should be kept to a minimum as follows:

- recommended drill speeds:  
2.0 mm Ø: approx. 800 rpm, 2.75 mm Ø: approx. 700 rpm, 3.0 mm Ø, 3.25 mm Ø and 3.4 mm Ø: approx. 500 rpm, 3.7 mm Ø, 3.9 mm Ø, 4.25 mm Ø, 4.9 mm Ø and 5.25 mm Ø: approx. 400 rpm
- intermittent drilling and cutting using sharp drills and cutters (use 10 times only)
- cool the instrument externally with sterile physiological saline (NaCl) or Ringer's solution.

**Important**

To determine the optimum implant length and prevent nerve injuries, check in advance that adequate bone is available (X-ray).

**Placing the implant after preparing the implant bone area****Removing the insertion post**

In order to remove the insertion post (for implant 3.6 mm brown; 4.0 mm green; 4.5 mm red; 5.5 mm blue), pull off the plastic implant carrier and loosen the vertical screw (blue) with the hexagonal screwdriver 1.2 mm. Dispose of the vertical screw. In cases with concomitant augmentation or where the primary stability is low, secure the implant against rotation with help of a counter wrench. Afterwards, the insertion aid can be lifted from the implant easily.

**Closing the soft tissue**

The cover / healing screw must be screwed onto the implant before the wound is closed with atraumatic suturing material.

**Exposure of the implant**

Once the individual healing period (minimum of three months for the mandible and six months for the maxilla) has elapsed, the implant should be exposed and the cover / healing screw removed. In case of augmentation techniques, the healing period must be prolonged accordingly. The interior of the implant is then cleaned and the correct height of gingiva former for the oral mucosa is screwed in. Make sure that the gingiva former fits properly! The restoration should not be fitted to the implant until the soft tissue has healed fully and shows no irritations.

**Impression taking****Open technique**

After removing the gingiva former, the impression aid of the appropriate diameter is screwed into the implant using a fixing screw. If an impression is subsequently taken with a customized impression tray, the impression aid is first coated with impression material. After the impression material has hardened and the fixing screw has been loosened, the impression together with the impression aid is removed and the implant is closed again with the gingiva former. A lab analog is carefully screwed to the impression post and then a model - preferably with a mucosa mask - is created. After the plaster has hardened, the impression aid or transfer abutment is removed and replaced by an appropriate prosthetic abutment.

**Closed technique**

Impression posts of the appropriate diameter are screwed onto the implant using a laboratory screw. Impressions are taken using a standard impression tray with the aid of transfer caps which are clicked onto the impression posts. The transfer caps remain in the impression and, after reinsertion of the impression posts, are fitted with lab analogs to produce a plaster model.

**Laboratory components**

When selecting the abutments, the angulation of the implant and the thickness of the gingiva must be taken into account. Those areas of the abutment which are in contact with the implant may not be blasted with abrasives or treated in any other manner.

**Documentation**

Each IMPLA implant is traceable by its lot number, which is noted on a self-adhesive label for internal user documentation. Any claims can only be asserted if this data is provided.

**Packaging and sterility**

IMPLA implants are packed sterile in a special clear container.

The chosen implant is removed from the sterile container with the insertion aid only directly before insertion and is inserted into the prepared bone cavity without touching the implant.

A special attachment on the bottom of the implant container holds the covering and healing screw which can be easily removed with the hex screwdriver SW 1.2 mm.

The outer cardboard box shows the user important data such as implant length, implant diameter, article no., lot and CE certification.

**General information**

- All implants are for single use only.
- Each implant is packed sterilely in a closed container, please note the expiration date of sterility.
- The implant is sterile (radiosterilization) as long as the packing remains unopened and undamaged. A subsequent sterilization by the user is not permitted by the manufacturer. A chemical sterilization of the implant (chemiclave) must not be carried out under any circumstances, as this could damage the implant surface. In case of damaged packaging, the implant is no longer sterile and may not be used anymore.
- Do not touch the implant.
- Store the implant in a dry place at room temperature. Avoid direct light incidence.
- System errors due to mixing up instruments and implants are to be avoided. Pay attention to the color codings and labelings.
- Thermal damages of the jaw bone by rotating instruments (such as pilot drills and extension drills) must be avoided.
- It is strongly recommended to perform implant surgery only after the patient consultation, after having the patient sign a consent form and after a medical check up. A surgery protocol is also recommended. Follow-up care with patient recalls to monitor the healing process and to support the patient's behavior after implantation is self-evident.

**Tightening torque (torque) of screw connections!**

10 Ncm horizontal screw connections

20 Ncm vertical screw connections

30 Ncm insertion of the implant, Locator Abutment screw connections

# IMPLA – Mini-Implants conetop / balltop

## Note:

Before performing implantology treatment, it is essential to be familiar with pertinent professional literature. It is absolutely necessary to have experience in surgery and in the assessment of the normal and pathological course of treatment.

**The „General notes on implantation“ (see page 15–16) must be considered before treatment.**

## Description:

The one-piece IMPLA Mini-Implant is made of pure grade 4 titanium in accordance with ISO 5832-2.

**Warning: If the implant breaks, it is possible that fragments remaining cannot be removed without damage to the surrounding bone tissue!**

**Implant threads:** The IMPLA Mini-Implant has primary self-cutting threads, allowing primary stable insertion into the jaw bone while preserving bone material. The screw design and resulting compression stability together with the pilot drill and the extension drill ensure sufficient primary stability even in spongy bone tissue.

**Implant surface:** The implant was designed with a structured surface in accordance with the implant's application as a permanent implant (see indications).

## Indications:

The IMPLA Mini Implant was developed for the following indications:

Mini implant „conetop“:

- Fixing full dentures
- Bar-supported dentures (on a minimum of 4 implants)
- Support implant for fixing bridges

Mini implant „balltop“:

- Fixing full dentures with help of an O-ring attachment

Contraindications:

### Absolute:

- Patients whose systemic condition does not allow implantation.

### Relative:

- Patients who are not capable of follow-up care for psychological or physical reasons and whose compliance is not anticipated
- Patients undergoing treatment for a malignant tumor
- Patients with uncontrolled diabetes
- Bruxism
- Heavy smokers
- Insufficiently available bone tissue and/or immediate vicinity to endangered structures (nerve, sinuses)

## Diagnostics / Patient information:

Diagnostics consist of a medical history, clinical examination, X-ray examination using OPG and, if necessary, CT. A medical check up is recommended. It is also advisable to perform biomechanical treatment planning to avoid overloading the bone, the implant or the prosthetic construction, taking into account the individual masticatory forces and occlusion. Implantations are elective procedures that require extensive patient education. Treatment alternatives and consequences of the implantations must be pointed out. In particular, the importance of oral hygiene and the costs of follow-up treatment must be pointed out. With regard to the type and scope of the information, reference is made to the relevant case law.

## Side effects:

- Even in a professional surgical procedure, potential horizontal and vertical bone loss must be taken into account as with all dental implants; the extent and kind of bone loss cannot be estimated in advance.
- Reversible or irreversible impairment or interaction is possible when implants are inserted in the vicinity of certain structures (e. g. adjacent teeth, maxillary sinus, or nerves).
- In rare cases, chronic pain can occur in association with dental implants.

To minimize the potential risk of infection, antibiotic prophylaxis can be recommended. In rare cases, allergies to the prosthetic components (immediate implantation) may prevent successful implantation.

**Procedure:**

The entire course of treatment has five stages:

1. Pilot drilling, 1.8 mm
2. Pilot drilling, 2.0 mm
3. Extension drilling, 2.5 mm
4. Initial implant insertion
5. Further insertion all the way to the implant's final position with help of the IMPLA screwdriver.

Point 1) Pilot drilling is performed at the optimal site (optionally also through the mucosa) with specification and control of the implant axis and length.

Points 2) and 3) The cavity is then finally prepared with the extension drill (diameter 2.5 mm, only for Mini conetop) according to the implant diameter 2.1 mm, 2.5 mm or 3 mm to be inserted. The distance to the definitive implants should not fall short of approx. 2 mm.

**The treatment must be adapted to the bone situation:**

**D1–D2 => 1.8 mm => 2.0 mm => 2.5 mm drill in the lower jaw**

**D3–D4 and upper jaw => 1.8 mm => 2.0 mm => 2.5 mm drill with decreased drilling depth**

Cooling with sterile saline solution is advisable.

**Important: Drill alternately at a speed of 500–600 rpm.**

The entire implant length should be predrilled with the 1.8 mm respectively 2.0 mm pilot drill.

Point 3) The preparation depth with the extension drill 2.5 mm should be adapted to the bone quality (see point 1).

Point 4) The packaging must not be opened until immediately before insertion of the implant. Remove the glass tube with the implant inside and remove the implant with the insertion aid. With help of the insertion aid, insert the implant into the cavity. Once an increasing torque is reached, the plastic insertion aid is removed.

Point 5) Use the IMPLA screwdriver to screw in the implant until the top edge of the structured implant area is level with the bone. If necessary, use the torque wrench to reach the final position.

**Important: Because of the thread configuration of the implant, a torque of min. 25 to max. 50 Ncm should be reached, depending on the indication.**

**Notes:**

- After implantation, the type of implant used and the LOT number must be recorded in the documentation (sticker).
- Implants must be stored in a closed container in a dry place. Before insertion, any contact with foreign substances must be avoided. The endosseous part of the implant must not be touched. In the event the implant, superstructure, or an accessory is swallowed, the location of the object must be determined and the appropriate medical steps are to be taken.
- Healing at the implant site must be monitored regularly, e. g. by x-ray. Implants that are infected or not osseointegrated must be removed promptly. The dentist must determine the time for removal.

**Scope of delivery:**

The IMPLA Mini-Implant is delivered sterile (radio sterilisation).

The implant is sealed in a film bag bearing a label for entering the implant data such as length and LOT number that can be put in the patient file. This film wrapping is the actual packaging.

The IMPLA Mini-Implant is intended to be used one time only!

If the packaging is damaged or not airtight, there is a risk of contamination with bacteria. The implant must not be used in this case.

# General notes on implantation

## Attention

All rotating instruments are subject to wear. Therefore, we cannot warrant the instruments' (pilot drills, extension drills, countersinks, cortical drills and thread cutters) functionality for more than 10 times of use, depending on the bone density. If the instruments are used more often, the manufacturer cannot guarantee a proper cutting efficiency and minimal temperature development any longer.

- Insufficient primary stability can be caused by poor preparation of the implant site or by insufficient bone quality. In this case, the loose implant should be removed and be replaced by an implant with a larger diameter. Implants should always be screwed into the bone all the way, even with their thread-free neck, when a covered healing process is considered (see above). Above all, a primary stability of the implant must be achieved.
- Implants should only be implanted when all prerequisites for a complication-free treatment and a longterm success are given. If there is any concern regarding the success of the implantation and/or the stability of the desired implant-based denture, a different solution without implantation should be favored.
- All rotating instruments can, in principle, be used several times. The user physician/expert alone decides under his/her own responsibility on the possibility of using the instrument again and the frequency of use on the basis of the respective case of use and the possible wear of the instruments. In case of doubt, the instruments should always be sorted out and replaced at an early stage.

## Care of instruments and parts

- Used instruments/parts must always be disinfected before cleaning.
- Damaged instruments/parts must be disinfected, cleaned, and disposed of.

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**Please consider our information on each material on page 29 – 36.**

## Disinfection

- Instruments/parts must be submerged in an adequate disinfection solution.
- Disinfection may be carried out manually or automatically. Please observe the manufacturer's instructions concerning dosing/concentration, residence time, temperature.
- In order to avoid damages, the instruments/parts may not have any contact during disinfection.
- It is highly recommended to wear protective gloves when working with contaminated instruments/parts (components).

## Cleaning

- Only use adequate cleaning agents.
- Always observe the manufacturer's instructions concerning concentration of the agent, temperature, etc.
- It is highly recommended to wear protective gloves (DIN EN 455) when cleaning contaminated instruments/parts manually.
- In case of automatic cleaning and disinfection, the instruments/parts may not have any contact in order to avoid damages.
- Only instruments/parts made of the same material should be cleaned together.

- In case of heavy contamination of the instruments/parts, ultrasound cleaning techniques are recommended.
- After cleaning, rinse the instruments/parts thoroughly with water and dry them immediately.
- Damaged instruments/parts must be removed and disposed of separately.

**Sterilization**

- All surgical instruments and parts must be sterilized after cleaning.

**Please consider our information on each material on page 29–36.**

**Storage**

- The sterilized instruments/parts must be stored in a dry and dust-free atmosphere at room temperature.
- After the expiration of the maximum storage period, the instruments/parts must be sterilized again.

Dose 25 kGy/min

# IMPLA Abutments\*

\*The term „abutments“ is used here to represent ball connector, adhesive base, cone connector, SAE cone adapter, IMPLA Direct base.

## **Important notes / warnings:**

IMPLA Abutments are a part of the IMPLA system and must only be used with the proper original IMPLA components. For the insertion of the abutments, only use the IMPLA instruments and tools. Otherwise, a damage of the original components cannot be ruled out.

## **Delivery:**

All IMPLA Abutments are supplied in unsterile packaging. Therefore, the components must be sterilized accordingly before use.

## **Material:**

IMPLA Abutments are manufactured from titanium grade 5 (Ti6Al4V ASTM 136).

IMPLA Scan abutments are made from PEEK (polyether ether keton).

## **Side-effects:**

Allergies and sensitivities towards the alloy cannot be ruled out completely and may occur in very rare cases.

## **Interaction:**

In case of occlusal contact, the presence of different alloy types in the patient's mouth may lead to galvanic reactions.

## **Safety advice:**

Metal dust and ceramic dust are hazardous to the health. When working on these materials, always wear a protective mask and protective goggles. Use an extraction system with a fine particle air filter commonly used in dental practices.

## **Indications:**

IMPLA ball abutments are used to anchor removable restorations to IMPLA implants.

The IMPLA adhesive base is used to fix CAD constructed and CAM produced restorations on IMPLA implants. This requires a suitable design software as well as material blanks and a milling machine to produce the prosthetic structures, which can then be bonded to the adhesive base.

Cone connectors are intended for the production of cementable individual crowns and bridges.

SAE cone adapters are used for the production of bar abutments and bridges which are screw-retained in a tension-free method.

The IMPLA Direct base is intended for the construction of bar abutments or screw-retained bridges using the sprue technique, casting technique or laser technique.

Scan abutments are only intended for temporary use in the patient's mouth for the time of optical scanning using the intraoral scanner.

## **Contra-indications:**

When selecting the patients, it is extremely important to ensure that they are able to maintain the necessary oral hygiene. In addition, the general rules and limitations in surgery and prosthetics for the selection and treatment of patients apply.

## **Procedure:**

If the abutment does not bear a screw thread, carry out the final fixation on the intended implant with the vertical screw standard, art. no. 636649, unless otherwise specified for the respective abutment type in separate instructions.

# LOCATOR® IMPLANT ABUTMENT SYSTEM

**IMPORTANT:** This document contains the most current instructions for use. Please read and retain.

**CLASSIFICATION:** Universal hinge, resilient attachment for endosseous implants.

## ATTACHMENT COMPONENT ORDER NUMBERS

Product numbers of Locator Implant Abutments vary with implant type, implant diameter, and tissue cuff height. It is necessary to have this information available to place your order.

## PARTS IDENTIFICATION

Implant Abutment  
 3.3/3.0 mm **No. 636060**  
 3.3/5.0 mm **No. 636061**  
 4.2/3.0 mm **No. 636062**  
 4.2/5.0 mm **No. 636063**  
 5.3/3.0 mm **No. 636064**  
 5.3/5.0 mm **No. 636065**  
 Implant Abutment Titanium alloy with TiN coating

Locator matrix set  
 3 pcs.  
**No. 636070**  
 Titanium Cap with black Low Density Polyethylene Male, Block Out Spacer, and Nylon Repl. Males (clear, pink, blue)

Locator plug connector clear, 4 pcs.  
**No. 636071**  
 Nylon (clear)  
 pull-off force: 2.270 g

Locator plug connector, pink slight retention, 4 pcs.  
**No. 636072**  
 Nylon (pink)  
 pull-off force: 1.360 g

Locator plug connector, blue extra-light retention, 4 pcs.  
**No. 636076**  
 Nylon (blue)  
 pull-off force: 680 g

Locator plug connector, green extended application range, 4 pcs.  
**No. 636073**  
 Nylon (green)  
 pull-off force: 1.820 g


Locator plug connector, red extended application range, slight retention, 4 pcs.  
**No. 636074**  
 Nylon (red)  
 pull-off force: 450 g

Locator impression post  
**No. 636067**  
 Aluminum housing with black Low Density Polyethylene Male

Locator lab implant (4 mm)  
**No. 636068**  
 Aluminium

Locator-Instrument, threepart (Male Removal Tool, Male Seating Tool and Abutment Driver)  
**No. 636066**

Locator Adapter with WS connection  
**No. 636075**

## INDICATIONS

The LOCATOR® Implant Attachment System is designed for use with overdentures or partial dentures retained in whole or in part, by endosseous implants in the mandible or maxilla.

## CONTRAINDICATIONS

Not appropriate where a totally rigid connection is required. Use on a single implant with divergence of greater than 20 degrees is not recommended.

## CAUTION

Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed dentist.

## SINGLE-USE DEVICES

**Locator Males:** The inadvertent re-use of Locator nylon males could cause loss of retention for the overdenture due to wear from previous use or damage during removal with the Locator Core Tool.

**Locator Abutments:** The inadvertent re-use of Locator abutments could contain patient contamination build-up and subsequent wear of the retention bands. This would result in the device to perform with improper fit and function which would result in loss of retention for the prosthesis.

## STERILIZATION

All components and instruments are supplied **NON-STERILE**.

Titanium abutments may be sterilized by autoclave or dry heat sterilization using the following parameters:

1. Autoclave sterilize using 121°C (250°F), (15 – 20 psi at sea level), for 20 minutes minimum.
2. Dry heat sterilize using 170°C (338°F) for 2 hours minimum.

Locator Core-Tools (**the disassembled state only**) may be sterilized by Autoclave or Dry Heat sterilization using the following parameters:

1. Autoclave sterilize using 121°C (250°F), 15 – 20 psi (at sea level), for 40 minutes minimum.
2. Dry heat sterilize using 170°C (338°F) for 2 hours minimum.

## FEATURES

1. **LOWEST VERTICAL HEIGHT:** The total height of the Locator Attachment (abutment plus male) is only 3.17 mm on an externally hexed implant.
2. **LOCATING DESIGN:** Self-locating design allows a patient to easily seat their overdenture without the need for accurate alignment of the attachment components.
3. **RETENTION INSIDE AND OUT:** The patented Dual Retention innovation provides the Locator Attachment with greater retention surface area than ever before available with other attachments. A combination of inside and outside retention ensures the longest lasting performance.
4. **ROTATIONAL PIVOTING ACTION:** The design of the pivoting Locator Male allows a resilient connection for the prosthesis without any resulting loss of retention. The retentive nylon male remains completely in contact with the abutment socket while its titanium denture cap has a full range of rotational movement over the male.
5. **USE WITH NON-PARALLEL IMPLANTS:** The Locator Replacement Males 636072, 636076 can be used to restore an implant with up to 10 degrees of divergence (20 degrees between implants). The Locator Extended Range Replacement Males 636073, 636074 can accommodate a divergent implant between 10 and 20 degrees (40 degrees between implants).

plug connector, red  
**No. 636074**

plug connector, blue  
**No. 636076**

plug connector, pink  
**No. 636072**

plug connector, green  
**No. 636073**

## PLACEMENT OF THE LOCATOR IMPLANT ABUTMENT

1. To select the proper Locator Implant Abutment, determine the type of implant and the diameter being used. Then measure the tissue thickness from the apical rim of the implant body to the crest of the gingiva at the highest side of the implant site. Choose the corresponding abutment tissue cuff height that exactly equals the tissue measurement, or is the next closest higher size available. The exact tissue cuff height of Locator abutment will position the proper 1.5 mm of working attachment above the surrounding gingival level (**which should not be submerged below the tissue**).
2. After the secondary gingival healing period is complete, remove the healing cuff according to instructions provided by the manufacturer of the implant system being used.
3. It is imperative that all bone and soft tissue be removed from the superior aspect of the implant body to guarantee complete seating of the Locator Implant Abutment.
4. A special gold plated Abutment Driver (contained in Locator Core Tool, 636066) is designed to engage the inside diameter of the Locator Abutment and thread it into the implant. (Fig. 1)
5. Final torque tightening of the Locator Abutment to prevent screw loosening is achieved using the 30 N cm Torque Wrench Kit.

**NOTE:** Use the IMPLA Torque wrench with either of these options to achieve the minimum seating force of 30 N cm that will help prevent screw loosening of the locator Implant Abutment.

Fig. 1

Abutment Driver/  
Locator Core Tool  
No. 636066

## LOCATOR MALE PLACEMENT BY THE DENTIST

1. Insertion of the proper Locator Implant Abutment at tissue level must be completed (see Section A-1) before beginning the procedure for placement of the Locator Male.
2. Place a White Block-Out Spacer (636070 Package) over the head of each Locator Abutment. (Fig. 2) The spacer is used to block out the area immediately surrounding the abutment. The space created will allow the full resilient function of the pivoting metal denture cap over the Locator Male.

**NOTE:** If the White Block-Out Spacer does not completely fill the space between the tissue and the metal denture cap, it is necessary to block out any remaining undercuts to prevent the added acrylic resin from locking the denture onto the abutment. This can be accomplished by stacking more Block-Out Spacers.


Fig. 2

White Block-out  
Spacer  
No. 636070

3. Insert a Locator Cap with Black Processing Male (636070 Package) into each Locator Implant Abutment, leaving the White Block-Out spacer beneath it (Fig. 3). The Black Processing Male will maintain the overdenture in the upper limit of its vertical resiliency during the processing procedure.
4. Prepare a recess in the denture to accommodate the protruding Locator Male. There must be no contact between the denture and the titanium cap. If the denture rests on the metal cap, excess pressure on the implant will result.
5. Use a light curing device to light cure bond the Locator Denture Cap Male into the denture (Fig. 4 and 5), or mix a permanent self-curing acrylic and place a small amount in the recess of the denture and around the metal cap of the Processing Cap Male.


Fig. 3

Processing Cap  
Nr. 636070

6. Insert the denture into position in the oral cavity. Guide the patient into occlusion, maintaining a proper relationship with the opposing arch. **Maintain the denture in a passive condition, without compression of the soft tissue, while the acrylic sets. Excessive occlusal pressure during the setting time may cause tissue recoil against the denture base and could contribute to dislodging and wear of the nylon males.**
7. After the acrylic resin has cured, remove the denture and discard the white spacer. Use a bur to remove excess acrylic, and polish the denture base before changing to the final male.


Fig. 4


Fig. 5

Acrylic Resin

8. Use the Locator Male Removal Tool (attached to the Locator Core Tool, Zest order 636066) to remove the Black Processing Male from the metal denture cap. The sharp circular edge on the end of the removal tool should be wedged tightly down into the very bottom of the plastic male so that it will catch the inside of the black plastic insert and pull it at an angle out of the metal housing (Fig. 6). To discard the nylon male from the tip on the Core Tool, point the tool down and away from you and tighten the Male Removal Tool clockwise back onto the Core Tool. This will activate the removal pin and dislodge the nylon male from the tip end of the Male Removal Tool.


Fig. 6  
Male Removal Tool/  
Locator Core Tool  
No. 636066

9. The Locator Male Seating Tool (contained in Locator Core Tool, 636066) is used to firmly push a Locator Replacement Male into the metal Denture Cap (Fig. 7). The replacement male must be placed securely in level with the rim of the cap (Fig. 8).


Fig. 7  
Male Seating Tool/  
Locator Core Tool  
No. 636066

**NOTE: The Replacement Male will not stay on the tool when it is turned upside down due to the varying sizes of males available. It is best to hold the denture with the base side down and snap the male into the metal denture cap.**

10. Instruct the patient in the path of insertion. Have the patient insert and remove the appliance several times.

Fig. 8

### LOCATOR MALE PLACEMENT BY THE LABORATORY

In the Operatory:

1. Insertion of the proper Locator Implant Abutment at tissue level must be completed (see Section A-1) before beginning the following impression procedure.
2. Place a Locator Impression Coping with Black Processing Male (636067) onto each Locator Abutment (Fig. 9).
3. Take an impression using a firm body impression material, exercising caution not to compress the soft tissue. The Locator Impression Coping is designed with minimum retention to be picked up with the impression material.
4. Snap a Locator Female Analog (636068 – 4 mm) onto each Impression Coping in the impression. The analog female must not fall off when turned upside-down with vibration (Fig. 10).


Fig. 9  
Impression Coping  
No. 8505

**NOTE:** An alternative reline impression technique using the patient's prosthesis is possible with use of the Locator Black Processing Cap Male (636070 Set). When the impression is withdrawn, the Processing Cap Male will remain on the abutment. Remove the Processing Cap Male from each abutment and snap it onto a Locator Female Analog. Reposition this assembly back into the impression making sure it is fully seated.

5. Pour the master cast. Upon separation, the Locator Female Analog is a part of the master cast replicating the position of the Locator Implant Abutment in the oral cavity.
6. Before waxing and processing the appliance, place a Locator Cap with Black Processing Male into each Female Analog in the master cast (Fig. 11). Make sure the male is fully seated.
7. Set the teeth and wax the appliance. Proceed with the processing technique of your choice through the boil-out step.
8. After the boil-out, remove the Processing Cap Male. Place a White Block-Out Spacer over the head of each Female Analog. The spacer is used to block out the immediate area surrounding the Locator Implant Abutment. The space created will allow the full resilient function of the pivoting metal denture cap over the Locator Male.


Fig. 10  
Female Analog  
Impression  
No.

9. Re-insert the Locator Black Processing Cap Male into each Female Analog, leaving the White Block-Out Spacer beneath it. The Black Processing Male will maintain the overdenture in the upper limit of its vertical resiliency during the processing procedure.


Fig. 11  
Denture Base  
Processing Model  
No.

10. Complete the processing and discard the white spacer. Avoid damage to the final male by polishing the denture base before changing to the final male.
11. Use the Locator Male 636066 Removal Tool attached to the Locator Core Tool to remove the Black Processing Male from the metal denture cap. The sharp circular edge on the end of the removal tool should be wedged tightly down into the very bottom of the plastic insert so that it will catch the inside of the black plastic insert and pull it at an angle out of the metal housing (Fig. 12).
12. The Locator Male Seating Tool is used to firmly push a Locator Replacement Male into the empty metal denture cap (Fig. 13). The replacement male must seat securely into place, level with the rim of the cap (Fig. 14).

**NOTE: The Replacement Male will not stay on the tool when it is turned upside down due to the varying sizes of males available. It is best to hold the denture with the base side down and snap the male into the metal denture cap.**


Fig. 12  
Male Removal Tool/  
Locator Core Tool  
No. 636066

#### HOW TO CHANGE THE LOCATOR MALE

1. The Locator Core Tool (636066), which contains a Locator Male Removal Tool and Locator Male Seating Tool, is used to remove the nylon male from the metal denture cap and replace it with another Locator Replacement Male.
2. Use the Male Removal Tool attached to the Locator Core Tool to remove the nylon male from the metal denture cap. The sharp circular edge on the end of the removal tool should be wedged tightly down into the very bottom of the plastic male so that it will catch the inside of the plastic insert and pull it at an angle out of the metal housing.
3. The Male Seating Tool is used to firmly push a Locator Replacement Male into the empty metal denture cap. The replacement male must it securely in place, level with the rim of the cap. Use of multiple Locator attachments (3 or more) in the same dental arch may require use of the 3.0 lbs. (light retention) pink colored Replacement Male No. 636072 or 1.5 lbs. (extra light retention) blue colored Replacement Male No. 636076, for easier removal of the prosthesis by the patient.

**NOTE: The Replacement Male will not stay on the tool when it is turned upside down due to the varying sizes of males available. It is best to hold the denture with the base side down and snap the male into the metal denture cap.**


Fig. 13  
Male Seating Tool/  
Locator Core Tool  
No. 636066


Fig. 14

#### PATIENT CARE

Good oral hygiene is vital to implant success. The Locator Implant Abutments must be thoroughly cleaned each day to prevent wear of the abutments due to buildup of abrasive plaque in the socket of the abutment. The use of a soft nylon bristle or end-tufted toothbrush, and superfloss to polish the abutments should be taught. A non-abrasive gel toothpaste, and an irrigation system is recommended to keep the socket of the Locator Abutment clean.

Patients should maintain a three to four month recall for cleaning and implant evaluation. The inside socket of the Locator Abutment and the sulcus area around the implant abutment are the primary areas of concern. Use plastic instruments for scaling the abutments. Do not use metal instruments which may create scratches on the abutment surface. Examine patients for signs of inflammation around the implant abutments, and for implant mobility. Use a 30 Ncm torque wrench to make sure the Locator Implant Abutment is tight before dismissal.

We continually strive to improve our products, and therefore, reserve the right to improve, modify or discontinue products and components at any time without notice or incurring obligation. Purchaser assumes all risks and liability resulting from the use of our products, whether used separately or in combination with other products not of our manufacture.

# HSL Abutment Connectors

## Description:

The IMPLA HSL Abutment Connector is used for the production of individualized, occlusally screw-fixed abutments in the cast-on technique. The abutment is composed of a base and a plastic tube (screw canal). The base is manufactured of a high-fusing cast-on alloy and the plastic tube burns out completely residue-free. During the cast-on procedure, the base is covered by the molt of a high-precious casting alloy or ceramic alloy. The metallic link is formed by diffusion of the components of both alloys.

## Important notes / warnings:

IMPLA HSL Abutment Connectors are a part of the IMPLA system and must only be used with the proper original IMPLA components (IMPLA implants). For the insertion of the abutments, use only the original IMPLA instruments and tools. Otherwise, a damage of the original components cannot be ruled out.

## Indications:

Single tooth gaps, partially or fully edentulous jaws, cement-fixed crowns or bridges, screw-fixed crowns.

## Contra-indications:

Do not fix inappropriate abutments on IMPLA implants. Ensure an accurate selection of the abutment version (diameter of the connector and implant geometry). Manipulating the connector geometry will lead to an imprecise fit which will render the abutment useless.

## Delivery:

All HSL Abutment Connectors are supplied in unsterile packaging. Therefore, the components must be sterilized by the practice staff accordingly before use.

## Procedure:

Modelling	The plastic tube can be trimmed and modelled as customary in dental technology. Ensure the preservation of the connector geometry to the HSL-base and keep the connector geometry of the abutment absolutely clean to avoid cast-on in the area.
Investing	Ensure that the investment remains bubble-free. Only use phosphate-bonded investment materials when casting high-precious metals.
Preheating	In order to allow for a complete cast-on, please adhere to the instructions of the preheating furnace as well as to the preheating temperatures and holding times prescribed by the manufacturer of the alloy.
Casting	Please adhere to the instructions of the casting machine. To prevent the IMPLA HSL Abutment Connectors from melting, please ensure that the casting temperature of the ceramic alloy you use is significantly lower than 1400°C (solidus temperature).
Cooling	We recommend letting the casting cylinder cool to room temperature to avoid tension in the cast object.
Divesting	To avoid a deformation of the object, especially of the connector geometry to the implant, only blast with beads at low pressure. Afterwards, the object can be stripped in an appropriate stripping solution.
Finishing	We recommend using ceramic grinders and hard-metal mills. Adhere to industrial safety regulations (protective goggles, extraction system, respiratory protection)!
Veneering	Please note the coefficient of thermal expansion of each veneering material to prevent cracking.

The abutments are fixed with the corresponding Aesthura vertical screw at 20 Ncm.

**Composition of the IMPLA HSL-Base:**

Gold	60 %
Platinum	24.9 %
Palladium	15 %
Iridium	< 1 %

**Technical data:**

Density	18.1 g/cm <sup>3</sup>
Melting interval	1350 – 1460°C
Vickers hardness HV5/30	150 – 205 – 230

**Side effects:**

Allergies and sensitivities towards the alloy cannot be ruled out completely and may occur in very rare cases.


**Interaction:**

In case of occlusal contact, the presence of different alloy types in the patient's mouth may lead to galvanic reactions.

# Multi-Unit Abutments

## Description:

The IMPLA Multi-Unit Abutment Concept facilitates improved usage of the available bone structure. By means of an angulated insertion of the terminal implants, the available bone is used in distal direction which improves the stability. With the according angulation, the Multi-Unit Abutments ensure a uniform insertion direction of the prosthetics to facilitate a tension-free screw fixation. For this purpose, the IMPLA Multi-Unit Abutment Concept contains straight and angulated components.

## Components Multi-Unit System

Abutments with the diameters 3.3 mm, 4.2 mm, 5.3 mm, each with gingiva heights of 1 mm and 3 mm, and angulations of 0°, 20° and 30° (not 5.3 mm).

## Accessories IMPLA Multi-Unit System

Impression system with impression aid and fixing screw as well as the according lab implant.

Gingiva sleeve for conditioning of the mucosa as well as the appropriate metal or plastic sleeve for the prosthetics with suitable screw for secondary screw fixation.

## Important notes / warnings:

IMPLA Multi-Unit Abutments are a part of the IMPLA system and must only be used with the proper original IMPLA components (IMPLA implants). For the insertion of the abutments, use only the IMPLA instruments and tools. Otherwise, a damage of the original components cannot be ruled out.

## Indications:

The IMPLA Multi-Unit Abutments are exclusively meant for a use with occlusal screw-fixed prosthetics (e. g. bridges, hybrid prostheses and cast bars).

A minimal implant length of 11.5 mm in the distal position is indicated for the use of IMPLA Multi-Unit Abutments.

The minimal number of implants for a complete reconstruction is four in the mandible and six in the maxilla.

A sufficient bone quality is a basic requirement for implantation, as it ensures the primary stability. The proper torque is between 30 and 40 Ncm.

The distal implants must be inserted with an angulation that corresponds to the angulation of the abutments inserted later, i.e. **0°, 20° or 30°**.

Presurgical planning as well as a comprehensive medical history and diagnosis (including radiological examination) are absolutely necessary.

In principle, the general rules and instructions for surgical and prosthetical procedures in dental implantology apply.

**Contra-indications:**

When selecting the patients, it is extremely important to ensure that they are able to maintain the necessary oral hygiene. In addition, the general rules and limitations in surgery and prosthetics for the selection and treatment of patients apply.

**Delivery:**

All components of the IMPLA Multi-Unit Abutment Concept are supplied in unsterile packaging. Therefore, the abutments must be sterilized accordingly before use.

**Procedure:**

Based upon presurgical planning and following the drill protocol for IMPLA implants, a pilot drilling for a mesial implant without angulation is performed first.

Subsequently, all other implants can be aligned to this first drilling with help of a drilling template. The drilling template improves the spacial orientation and presets the angulation of the distal implants. The drill is guided through a sleeve in the drilling template into the bone in the exact angle intended.

After the cavities have been prepared in accordance with the IMPLA drilling protocol, the implants can be inserted.

The prosthetic restoration can be carried out either in two stages or by immediate load.

In case of a two-stage restoration, we recommend to cover the implant with a surgical cover screw (covered healing) or with a gingiva former for open healing.

A primary stability of 30–40 Ncm is absolutely necessary for the immediate load approach.

Impressions are taken with help of the according impression components of the IMPLA Multi-Unit Abutment Concept.

Depending on the previously selected implant angulation (**0°, 20° or 30°**), the mesially positioned implants are then restored with a straight abutment and the distally positioned implants with an angulated abutment. The abutments are fixed with the appropriate vertical screw with a torque of 20 Ncm.

Finally, the prosthetic construction is set upon the abutments tension-free and is fixed with a prosthetic screw with a torque of 15 Ncm.

**Material:**

The IMPLA Multi-Unit Abutments are made of titanium grade 5 (Ti6Al4V ASTM 136).

**Side effects:**

Allergies and sensitivities towards the alloy cannot be ruled out completely and may occur in very rare cases.

**Interaction:**

In case of occlusal contact, the presence of different alloy types in the patient's mouth may lead to galvanic reactions.

# IMPLA PreFace Abutments

## **Description:**

The IMPLA PreFace Abutment is finished with an appropriate milling system with help of the original PreFace Abutment holder. The PreFace Abutment holder serves to fix the PreFace Abutment inside the milling system in the appropriate working position. The PreFace abutments are held in the PreFace abutment holder on the front side, making sure that they are seated exactly in the end position. The PreFace abutment holders are available for a wide range of milling machines.

## **Important notes / warnings:**

IMPLA PreFace abutments are part of the IMPLA system and must only be used with original components of the system (IMPLA implants). Process the PreFace abutments only with the corresponding suitable abutment holders.

## **Indications:**

The IMPLA PreFace Abutments are used for the production of individual titanium abutments on IMPLA implants in CAD/CAM milling machines. The individual abutments can be combined with crowns or bridges to reconstruct function and aesthetics.

## **Contraindications:**

The IMPLA PreFace Abutments must only be combined with IMPLA implants. Manipulating the connector geometry will lead to an imprecise fit which will render the abutment useless.

## **Delivery:**

All components of the IMPLA PreFace Abutment concept are supplied in unsterile packaging. Therefore, the components must be sterilized accordingly before use by the practice staff.

## **Procedure:**

The shape of the PreFace Abutment can be adjusted individually to the present anatomic conditions. The production of the PreFace Abutment is done digitally with an appropriate CAD software. The basic rules must be adhered to.

Design the emergence profile gently. Any further expansion of the gingiva than the expansion obtained with the gingiva former should always be discussed with the attending dentist beforehand. The gingiva can be expanded gradually by individually designed temporary abutments.

In general: The circular notch should be positioned slightly towards the subgingival direction in the vestibular area and orally on the gingiva line. This way, the cement between coping and abutment is easy to remove and the transition is not visible from vestibular direction. The minimal wall thickness of the abutments should not fall below 0.5 mm (screw hole to the outer surface of the abutment). The diameter of the emergence profile should not exceed 13 mm, while the height of the emergence profile should not exceed 6 mm. The total length of the PreFace Abutment should not exceed 15 mm, the maximum angle of the PreFace Abutment should not exceed 30°. The taper of the PreFace Abutment should not exceed 6°. The PreFace Abutment should be rounded occlusally and should taper off. Avoid sharp edges.

The recommended torque for the fixation of IMPLA PreFace Abutments is 20 Ncm.

## **Finishing:**

In case that further work steps are necessary after shaping the PreFace Abutment in the CAD/CAM milling machine, please note the following:

To mill PreFace Abutments, use appropriate, flawless mills with only little pressure. The minimal wall thickness must not fall below 0.5 mm. Avoid ridges and edges. We recommend to fix the PreFace Abutments in the lab analog to protect the connection geometry.

## **Scan body:**

A scan abutment is used to determine the exact position of the implant during the scanning procedure. The scan abutment's diameter must correspond to the diameter of the selected implant.

\* PreFace is a registered trademark of Medentika GmbH

# IMPLA Aesthura Abutments

## **Description:**

The IMPLA Aesthura abutment is used to fabricate occlusally screw-retained restorations. It is screwed onto the intended implants with the vertical screw, art. no. 638665.

## **Important information / warning:**

The Aesthura Abutment must not be altered by measures of any kind.

IMPLA Aesthura Abutments are part of the IMPLA system and must only be used with original components of the system (IMPLA implants). For the insertion of the abutments, use only the IMPLA instruments and tools. Otherwise, a damage of the original components cannot be ruled out.

## **Indications:**

Individual gaps, partially or completely edentulous jaws, screw-fixed individual tooth restorations.

## **Contra-indications:**

When selecting the patients, it is extremely important to ensure that they are able to maintain the necessary oral hygiene. In addition, the general rules and limitations in surgery and prosthetics for the selection and treatment of patients apply.

## **Delivery:**

All IMPLA Abutments are supplied in unsterile packaging. Therefore, the components must be sterilized accordingly before use by the practice staff.

## **Material:**

The IMPLA Aesthura Abutments are made of grade 5 titanium (Ti6Al4V ASTM 136).

## **Side effects:**

Allergies and sensitivities towards the alloy cannot be ruled out completely and may occur in very rare cases.

## **Interaction:**

In case of occlusal contact, the presence of different alloy types in the patient's mouth may lead to galvanic reactions.

## **Procedure:**

The abutments are fixed with the corresponding Aesthura vertical screw at 20 Ncm.

# Disinfection, cleaning and sterilization of surgical and prosthetic tools and instruments of the IMPLA system

Following surgical and prosthetic instruments/tools of the IMPLA system are composed of surgical steel/ titanium/PEEK/PP and are to be treated accordingly:

X-ray spheres, drills, drill extensions, countersinks, bone thread cutters, ratchet, screw drivers, IPK, depth stop, gingiva cutter, contra angle, flat wrench, sleeve holder, tray, carrier and guide sleeves, depth gauge, tweezers, sealing and healing cap, gingiva former, abutments, scan abutments, impression aid, horizontal rosette, vertical screw, bar screw, bone screws

<b>Please observe!</b>	Neutralize and dry the instruments thoroughly in order to avoid precipitation on the instruments to be sterilized. It is recommended to autoclave unmixed materials only.
Disinfection	Always wear protective gloves when disinfecting components suited to this treatment. Avoid drying of sticking residues of blood and/or tissue. Always use an appropriate disinfection bath and dose according to the manufacturer's instruction.
Cleaning	The components can be cleaned with all cleaning agents suitable for reprocessing surgical instruments and are resistant to any cleaning process, e.g. ultrasonic bath. Please avoid the use of metal brushes and abrasives. After cleaning, subsequently clean pH-neutrally and rinse off the cleaner and particle residues with demineralized water. Dry immediately.
Sterilization	The components are resistant to any kind of sterilization process. Make sure that all deposits have been removed by cleaning. The parts may not have any contact during sterilization. The temperature should not exceed 176°C (248°F). Please observe the instructions of the autoclave manufacturer.

### Instructions for the practice:

1. Do not let surgical residues (blood, tissue, etc.) dry on the instruments/parts. Submerge the instruments/parts in disinfective agent directly after the surgery.
2. Clean the instruments/parts only with a nylon brush. (Do not use metal brushes, such as brass brushes, etc.)
3. Rinse off disinfectant residue thoroughly with clear water.
4. Dry instruments/parts directly after cleaning and rinsing.
5. Never sterilize instruments made of different materials (e.g. ratchet) in the same cycle!
6. Please observe the instructions for use of the cleaning agents, as suitability for all instruments is limited (see material groups).

### 1 Use of the torque ratchet:

The torque ratchet is build up of the ratchet head (Fig. 1-b), the ratchet handle (Fig. 1-c) and the setting screw (Fig. 1-e). Inside the ratchet head, there is a tool fitting (Fig. 1-a).


Fig. 1

The ratchet head is lockable to lock the torque limit (see 1.1 below). If the ratchet head is not locked, the torque can be adjusted infinitely between 10 and 30 Ncm with the setting screw (see 1.2).

Before its first use and after each use in a patient case, the ratchet must be disassembled, cleaned and dried. **It must be sterilized directly before each use** (see section 3).


Fig. 2


Fig. 3

Fig. 4


Fig. 5

#### 1.1 Adjusting the ratchet

In the factory setting, the ratchet head is not locked. If you would like to lock it, turn the setting screw to less than 10 Ncm (Fig. 4), to relieve all pressure on the spring.

Subsequently, pull the handle from the ratchet head (Fig. 2). Then, turn the handle clockwise by 90° and let it lock into place in the ratchet head. Finally, preset a torque of 20 Ncm to ensure that the ratchet remains locked.

**Attention:** In the fixed condition, the maximum load is 100 Ncm!

#### 1.2 Setting the torque function

First, ensure that the ratchet is not locked. Turn the setting screw until you have reached the necessary torque. The correct torque is set once the according marking ring at the end of the handle is just still visible (Fig. 4). The ratchet head will now bend once the preset torque has been reached (Fig. 5). Afterwards, do not keep on turning, as this could lead to exceeding the preset torque significantly.

#### Important notice:

Whenever the ratchet is not in use, turn back the setting screw to relieve all pressure from the spring.


#### 1.3 Change tool fitting

(See 2.1 [1.] and 2.2 [5.]

## 2 Maintenance of the torque ratchet

### 2.1 Disassembly

1. The tool fitting (Fig. 6-a) is removed from the ratchet head (Fig. 6-b) by pulling on the two pin ends (arrows) towards the ratchet handle. Please note that the tool fitting can only be pulled from the ratchet head towards one side. **Do not use force.**
2. Completely screw off the setting screw (Fig. 6-e) from the ratchet handle (Fig. 6-c).
3. Pull off the ratchet handle (Fig. 6-c) from the ratchet head (Fig. 6-b) and remove the spring.
4. Disinfect, clean and dry all parts (see section 3).  
Afterwards and **before the sterilization**, the ratchet is reassembled.

### 2.2 Reassembly

1. Before reassembling the ratchet, all parts must be disinfected, cleaned and dried completely. Cover the areas marked with ● with oil (see 3.4).
2. Push the ratchet handle (Fig. 6-c) over the threaded bolt into the ratchet head (Fig. 6-b).
3. Push the spring (Fig. 6-d) into the ratchet handle (Fig. 6-c).
4. Screw the setting screw into the ratchet handle (Fig. 6-e) until the resistance of the spring becomes evident.
5. The tool fitting (Fig. 6-a) is set back into the ratchet head (Fig. 6-b) by pulling the two pin ends (arrows) towards the ratchet handle. Please consider that the tool fitting can only be set into the ratchet head from one side. **Do not use force.**

#### Important information:

Always relieve pressure on the spring when ratchet is not in use. Screw off the setting screw until the pressure on the spring is relieved completely.

## 3 Reconditioning

### 3.1 Disinfection

Before first use and immediately after each use, all ratchet parts must be disinfected in a suitable disinfectant. Adhere to the instructions of the disinfectant's manufacturer.

Chemical as well as thermal disinfection techniques are suitable, if they come up to the requirements of medical instrument disinfection.

### 3.2 Cleaning

Immediately after disinfection, all parts of the ratchet must be cleaned with a suitable cleaning agent. Adhere to the instructions of the cleaning agent's manufacturer. Otherwise, impurities will dry off onto the instrument and make cleaning difficult.

In order to avoid corrosion and discoloration, closely adhere to the instructions on concentration, submergence time and temperature given by the manufacturer of the cleaning agent.

Remove any deposits, impurities, etc. thoroughly with a soft brush (do not use metal brushes, scouring agents, wire wool or too much pressure).

#### Please note:

When disinfecting or cleaning the device with an ultrasound unit, please adhere to the instructions of the unit's manufacturer.

### 3.3 Rinsing and drying

Directly after disinfection, cleaning or ultrasound treatment, the individual parts must be rinsed carefully with clear desalinated or distilled water (do not use tap water). Afterwards, dry all parts thoroughly to avoid water stains.

If the individual parts are not rinsed, residue is not removed and/or the parts are not dried completely, this could lead to discoloration or corrosion of the device during sterilization.

### 3.4 Care and checking

Parts that have joints must be treated with a sterilizable oil spray (paraffin-based) after all cleaning and disinfection procedures. Such a spray will reduce metal-on-metal-friction, avoids fretting corrosion and reduces the danger of oxidation during heavy use.

### 3.5 Sterilization

Only sterilize carefully cleaned, dried devices without any type of corrosion.

Badly cleaned and insufficiently dried parts bear the danger of burn-in or corrosion damage.

If corroded devices are sterilized with immaculate instruments, flash rust can cause damages on the immaculate instruments.

Please note the instructions of the manufacturer.

#### Note:

Only sterilize the torque ratchet directly before use.


### 3.6 Storage

#### Important information regarding the storage of the ratchet:

The spring of the ratchet should be fully relieved immediately after use to ensure proper use in the long term. Only store the ratchet in a relieved state.

## Cleaning and sterilization of the ratchet

Screw off turn-knob, remove spring and handle.

To clean, pull back the interlock to unlock the receptacle

Sterilization of the ratched only in disassembled state, without the spring.  
Attention: Autoclave spring separately!

# Application and safety instructions

All rotating instruments and parts for medical use (dental implantology, surgery) were developed for their specific application. For that reason, incorrect use may cause premature wear of the instruments and endanger the patients.

Rotating instruments for the medical field must only be used by physicians and other experts who are intensively familiar with the use of these instruments due to their training and experience.

## Correct application

- Take care to use only technically and hygienically flawless, well-maintained and cleaned contra-angle handpieces.
- The instruments must be inserted correctly.
- The instrument must revolve at operating speed before being placed on the object (see pages 4, 5, 10, 14).
- Avoid wedging or levering the instrument as this may increase the risk of fracture.
- Thermal damage through the use of rotating instruments must be avoided (we recommend the slow-speed mode with sufficient cooling). Please see Surgical procedure or Procedure on pages 4, 5, 8, 10 and 14 of this manual.

## Pressing forces

- Excessive pressing forces must be avoided at all times. Excessive pressing forces may damage the working sections of the instruments or fracture their blades. This also results in an increased generation of heat.

## Cooling

- To prevent excessive heat generation during preparation, provide for adequate air or water cooling (at least 50 ml/min).
- Inadequate cooling may lead to irreversible damages of the bone and/or tissue.

# Manufacturer's information on rotating surgical instruments

## Product:

The present manufacturer's information applies to all instruments that are used for surgical application. These are rotating instruments made of stainless steel.

**Instruments delivered in non-sterile condition must be prepared prior to first use.**

## Limited number of reprocessing cycles:

The end of product life is basically determined by wear and damage from use. If necessary, known restrictions on the frequency of use of instruments must be observed (see pages 5, 10, 15). Frequent reprocessing does not affect the performance of these instruments.

## Work station:

Hygienic precautions according to the provisions valid in your country.

## Storage and transport:

Place instruments in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). It is recommended to process the instruments within one hour of use at the very latest. The instruments should be transported to the reprocessing site in the disinfection tank. Endodontic instruments can be transported in a special interim support equipped with a foam insert drenched in disinfecting solution.

# Standardized manual reprocessing for rotating instruments

## Equipment used:

- Nylon brush
- Suitable detergent/disinfectant for rotary instruments with proven disinfecting effect
- Ultrasonic bath

## Reprocessing:

1. Rinse off surface contamination under running water. Remove stubborn contamination with a nylon brush under running water, turning the instrument constantly.
2. Rinse instrument thoroughly with running water.
3. Place the instruments in a suitable sieve into the ultrasonic device filled with detergent/disinfectant.
4. Observe the instructions of the manufacturer regarding concentration and immersion time for chemical cleaning/disinfection in the ultrasonic device. Be sure to observe the full correct immersion time which does not start until the last instrument has been placed into the ultrasonic device. Attention: Do not exceed 45°C (113°F) (risk of protein coagulation)!
5. On completion of the immersion time, rinse instruments thoroughly with suitable water (preferably with demineralised water to avoid residues).
6. Dry instruments (according to RKI recommendations, preferably with compressed air).
7. Visual examination with a suitable magnifying device to ensure that the instrument is clean and undamaged (experience has shown that a magnification factor of 8 permits a visual examination). If there are still visible residues of contamination, repeat the cleaning and chemical disinfecting process until no visible contamination is left.

**Control and functional test:**

Instruments showing the following defects are to be discarded immediately:

- Blunt and chipped blades
- Deformations (e. g. bent/twisted/fractured instruments)
- Corroded surfaces

**Packing:**

Make sure that the packaging is suitable for the instrument and the chosen method of sterilization. Single pack: The packaging must be large enough to ensure that there is no pressure on the seal. In the set: Place instruments onto the tray provided or onto universal sterilization trays. The instruments must be protected. Use an appropriate method to pack the tray. Instruments with limited times of use are to be marked accordingly.

**Sterilization:**

Steam sterilization using a vacuum process at 134°C (273°F) in a device that complies with the provisions of DIN EN 13060; with validated processes.

- Fractionated pre-vacuum (type B) respectively simplified pre-vacuum (type S).
- Sterilization temperature: 134°C (273°F)
- Hold time: 5 minutes minimum (full cycle)
- Drying time: 10 minutes minimum

In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer.

**Transport and storage:**

The packed sterile goods must be protected from dust, moisture and recontamination during transport and storage.

**General note:**

Observe the legal provisions regarding the reprocessing of medical products valid in your country (e. g. [https://www.rki.de/EN/Home/homepage\\_node.html](https://www.rki.de/EN/Home/homepage_node.html)).

The manufacturer confirms that the above detailed reprocessing methods are **suitable** for preparing the above named instrument group to enable their reuse. The person carrying out the reprocessing is responsible for ensuring that the applied method is carried out with appropriate equipment, materials and trained personnel at the reprocessing site and that it actually achieves the desired result. To guarantee this, routine controls of the validated mechanical and/or manual preparation methods are necessary. Any deviation from the above detailed process (e. g. use of different chemicals) must be carefully checked by the operator to ensure effectiveness and to avoid possible adverse consequences.

# Information on storage, disinfection, cleaning procedures and sterilization of rotating instruments

- Rotating instruments must be disinfected, cleaned and – whenever the instrument material permits it and it seems suitable – sterilized prior to first use on patients and immediately after use.
- Used instruments must always be disinfected before cleaning.
- Before first use, the rotating instrument should be kept in its original packaging and be stored at room temperature, protected against dust and moisture.
- All instruments and parts are packaged non-sterile and must be sterilized before use according to their application. Cleaning and sterilisation not in accordance with the instruction may lead to microbiological contamination and infection of the patient.
- All materials can be sterilized in an autoclave according to DIN ISO 17665 „Sterilization of health care products - Moist heat“ (e. g. at 134°C [273°F] with a minimum holding time of 3 min. or 121°C [249°F]/15 min.). Please be sure to follow the manufacturer’s operating instructions when sterilizing.
- Before sterilization, the instruments must be cleaned and disinfected properly.
- Disinfection may be carried out manually or by machine.
- In case of severe contamination of the instrument, ultrasound cleaning techniques are recommended.
- After cleaning rinse the instruments thoroughly with water and dry them immediately.
- Contact between the instruments during cleaning and disinfection procedures may cause damage and should be avoided.
- It is highly recommended to wear protective gloves when working with contaminated instruments.
- Thermal disinfection procedures are generally not suitable for rotating instruments.
- Avoid sterilization temperatures exceeding 180°C (356°F). Exceeding this temperature reduces the hardness of the working part and reduces the instrument’s service life.
- All sterilized instrument containers should be labeled with an indicator strip bearing the date of the sterilization as well as the expiry date.
- Non-rustproof instruments may corrode. This may cause discolorations and reduce the service life. Therefore, it is recommended to use anti-corrosion disinfectants and cleaning agents in these cases.
- The instructions for use, reaction time, and suitability of disinfectants and cleaning agents for certain types of instruments are provided by the manufacturer's instructions. They must be observed.
- Contact of the instruments with hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) must be avoided.
- The instruments should be stored in appropriate, hygienically well-kept containers. This also applies to sterilized instruments. When stored, the instruments should be protected against dust, humidity, and recontamination. Observe the maximum storage period.
- Broken-out cutting edges cause vibrations and high contact forces and thus lead to broken-out preparation edges and rough surfaces. They should be discarded immediately.
- Bent or non-concentric instruments should be discarded immediately.