# Capo Slow Flow

(hochviskoses Flowable Composite)

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

CE 0297

Schütz Dental GmbH  
Dieselstr. 5-6 • 61191 Rosbach/Germany  
Phone: +49 (0) 6003 814-0 • Fax: +49 (0) 6003 814-106  
info@schuetz-dental.de • www.schuetz-dental.com

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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### Instructions for use

Capo Slow Flow is a light-curing, flowable, highly radiopaque (210% Al) composite with a high viscosity. The guidelines and requirements of DIN EN ISO 4049 apply. Capo Slow Flow Composite is available in syringes and compules. The compules are for single use only. Please do not reuse them, as this makes it impossible to rule out contamination and germ formation.

**Composition:**  
Glass powder, diurethane dimethacrylate, tetramethylene dimethacrylate, silicon dioxide

**Total filler:**  
77 % by weight (57 % by volume) anorganic filler (0.005 – 40 µm)

### Indications:

- Fissure sealing
- Extensive fissure sealing on molars and premolars
- Fillings in Black's class V (cervical caries, root erosions, v-shaped lesions)
- Minimally invasive fillings in Black's classes I, II and III
- Corrections of enamel defects
- Blocking out undercuts
- Small corrections of shape and color of the enamel

### Procedure:

**Pre-treatment**  
Before starting the treatment, clean the hard tooth substance with a fluoride-free polishing paste. Select the correct color with help of the Vita® shade guide while the tooth is still moist.

- Cavity preparation**  
Prepare the cavity according to the general rules of the adhesive technique. Use a procedure that is gentle on the hard tooth substance. Bevel the incisal edges on anterior teeth. Do not bevel the edges on posterior teeth. Avoid

irregular edges. Next, clean the cavity with water spray, remove all residue and dry. It is necessary to dry the site completely. We recommend the use of a kofferdam.

- Pulp protection/liner**  
The use of a liner can be foregone, if an enamel dentin adhesive is used. In case of very deep cavities that come close to the pulp, cover the according surfaces with a calcium hydroxide compound.
- Design of the approximal contacts**  
For cavities in approximal areas, insert and fix a transparent matrix.
- Adhesive system**  
Etch (e.g. with Capo Etch) and bond (e.g. with Capo Bond) according to the manufacturer's instructions.
- Application of Capo Slow Flow syringes**  
Fill Capo Slow Flow in thin layers (max. 2 mm) directly into the cavity. For hygiene reasons, do not use the application tips more than one time. Avoid introducing air bubbles into the material during application with the tip. Ensure that the prepared hard tooth substance is well-coated.
- Application of Capo Slow Flow compules**  
Insert the compule into the dispenser. Remove the sealing cap. Place the compule in the correct angle towards the cavity. Inject the material into the cavity. Apply slow and steady pressure to the compule. Do not use excessive force! Max. layer thickness 2 mm. To remove the compule from the dispenser after use, retract the plunger. Next, remove the compule. Please note: For hygiene reasons, compules are intended for single use only.
- Polymerization**  
Polymerize each layer with a commercially available polymerization unit for 40 seconds, holding the light guide as close as possible to the surface. During polymerization, a dispersion layer forms on the material surface. This layer must not be touched or removed. If forms a chemical bond with the following layer of composite applied.
- Finishing**  
Capo Slow Flow is ready to be finished and polished directly after polymerization. Suitable tools are finishing diamonds, flexible wheels, silicone polishers and polishing brushes.

### Special notes:

- In case of extensive restorations, the surgical light should be moved away from the working area temporarily to avoid premature curing of the composite. Alternatively, the material can be covered with a light-tight foil.
- For hygiene reasons, do not use the curved application tips included with the product more than once. Dispose of the tips directly after use.
- Use a light polymerization system with an emission range of 350 – 500 nm to polymerize the material. The required physical properties are only reached if the polymerization light functions properly. Therefore, it is necessary to check the light intensity regularly according to the manufacturer's instructions. Light intensity for polymerization  $\geq 650 \text{ mW/cm}^2$   
Wave length for polymerization 350 – 500 nm  
Polymerization time 40 sec.

### Hazard statement / precautionary statement: Contains tetramethylene dimethacrylate, di-urethane dimethacrylate

**Warning:** May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects. Avoid breathing vapours / spray. Avoid release to the environment. Wear protective gloves. If skin irritation or rash occurs: Get medical advice / attention.

### Contra-indications

If a patient has known hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth before using the product.

### Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us. In case of exposed dentin in a cavity, ensure a proper protection of the pulp to avoid possible pulp reactions (e.g. line the cavity with a calcium hydroxide compound).

### Cross-reactions with other substances:

Phenolic substances (e.g. eugenol) inhibit polymerization. Therefore, these materials must not be used as liners.

### Disinfection / Protection from cross-contamination:

Place the functioning syringe with attached application tip into a suitably shaped barrier sheath; pierce end of sheath with cannula, exposing the cannula for use. Using a barrier sheath facilitates cleaning and disinfection of the syringe between patients. After use of sheathed syringe, remove delivery tip and sheath by grasping on the hub of the delivery tip through the sheath, twist and remove tip along with sheath. Discard used tip and sheath in appropriate waste stream. Replace syringe storage cap. Disinfect – After removal of the delivery tip and sheath, disinfect this product using an intermediate level disinfection process (liquid contact) as recommended by the Centers for Disease Control and endorsed by the American Dental Association, Guidelines for Infection Control in Dental Health-Care Settings – 2003 (Vol. 52; Nr. RR-17), Centers for Disease Control and Prevention (USA).

### Storage:

Store at temperatures between 10 °C and 25 °C (50 °F to 77 °F). Screw the cap back onto the syringe tightly after each use. Let the material reach room temperature before use. Withdraw the plunger slightly after use to keep the outlet from becoming plugged. Do not use after the expiration date (see label on the syringe). Only for use in dentistry. Keep out of children's reach. This product has been developed for the specific use illustrated above. Only process as described in these instructions. The manufacturer will not be held liable for any damages that result from improper use or improper processing.

\* Vita is a registered trademark of Vita® Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen.

### Troubleshooting

Trouble	Cause	Remedy
Composite does not cure	Luminous intensity of the polymerization unit insufficient Emitted spectral range of the polymerization unit insufficient	Check luminous intensity; replace light source, if necessary Consult manufacturer of polymerization unit; recommended spectral range: 350 – 500 nm
Composite seems to be too hard / firm inside the syringe	Material was stored at temperatures below 10 °C for a longer period of time Syringe was not closed tightly which caused part of the material to cure	Let composite reach room temperature before use Close syringe correctly with the cap after each use
Composite does not cure sufficiently	Layer thickness per polymerization cycle too high	Keep to max. layer thickness of 2 mm
Restoration seems too yellow when compared to color reference	Insufficient polymerization of the composite layers	Repeat polymerization cycle several times, for a minimum of 40 sec.

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