FuturaBasic Hot

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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Instructions for use

Composition:

- Powder: pearl polymer made of polymethyl methacrylate, pigments, initiators
- Liquid: methyl methacrylate, tetramethylene dimethacrylate, initiators, stabilizers
- FuturaBasic Hot is free of cadmium.

Indications:

- For the fabrication of partial and full dentures using the flask technique.
- For the fabrication of splints
 - bite splints
 - reflex splints / occlusal splints

Contra-indications / interaction with other agents:

Do not use in case of a known allergy against one of the components.

Application / processing:

Preparatory steps:

Produce plaster model and wax-up as usual. The prefabricated teeth set up in the wax are fixed with plaster or silicone (stent/flask). Subsequently, the wax is eliminated as usual. Roughen the acrylic teeth, provide them with undercuts and coat them with monomer or HS Cross-Liquid (Schütz Dental) to ensure that they are securely retained. The surfaces of the plaster models are insulated with the alginate separating agent FuturaSep Plus (Schütz Dental) according to the manufacturer's instructions.

Mixing ratio:

2.5 parts powder by weight, 1 part liquid by weight approx. **25 g powder** and **10 g liquid** are required for a full denture. If desired, veins can be added. Please note the separate manufacturer's information and instructions.

Production of acrylic dentures:

After adding the powder to the liquid, stir FuturaBasic Hot thoroughly with a spatula for 30 seconds in a mixing cup to create a homogeneous mixture.

Allow a doughing time of ca. 10 min. Afterwards, the acrylic is ready for use. Knead the acrylic dough for ca. 30 seconds. The processing time is ca. 30 - 40 minutes at an ambient temperature of ca. 21°C (ca. 69°F). Once FuturaBasic Hot becomes viscoelastic, pack it into the prepared flask and press. Hold the pressure for ca. 15 minutes. Secure the flask in a clamp for polymerisation.

Polymerize FuturaBasic Hot in the closed flask with the clamp for 20 minutes in boiling water (or in a pressure pot at 95°C). Allow the flask to cool to room temperature in the clamp and then deflask in the usual way.

Finishing the material:

After polymerization, the restoration is lifted from the model and it is finished with appropriate tools.

Silicone polishers (wheel, lentil, roll) as well as hard metal cutters or diamond grinding tools are suited for the finishing the material. Prepolish with sandpaper and pumice.

Polishing:

Polish the material with goat hair brushes, polishing paste and soft woolen wheels. A careful surface treatment and polish are a precondition for an optimal result.

Repairs:

Cold curing resins (e. g. FuturaGen, FuturaPress LT, FuturaSelf or FuturaPress N) can be used for repairs. Roughen the fracture surfaces, clean them and coat them with monomer before applying the acrylic dough. The polymerization time is 15 min. in a pressure polymerization system, at a pressure of 2 - 5 bar and a temperature of 45°C (113°F).

Production of splints:

Preparatory steps:

Production of a plaster model. Blocking out of undercuts with suitable materials. Prepare a wax model of the desired splint. Use plaster to produce a counter model and remove wax. Separate plaster model with suitable materials (e. g. FuturaSep Plus, Schütz Dental) according to the manufacturer's instructions.

Mixing ratio:

2.5 parts powder by weight, 1 part liquid by weight

Production of splints:

After adding the powder to the liquid, stir FuturaBasic Hot thoroughly with a spatula for 30 seconds in a mixing cup to create a homogeneous mixture. Allow a doughing time of ca. 10 min. Afterwards, the acrylic is ready for use. Add the acrylic to the model as required, close the flask. Secure the flask in a clamp and press. Hold the pressure for ca. 15 minutes. The flask is secured in a clamp for polymerization.

Polymerize FuturaBasic Hot in the closed flask with clamp for 20 min. in boiling water (or in a pressure pot at 95 °C [203 °F]). Let the flask with the clamp cool down to room temperature and deflask as usual.

Finishing the material:

After polymerization, the restoration is lifted from the model and it is finished with appropriate tools.

Silicone polishers (wheel, lentil, roll) as well as hard metal cutters or diamond grinding tools are suited for the finishing the material.

Repolish with sandpaper and pumice.

Polishing:

Polish the material with goat hair brushes, polishing paste and soft woolen wheels. A careful surface treatment and polish are a precondition for an optimal result.

Repairs:

Cold curing resins (e. g. FuturaGen, FuturaPress LT, FuturaSelf or FuturaPress N) can be used for repairs. Roughen the fracture surfaces, clean them and coat them with monomer before applying the acrylic dough. The polymerization time is 15 min. in a pressure polymerization system, at a pressure of 2 - 5 bar and a temperature of 45 °C (113 °F).

Troubleshooting / FAQ liste:

Problem	Cause	Remedy
acrylic does not cure properly	incorrect mixing ratio	<ul style="list-style-type: none"> – observe instructions for use – observe mixing ratio
acrylic shows a whitish tint along the bottom	<ul style="list-style-type: none"> – incorrect or improper insulation – model was not soaked in water 	<ul style="list-style-type: none"> – insulate properly – check insulation agent, if necessary – soak model in lukewarm water for ca. 10 min.
acrylic forms bubbles along the bottom	– model was not soaked in water	soak model in lukewarm water for ca. 10 min.
acrylic surface white / spongy	too much time passed until the acrylic was put in the pressure pot	observe processing times
color deviations of repairs	incorrect mixing ratio	observe mixing ratio
Kacrylic fractures when lifted	improperly insulated	check insulation
Kacrylic / prosthesis shows whitish streaks	<ul style="list-style-type: none"> – improperly mixed – powder content too high 	<ul style="list-style-type: none"> – mix acrylic dough thoroughly with a spatula – observe mixing ratio
splint fractures	design too thin / unsuitable	ensure a sufficient layer thickness

Storage:

Storage temperature at 10 - 25 °C / 50 - 77 °F.
Close container carefully.

Shelf-life:

The maximum shelf-life is printed on the label of each packaging. Do not use after the expiration date.

Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-indications / interactions:

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth before using the product.

Note:

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

Disposal:

Powder and liquid must be disposed of as hazardous waste.

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