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Interaction with other substances

As phenolic substances (such as eugenol) inhibit polymerization, do not use cavity liners (such as zinc-oxide eugenol cements) containing such substances.

Hazard statement**Contains tetramethylene dimethacrylate, aliphatic di-urethane dimethacrylate**

Warning: May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects.

Avoid breathing vapours/spray. Avoid release to the environment. Wear protective gloves. If skin irritation or rash occurs: Get medical advice/attention.

Application**Preparatory measures**

Before commencing the treatment, clean the tooth with non-fluoride polishing paste. Use a Vita* shade guide to select the shade while the tooth is still moist.

1. Cavity preparation

Prepare the cavity minimally invasively as generally required for adhesive techniques. All enamel margins in the anterior region must be bevelled. Do not bevel the margins in the posterior region and avoid slice preparations. Spray the cavity with water to clean it, remove all debris and dry it. The cavity must be isolated. It is advisable to place a rubber dam.

2. Pulp protection/Cavity liner

If an enamel-dentine adhesive is used, no cavity liner is required. In very deep cavities those areas in close proximity to the pulp must be coated with a calcium hydroxide material.

3. Approximal contact areas

When filling cavities with approximal sections, place a transparent matrix and fix it in place.

4. Etching and Bonding

Etch (e.g. CapoEtch) and bond (e.g. Capo Bond) according to the manufacturer's instructions.

5. Placing the composite in the cavity

"Fillis"

Insert the "Fillis" in the dispenser and remove the sealing cap. Exert pressure slowly and uniformly to dispense the required amount of filling material from the "Fillis" directly into the cavity, layer-by-layer. Contour the composite with standard metal instruments.

Please note: For hygienic reasons "Fillis" are for single-use only.

The layers must not be more than 2 mm thick. Due to the effect of the oxygen in the air, a thin smear layer of unpolymerized material remains on the surface of each layer. This bonds the layers chemically and must not be touched or contaminated with moisture.

Curing

All shades must be light cured for 40 seconds per layer using a commercially available light curing unit.

The light guide must be held as close possible to the surface of the filling. Multiple surface fillings should be light cured from each side.

6. Trimming

ReVeneer Composite can be trimmed and polished immediately after curing using finishing diamonds, flexible disks, silicone polishers and polishing brushes.

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Check the occlusion and articulation and spot grind to eliminate high spots or undesirable paths of articulation from the surface of the filling.

Indirect method**Cavity preparation**

The cavity should be prepared as minimally invasively as possible with only slightly diverging sides. To prevent the

**ReVeneer
Composite****Instructions for use****ReVeneer Composite**

ReVeneer Composite is a light curing, hybrid composite containing an ultrafine, radiopaque glass filler and is indicated for placing fillings using adhesive techniques. It can be polished to a high lustre.

Due to the ultra-fine particle filler, extremely homogeneous restorations can be placed which are easily polished to a high lustre. The chameleon effect matches the shade of the filling perfectly to the tooth structure. The guidelines in EN ISO 4049 have been complied with.

ReVeneer Composite is supplied in **single-use**, 0.3 g "Fillis". The Fillis must not be used more than one time, as a contamination of the material and a transfer of germs cannot be ruled out. For indications 2-6, use ReVeneer Composite in combination with the CapoEtch etching gel and CapoBond components.

Composition:**Monomer matrix:**

Aliphatic diurethane dimethacrylate, tetramethylene dimethacrylate, isopropylidene-bis[2(3)-hydroxy-3(2)-(4-phenoxy)propyl]bismethacrylate

Total filler:

75% by weight (52% by volume)
glass filler (mean particle size: 0.7 µm),
pyrogenic silicic acid (mean particle size: 0.04 µm)

Indications

- Intra and extra oral repair of fractured porcelain veneers (see 6. and 7.)
- Direct anterior and posterior restorations in Black's class I, II, III, IV and V cavities.
- Indirect restorations such as inlays, onlays and laminate veneers
- Extended fissure sealing in molars and premolars
- Cores
- Splinting mobile teeth
- Adjusting the contours and shade to improve aesthetics

Contra-indications / Interactions

If a patient has known allergens against or hypersensitivity towards a component of this product, we recommend not to use it to do so only under strict medical supervision. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth before using the product.

Unpolymerized composite may cause skin allergies. The user must take adequate precautions. In case of irritation or allergy due to one of the constituents listed under "Composition", do not use this material.

Placing the inlay, onlay or laminate veneer

Remove the temporary restoration and clean the cavity. Place a rubber dam before cleaning and drying the prepared

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*Vita è un marchio registrato della Vita® Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen