


Tizian Blank translucent

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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




Instructions for Use

Material Properties

Tizian Blank translucent discs are made of Zirconium Oxide (Y-TZP ZrO_2) for dental applications.

This material is specifically made for manufacturing of permanent dental prosthesis.

Applications include: Anterior, posterior crown and bridges, conical telescopic copings.

After completion of the specified final sintering all Tizian Blank translucent zirconia meets the requirements of EN ISO 6872.

Technical Data

Composition

Material / Product	Tizian Blank translucent
Zirconium oxide ZrO_2	$\geq 89\%$
Yttrium oxide Y_2O_3	4 – 6%
Hafnium oxide HfO_2	$\leq 5\%$
Aluminium oxide Al_2O_3	$< 1\%$
Chemical solubility [$\mu\text{g}/\text{cm}^2$]	< 100

Properties

Material / Product	Tizian Blank translucent
Linear thermal expansion / CTE [10^{-6} K^{-1}]	$10,1 \leq \text{CTE}^* \leq 11,1$ (*span 25 – 500° C)
Biaxial flexural strength [MPa] ¹	1270
Fracture Toughness [MPa*m ^{1/2}] ¹	≥ 5
Translucency [1-CR]*100	42% ²
Typ / Class	Typ II / Class 5

¹ Typical values acc. to EN ISO 6872 (polished specimen)

² Pre-shaded discs show lower values

Warranty / Storage

Warranty

- Technical information and user recommendations, whether given orally or in writing, as well as practical training are guidelines.
- Sintering ovens vary in their performance. It is extremely important that furnaces be calibrated on a regular basis to achieve optimum results. Follow the manufacturer's recommended calibration instructions.
- Our products are subject to continuous development and improvements. We will advise you of these changes.
- We strive to provide the best quality products. Upon receipt inspect the product for any visual defects prior to milling.
- After discs and blocks have been partially milled, complaint claims will be voided.

Storage

Store all Tizian Blank translucent zirconia in its original packaging in a dry environment at normal room temperature.

Application / Design and Milling

Indications For Use

Material / Indication	Tizian Blank translucent
Single unit restorations (anterior / posterior)	✓
3-unit bridges (anterior / posterior)	✓
Multi-unit bridges (anterior / posterior)	✓

Contraindications

- Insufficient tooth structure reduction.
- Insufficient tooth structure for proper adhesion and force distribution.
- Insufficient oral hygiene.
- Insufficient interproximal space for sufficient joints in bridges.
- Known allergies.
- Known incompatibilities to product composition.

Design and Milling

- Follow instructions for CAD/CAM software to scan and design restorations.
- Milling systems need to be calibrated for best results. All systems are not alike and can produce adverse results if the minimum thickness is not followed.
- For bridges, always design auxiliary supports to prevent warping during sintering.

Infiltration

Infiltration of Tizian Blank translucent restorations before sintering

- Please follow the respective Instructions For Use of the supplier.

WARNING!

- User must take precautions when handling green state zirconia. Always work in a well-ventilated environment.
- Use synthetic gloves when handling zirconia and pre-stain liquids.
- Pre-stain liquids should not come in contact with skin.
- It is recommended to wear gloves, safety goggles and suitable protective clothing.
- Use appropriate vacuum in a well-ventilated area to capture and contain dust.

Sintering Recommendations

Sintering ovens vary in their performance. It is extremely important that furnaces be calibrated on a regular basis to achieve optimum results.

	Number of units	Duration (h)	Phase	Temperature		Heating rate / Cooling rate		Holding time (min.)
				°C	°F	°C/min.	°F/min.	
Tizian Blank translucent								
Standard	1 – 5	- 3.7	Phase 1	20 – 1300	68 – 2372	30	54	30
			Phase 2	1300 – 1530	2372 – 2786	40	72	60
			Phase 3	1530 – 900	2786 – 1652	15	27	–
			Phase 4	900 – 80	1652 – 176	20	36	–
	5 – 10	- 5.2	Phase 1	20 – 1300	68 – 2372	30	54	60
			Phase 2	1300 – 1530	2372 – 2786	40	72	120
			Phase 3	1530 – 900	2786 – 1652	15	27	–
			Phase 4	900 – 80	1652 – 176	20	36	–
Long	1 – 20	10.7	Phase 1	20 – 900	68 – 1652	10	18	10
			Phase 2	900 – 1530	1652 – 2786	3	5.4	150
			Phase 3	1530 – 80	2786 – 176	8	14.4	–
	> 21	11.7	Phase 1	20 – 900	68 – 1652	10	18	10
			Phase 2	900 – 1530	1652 – 2786	3	5.4	210
			Phase 3	1530 – 80	2786 – 176	8	14.4	–
Over-night	unlimited	-14.3	Phase 1	20 – 250	68 – 482	2	3.6	–
			Phase 2	250 – 1530	482 – 2786	4	7.2	240
			Phase 3	1530 – 80	2786 – 176	8	14.4	–

- Place objects to be sintered on beads in sintering tray.
- Space objects in tray to allow for convection heat.
- Fired objects will have a slight luster.

**WARNING!**

- Sintering furnaces must be located in a fireproof well-ventilated area.
- Slow cooling is essential to the final outcome; do not cool too fast.
- Opening the furnace too early may cause zirconia to crack.

Sandblasting / Post Processing**Frame Fitting**

After final sintering, the zirconia restorations can be fit and shaped using suitable diamond grinding points. Use a water cooled lab turbine to prevent fractures.

Margins can be thinned using a soft rubber abrasive wheels, especially designed for such use.

Sandblasting

After performing any adjustments, the object should be lightly sandblasted with pure white 50 µm corundum (aluminum oxide) at approximately 2.5 bars.

Re-sintering

After sandblasting and steam cleaning, the objects should be re-sintered in a porcelain furnace to seal any micro fractures that may have developed during grinding.

Raise temperature at 40 °C/min. to 1000 °C. Hold in air for 5 minutes. Slow cool to room temperature.

Restorations are now ready for veneering, staining, and glazing.

**WARNING!**

- Any grinding performed on sintered zirconia should be carried out in well-ventilated areas.
- Do not inhale particle dust.
- Use appropriate vacuum units to capture dust.
- Use safety glasses when grinding and sandblasting.
- Sandblast only in approved units with vacuum.

Veneering / Staining and Glazing**Veneering**

- A thin wash of bonding porcelain should be applied to the veneering surface and fired.
- Apply zirconia veneering porcelain as required.
- Follow manufacturer's recommendations for firing parameters.
- Follow technical information for coefficient of thermal expansion for zirconia as well as veneering porcelain coefficient.

Staining and Glazing

- Stain and glaze in thin layers to preferred luster.
- Use stains and glazes designed for use with zirconia.
- Use manufacturer's recommendations for firing parameters.

Post Processing by Dentist

When occlusal and proximal adjustments are required by the dentist it is recommended that fine diamond grinding points be used. The restoration should be cooled during the grinding process. Diamond grit size should be approximately 40 microns.

After grinding, smooth areas with a rubber wheel and polish with 10 micron diamond polishing paste.

Note that if the restoration is not sufficiently polished, the antagonist may experience abrasion over time.

Cementation

Conventional Cementation

The inherent properties of Tizian Blank translucent zirconia give it maximum strength and stability. Therefore conventional fixation with zinc – oxide phosphate or glass ionomer cement is possible in most cases. Lightly sandblast internal of restoration with pure white 50 micron Aluminum Oxide and steam clean prior to cementation.



When applying the conventional cementation technique, it is important to observe the correct requirements of abutment retention.

Adhesive Fixation

For adhesive fixation, we recommend the bonding composite. These adhesive cements will create an excellent bond between tooth structure and the zirconium-oxide frame material.

Zirconia Fixation as a Provisional

Although not recommended, if a restoration needs to be placed temporarily, care must be taken during removal as frames can be subject to damage.

CE 0297

SCHÜTZ DENTAL
Micerium Group

Schütz Dental GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany
Telefon +49 (0) 6003 814-0 · Telefax +49 (0) 6003 814-906
info@schuetz-dental.de · www.schuetz-dental.de