dialog Pastenopaker

!

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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Instructions for use

Product designation
Light-curing paste opaquer for the crown and bridge technique with veneering composites.

General product description
Light-curing one-component material for color coverage of framework retentions made of conventional dental metals and alloys. The material is light-curing and it is polymerized in a light-curing system.

Composition
Silicon dioxide, glass powder, urethane dimethacrylate, pigments, initiators.
Filler content: 62% by weight; anorganic fillers (0.005 - 3.0 µm)

Indications
For covering framework materials in the crown and bridge technique (conventional dental metals/dental alloys), for new restorations or repairs.

Contraindication
Do not use in case of a known allergy against one of the components.

Interaction with other agents
Phenolic substances (e.g. Eugenol) inhibit polymerization. Therefore, do not use any materials containing such substances.

Application / processing
Working range of the materials: 1- 3 min., depending on the lighting conditions.

Preparatory steps
The frameworks are designed, cast or milled in CAD / CAM, and processed in accordance with the applicable guidelines of dental technology. Conventional dental metals and alloys, PMMA, PEEK and zirconium dioxide can be used as framework materials.

Troubleshooting / FAQ list
– The framework surfaces must be free of grease and polishing residue.
– Bond opaque must not be mixed with opaquer liquid or powder.
– Opaquer layers which are too thick prevent an optimal polymerization and thus weaken the bond.
– Extending the polymerization time does not result in a greater polymerization depth.
– Do not use intermediate polymerization systems for polymerization.
– Opaquers must only be applied in very thin layers.

Problem	Cause	Remedy
does not cure	applied layer too thick <ul style="list-style-type: none">insufficient polymerization	only apply wash-like in several layers <ul style="list-style-type: none">observer polymerization times <ul style="list-style-type: none">– check lamp, replace lamp if necessary – extending the polymerization time does not result in a greater polymerization depth – do not use intermediate polymerization systems
opaquer chips	opaquer not polymerized correctly / applied layer too thick	only apply wash-like respectively in very thin layers <ul style="list-style-type: none">incorrect polymerization times / check lamp / replace lamp if necessary
material is too hard inside the syringe	syringe stored at temperatures under 10 °C (50 °F)	store material at room temperature or let material reach room temperature in time for application
greasy surface	insufficient polymerization	observe polymerization times <ul style="list-style-type: none">check system / regular system maintenance

Guarantee
Our processing instructions, regardless of whether they are issued verbally, in writing or in the form of a practical demonstration, are based upon our own experience and can therefore be regarded solely as guidelines. Our products are subject to continuous development. We therefore reserve the right to make possible changes.

General notes
Use an appropriate light polymerization system with an emission spectrum of at least 310 - 500 nm. The necessary physical properties will only be achieved with a flawless lamp. Therefore, it is necessary to check the light intensity regularly in accordance with the manufacturer’s instructions. Do not use any intermediate polymerization systems.

Storage
Storage temperature 10 -25 °C / 50 -77 °F.
Close syringe well after each use.

Shelf-life
The maximum shelf-life is printed on the label of each packaging. Do not use after the expiration date.

Side-effects
With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (e.g. allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-indications / interactions
If a patient has known hyper-sensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient’s mouth before using the product.

Note
Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

Disposal
After pretreatment, dispose of in an approved hazardous waste recycling facility in compliance with the hazardous waste regulations. Do not dispose in household waste. Do not allow to enter the sewage system.

All materials are to be processed and prepared for further processing (sand-blasting, conditioning) in accordance with the manufacturer’s instructions. The frameworks must be conditioned before the opaquer is used. The method of framework conditioning depends on the framework material used. Please observe the manufacturer’s instructions.

Bonding agents and framework materials
Sebond Smart serves as a bonding agent between acrylics / composites and metals / alloys.
Sebond Implant serves as a bonding agent between acrylics / composites and zirconium dioxide.
Bonding (e.g. dialog Bonding Fluid, Schütz Dental) serves as a bonding agent between acrylics/composites and PEEK or PMMA materials.

Procedure
Apply dialog paste opaquer thinly (wash-like) with a tight shorthaired brush. In order to reach a complete coverage, apply at least two layers. An insufficient opaque layer will negatively influence the color effect. If necessary, optimize the flow behavior of dialog paste opaquer by stirring it on a mixing block. Make sure to apply an extra thin coat along the retentions. Polymerize each layer separately (see table „polymerization times“).

Subsequent worksteps
Finally, the restoration is designed with a conventional veneering composite. Please observe the manufacturer’s instructions.

Repairs
All repairs are to be carried out extraorally. Prepare the surfaces to be repaired as described under „Preparatory steps“. Follow the procedures described under „Procedure“ and „Subsequent worksteps“.

Light curing system	Time
Spektra LED	1 min.
Spektra 2000	3 min.
HiLite / UniXS	90 sec.
Labolight LV-II / III	1 min.
Solidilite	1 min.

Problem	Cause	Remedy
does not cure	applied layer too thick <ul style="list-style-type: none">insufficient polymerization	only apply wash-like in several layers <ul style="list-style-type: none">observer polymerization times <ul style="list-style-type: none">– check lamp, replace lamp if necessary – extending the polymerization time does not result in a greater polymerization depth – do not use intermediate polymerization systems
opaquer chips	opaquer not polymerized correctly / applied layer too thick	only apply wash-like respectively in very thin layers <ul style="list-style-type: none">incorrect polymerization times / check lamp / replace lamp if necessary
material is too hard inside the syringe	syringe stored at temperatures under 10 °C (50 °F)	store material at room temperature or let material reach room temperature in time for application
greasy surface	insufficient polymerization	observe polymerization times <ul style="list-style-type: none">check system / regular system maintenance

Color according to Vita®*	Paste opaque		
	Cervical	Dentin	Incisal
A 1	P-H 1	P-D 1	P-S 1
A2		P-D 2	
A3	P-H 2	P-D 3	
A3,5		P-D 4	
A4		P-D 4	
B1	P-H 1	P-D 5	
B2		P-D 6	
B3		P-D 7	
B4		P-D 7	
C1	P-H 3	P-D 8	
C2		P-D 9	
C3		P-D 9	
C4		P-D 8	
D2		P-D 10	
D3	P-H 3	P-D 2	P-D 4
D4		P-D 4	

*Vita is a registered trademark of Vita® Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.

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*Vita es una marca registrada de la Vita® Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen.