

# MicroClean

Ask for  
an offer  
now!



Let the MicroClean system do the work –  
and you make the money!

**Hygienic cleaning of prostheses, crowns, bridges and many other types of restorations.**

- **Fast:** Dual cleaning system (mechanical and chemical cleaning in one step)
- **Gentle:** Specifically rounded cleaning needles protect your dental work.
- **Compact:** Fits in every dental practice and every lab.
- **Health promoting:** Clean old prostheses from calculus and plaque before getting started e. g. with the repair.



**SCHÜTZ DENTAL**  
Micerium Group

... suitable for  
orthodontics as well!

# MicroClean

Finally, no more manual cleaning of disgusting prostheses!

Even the toughest stains such as plaque and calculus are removed easily and automatically by MicroClean, including hard-to-reach areas.

Worn prostheses or new prosthetic constructions are cleaned and polished almost effortlessly – down to the finest interdental spaces, fissures and telescopic tips.

## Case study:



[www.sdent.eu/microclean-case-study](http://www.sdent.eu/microclean-case-study)



Product description	art. no.
<b>MicroClean, Denture Cleaning System</b>	<b>676410</b>
Needles short, 150 g	678866
Needles long, 150 g	678865
Cleaning fluid green (acidic), 500 ml	678857
Cleaning fluid orange (universal), 500 ml	678855
Plastic container, 500 ml	678820
Tweezers, 18.00 cm	678852
Holding magnet complete	678849
Measuring cylinder, 50 ml	676452
Stand-up display	678860

### Technical data:

Voltage:	230 V~
Power output:	200 VA
Max. revolutions per minute:	1,450 rpm
Acoustic level:	< 65 dB (A)
Dimensions:	185 x 175 x 195 mm
Weight:	7 kg



Further questions? We are here for you:

Telephone: +49 (0) 6003 814-365

E-mail: [export@schuetz-dental.de](mailto:export@schuetz-dental.de)

Pictures may vary from actual product appearance.

**SCHÜTZ DENTAL**  
Micerium Group