Printing date 05.02.2016

\*

\*

\*

Version: 5

Revision: 05.02.2016

SECTION 1: Identification of the substance/mixture and of the company/undertaking
· 1.1 Product identifier
· Trade name: <u>ReVeneer Base</u>
<ul> <li>Article number: 232210-1</li> <li>1.2 Relevant identified uses of the substance or mixture and uses advised against No further relevant information available.</li> <li>Application of the substance / the mixture Auxiliary for manufacture of dental prothesis</li> </ul>
<ul> <li>1.3 Details of the supplier of the safety data sheet</li> <li>Manufacturer/Supplier: Schütz Dental GmbH, Dieselstrasse 5-6, D-61191 Rosbach v.d.H. (Germany) Tel.: +49 (0)6003/814-0 Fax: +49 (0)6003/814-900 www.schuetz-dental.de; e-mail: info@schuetz-dental.de</li> </ul>
<ul> <li>Further information obtainable from: F&amp;E, Tel. +49 (0)6003/814-100</li> <li>1.4 Emergency telephone number: +49 (0) 6003 8140 Schütz Dental (8:00 - 17:00 Uhr) or +49 (0) 6131 19240 Poison Information Center, University Mainz (24 h)</li> </ul>
SECTION 2: Hazards identification
<ul> <li>• 2.1 Classification of the substance or mixture</li> <li>• Classification according to Regulation (EC) No 1272/2008 The product is not classified according to the CLP regulation.</li> </ul>
<ul> <li>2.2 Label elements</li> <li>Labelling according to Regulation (EC) No 1272/2008 Void</li> <li>Hazard pictograms Void</li> <li>Signal word Void</li> <li>Hazard statements Void</li> <li>Additional information: Contains cristobalite which may damage in a solid form the lung by inhalation. Because of liquid/pasty form of the product this danger is not given. Safety data sheet available on request.</li> <li>Labelling of packages where the contents do not exceed 125 ml</li> <li>Hazard pictograms Void</li> <li>Signal word Void</li> <li>Hazard statements Void</li> <li>2.3 Other hazards</li> <li>Results of PBT and vPvB assessment</li> <li>PBT: Not applicable.</li> </ul>
SECTION 3: Composition/information on ingredients
• <b>3.2 Chemical characterisation: Mixtures</b> • <b>Description:</b> Liquid based on methacryl acid ester, containing an activator.
• <b>Dangerous components:</b> Void • <b>Additional information</b> For the wording of the listed risk phrases refer to section 16.
SECTION 4: First aid measures
<ul> <li>4.1 Description of first aid measures</li> <li>General information Immediately remove any clothing soiled by the product.</li> <li>After inhalation Supply fresh air; consult doctor in case of complaints.</li> <li>After skin contact Immediately wash with water and soap and rinse thoroughly. If skin irritation continues, consult a doctor.</li> </ul>
(Contd. on page 2) GB

Printing date 05.02.2016

Version: 5

Revision: 05.02.2016

Trade name: ReVeneer Base

- (Contd. of page 1) • After eye contact Rinse opened eye for several minutes under running water. Then consult a doctor. • After swallowing
- *Rinse out mouth and then drink plenty of water. Seek medical treatment.*
- 4.2 Most important symptoms and effects, both acute and delayed No further relevant information available.
- $\cdot$  4.3 Indication of any immediate medical attention and special treatment needed
- No further relevant information available.

### **SECTION 5: Firefighting measures**

- · 5.1 Extinguishing media
- · Suitable extinguishing agents
- CO2, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- · For safety reasons unsuitable extinguishing agents Water with full jet.
- 5.2 Special hazards arising from the substance or mixture No further relevant information available.
- · 5.3 Advice for firefighters
- · Protective equipment:
- No special measures required.

Do not inhale explosion gases or combustion gases.

#### **SECTION 6: Accidental release measures**

- 6.1 Personal precautions, protective equipment and emergency procedures Ensure adequate ventilation Keep people at a distance and stay on the windward side.
- 6.2 Environmental precautions: Do not allow to enter sewers/ surface or ground water.
- 6.3 Methods and material for containment and cleaning up:
- Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- · 6.4 Reference to other sections
- No dangerous substances are released.
- See Section 7 for information on safe handling
- See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

### **SECTION 7: Handling and storage**

- 7.1 Precautions for safe handling Ensure good ventilation/exhaustion at the workplace.
- · Information about fire and explosion protection: Keep ignition sources away Do not smoke.
- · 7.2 Conditions for safe storage, including any incompatibilities
- · Storage
- *Requirements to be met by storerooms and receptacles: Close container immediately after use and protect product from light incidence. Storage between 10 °C and 25 °C.*
- · Information about storage in one common storage facility: Not required.
- · Further information about storage conditions: None.
- 7.3 Specific end use(s) No further relevant information available.

### SECTION 8: Exposure controls/personal protection

- Additional information about design of technical facilities: No further data; see item 7.
- · 8.1 Control parameters
- · Ingredients with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

(Contd. on page 3)

в-

Printing date 05.02.2016

Version: 5

Revision: 05.02.2016

(Contd. of page 2)

Trade name: ReVeneer Base

• Additional information: The lists valid during the making were used as basis.

· 8.2 Exposure controls

· Personal protective equipment

· General protective and hygienic measures Immediately remove all soiled and contaminated clothing

• Respiratory protection: Not required.

· Protection of hands: Protective gloves.

· Material of gloves Butyl rubber, BR

• Penetration time of glove material

0,1 - 0,2 mm

Penetration time: 5 min.

· Eye protection: Safety glasses

9.1 Information on basic physical a	and chemical properties	
General Information		
Appearance:		
Form:	Pasty Light group	
Colour: Odour:	Light green Weak, characteristic	
Odour threshold:	Not determined.	
pH-value:	Not determined.	
Change in condition		
Melting point/Melting range:	undetermined	
Boiling point/Boiling range:	251 °C	
Flash point:	151 °C	
Flammability (solid, gaseous)	Not applicable.	
Ignition temperature:	> 370 °C	
Decomposition temperature:	Not determined.	
Self-igniting:	Product is not selfigniting.	
Danger of explosion:	Product does not present an explosion hazard.	
Explosion limits:		
Lower:	Not determined.	
Upper:	Not determined.	
Vapour pressure:	Not determined.	
Density:	Not determined	
Relative density	Not determined.	
Vapour density	Not determined.	
Evaporation rate	Not determined.	
Solubility in / Miscibility with		
Water:	Not miscible or difficult to mix	
Partition coefficient (n-octanol/wat	er): Not determined.	
Viscosity:		
dynamic:	Not determined.	
kinematic:	Not determined.	
Solvent content:		
Organic solvents:	0,0 %	

Printing date 05.02.2016

Version: 5

Revision: 05.02.2016

(Contd. of page 3)

Trade name: ReVeneer Base

• 9.2 Other information

No further relevant information available.

#### **SECTION 10: Stability and reactivity**

- · 10.1 Reactivity No further relevant information available.
- · 10.2 Chemical stability
- Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- · 10.3 Possibility of hazardous reactions Exothermic polymerization
- · 10.4 Conditions to avoid No further relevant information available.
- · 10.5 Incompatible materials:

In presence of radical formers (e. g. peroxides), deoxidizing substances, and/or heavy metal ions, polymerization with heat release is possible.

· 10.6 Hazardous decomposition products: No dangerous decomposition products known

#### **SECTION 11: Toxicological information**

· 11.1 Information on toxicological effects

- Acute toxicity Based on available data, the classification criteria are not met.
- · Primary irritant effect:
- Skin corrosion/irritation Based on available data, the classification criteria are not met.
- · Serious eye damage/irritation Based on available data, the classification criteria are not met.
- · Respiratory or skin sensitisation Based on available data, the classification criteria are not met.
- · CMR effects (carcinogenity, mutagenicity and toxicity for reproduction)
- · Germ cell mutagenicity Based on available data, the classification criteria are not met.
- · Carcinogenicity Based on available data, the classification criteria are not met.
- *Reproductive toxicity Based on available data, the classification criteria are not met.*
- · STOT-single exposure Based on available data, the classification criteria are not met.
- · STOT-repeated exposure Based on available data, the classification criteria are not met.
- · Aspiration hazard Based on available data, the classification criteria are not met.

### **SECTION 12: Ecological information**

- · 12.1 Toxicity
- · Aquatic toxicity: No further relevant information available.
- · 12.2 Persistence and degradability No further relevant information available.
- · 12.3 Bioaccumulative potential No further relevant information available.
- · 12.4 Mobility in soil No further relevant information available.
- · Additional ecological information:
- · General notes:

Water hazard class 2 (German Regulation) (Self-assessment): hazardous for water. Do not allow product to reach ground water, water course or sewage system. Danger to drinking water if even small quantities leak into the ground.

- · 12.5 Results of PBT and vPvB assessment
- **PBT:** Not applicable.
- · vPvB: Not applicable.
- · 12.6 Other adverse effects No further relevant information available.

#### **SECTION 13: Disposal considerations**

· 13.1 Waste treatment methods

· Recommendation

After prior treatment product has to be disposed of in an incinerator for hazardous waste adhering to the regulations pertaining to the disposal of particularly hazardous waste.

(Contd. on page 5)

Printing	date	05.	02.	2016
----------	------	-----	-----	------

Version: 5

Revision: 05.02.2016

(Contd. of page 4)

Trade name: ReVeneer Base

· European waste catalogue

08 01 12 waste paint and varnish other than those mentioned in 08 01 11

· Uncleaned packaging:

· Recommendation:

Packagings that may not be cleansed are to be disposed of in the same manner as the product.

· Recommended cleansing agents: Alcohol

SECTION 14: Transport informa	tion	
· 14.1 UN-Number · ADR, ADN, IMDG, IATA	Void	
· 14.2 UN proper shipping name · ADR, ADN, IMDG, IATA	Void	
· 14.3 Transport hazard class(es)		
· ADR, ADN, IMDG, IATA · Class	Void	
· 14.4 Packing group · ADR, IMDG, IATA	Void	
· 14.5 Environmental hazards: · Marine pollutant:	No	
· 14.6 Special precautions for user	Not applicable.	
· 14.7 Transport in bulk according to Ann Marpol and the IBC Code	<b>ex II of</b> Not applicable.	
· UN ''Model Regulation'':	Void	

## **SECTION 15: Regulatory information**

· 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

· Directive 2012/18/EU

· Named dangerous substances - ANNEX I None of the ingredients is listed.

· National regulations

- · Waterhazard class: Water hazard class 2 (Self-assessment): hazardous for water.
- 15.2 Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

### **SECTION 16: Other information**

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

• Department issuing MSDS: F & E

· Contact: Dr. U. Krichbaum

• Abbreviations and acronyms:

RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail)

ICAO: International Civil Aviation Organisation

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

(Contd. on page 6)

GB

Printing date 05.02.2016

Version: 5

Revision: 05.02.2016

(Contd. of page 5)

#### Trade name: ReVeneer Base

ELINCS: European List of Notified Chemical Substances CAS: Chemical Abstracts Service (division of the American Chemical Society) PBT: Persistent, Bioaccumulative and Toxic vPvB: very Persistent and very Bioaccumulative • \* Data compared to the previous version altered.

GB