

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**Tizian Blank PMMA Splint**

Revision date: 29.08.2019

Product code: 214

Page 1 of 11

**SECTION 1: Identification of the substance/mixture and of the company/undertaking****1.1. Product identifier**

Tizian Blank PMMA Splint

**1.2. Relevant identified uses of the substance or mixture and uses advised against****Use of the substance/mixture**

The product is intended for professional use.

**Uses advised against**

Do not use for private purposes (household).

**1.3. Details of the supplier of the safety data sheet**

Company name:	Schütz Dental GmbH	
Street:	Dieselstr. 5-6	
Place:	D-61191 Rosbach / Germany	
Telephone:	+49 (0) 6003-814-0	Telefax: +49 (0) 6003-814-906
e-mail:	info@schuetz-dental.de	
Internet:	https://www.schuetz-dental.de/de/	

**1.4. Emergency telephone number:**

+49-(0)551-19240 (Giftinformationszentrum-Nord)

**Further Information**

none

**SECTION 2: Hazards identification****2.1. Classification of the substance or mixture****Regulation (EC) No. 1272/2008**

This mixture is not classified as hazardous in accordance with Regulation (EC) No. 1272/2008.

**2.2. Label elements****Regulation (EC) No. 1272/2008****Special labelling of certain mixtures**

EUH208	Contains methyl 2-methylprop-2-enoate; methyl 2-methylpropenoate; methyl methacrylate, dibenzoyl peroxide; benzoyl peroxide. May produce an allergic reaction.
EUH210	Safety data sheet available on request.

**Additional advice on labelling**

There is no requirement for the product to be specially labelled according to EC directives or the corresponding national laws.

**2.3. Other hazards**

Harmful dust is produced during dry-state pulverisation. Inhalation of dust may cause irritation of the respiratory system.

**SECTION 3: Composition/information on ingredients****3.2. Mixtures****Chemical characterization**

POLYMETHYL METHACRYLATE

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

## Tizian Blank PMMA Splint

Revision date: 29.08.2019

Product code: 214

Page 2 of 11

### Hazardous components

CAS No	Chemical name			Quantity
	EC No	Index No	REACH No	
	Classification according to Regulation (EC) No. 1272/2008 [CLP]			
	PMMA			> 98 %
80-62-6	methyl 2-methylprop-2-enoate; methyl 2-methylpropenoate; methyl methacrylate			< 1 %
	201-297-1	607-035-00-6		
	Flam. Liq. 2, STOT SE 3, Skin Irrit. 2, Skin Sens. 1; H225 H335 H315 H317			
94-36-0	dibenzoyl peroxide; benzoyl peroxide			< 1 %
	202-327-6	617-008-00-0		
	Org. Perox. B, Eye Irrit. 2, Skin Sens. 1; H241 H319 H317			

Full text of H and EUH statements: see section 16.

### Further Information

May cause sensitisation by skin contact.

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

#### General information

No special measures are necessary.

#### After inhalation

Inhalation of dust may cause irritation of the respiratory system. In case of allergic symptoms, especially in the breathing area, seek medical advice immediately. Provide fresh air. In all cases of doubt, or when symptoms persist, seek medical advice.

#### After contact with skin

After contact with skin, wash immediately with plenty of water and soap. In case of skin irritation, consult a physician.

#### After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water. In case of eye irritation consult an ophthalmologist.

#### After ingestion

No known symptoms to date. Observe risk of aspiration if vomiting occurs. When in doubt or if symptoms are observed, get medical advice.

### 4.2. Most important symptoms and effects, both acute and delayed

Allergic reactions

### 4.3. Indication of any immediate medical attention and special treatment needed

When in doubt or if symptoms are observed, get medical advice.

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

#### Suitable extinguishing media

Water mist, Foam, Carbon dioxide (CO<sub>2</sub>), Extinguishing powder

#### Unsuitable extinguishing media

High power water jet

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

**Tizian Blank PMMA Splint**

Revision date: 29.08.2019

Product code: 214

Page 3 of 11

**5.2. Special hazards arising from the substance or mixture**

This material is combustible, but will not ignite readily. Exothermal decomposition with formation of: Carbon dioxide (CO<sub>2</sub>), Carbon monoxide.

**5.3. Advice for firefighters**

Wear a self-contained breathing apparatus and chemical protective clothing. Co-ordinate fire-fighting measures to the fire surroundings. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water. Use caution when applying carbon dioxide in confined spaces. Carbon dioxide can displace oxygen.

**Additional information**

Do not inhale explosion and combustion gases. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

**SECTION 6: Accidental release measures****6.1. Personal precautions, protective equipment and emergency procedures**

Avoid dust formation. Wear breathing apparatus if exposed to vapours/dusts/aerosols. Use personal protection equipment.

**6.2. Environmental precautions**

No special environmental measures are necessary. Avoid release to the environment.

**6.3. Methods and material for containment and cleaning up**

Wet clean or vacuum up solids. Use approved industrial vacuum cleaner for removal. Do not use a dry brush as dust clouds or static can be created. Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

**6.4. Reference to other sections**

Safe handling: see section 7. Personal protection equipment: see section 8. Disposal: see section 13.

**SECTION 7: Handling and storage****7.1. Precautions for safe handling****Advice on safe handling**

Measures to prevent aerosol and dust generation. Do not breathe dust. Dust should be exhausted directly at the point of origin. If local exhaust ventilation is not possible or not sufficient, the entire working area should be ventilated by technical means.

**Advice on protection against fire and explosion**

No special fire protection measures are necessary.

**Further information on handling**

Dust should be exhausted directly at the point of origin.

**7.2. Conditions for safe storage, including any incompatibilities****Requirements for storage rooms and vessels**

No special measures are necessary.

**Advice on storage compatibility**

To follow: Storage class

**Further information on storage conditions**

No special measures are necessary.

**7.3. Specific end use(s)**

Specific use(s): Observe instructions for use.

**SECTION 8: Exposure controls/personal protection****8.1. Control parameters**

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

## Tizian Blank PMMA Splint

Revision date: 29.08.2019

Product code: 214

Page 4 of 11

## Exposure limits (EH40)

CAS No	Substance	ppm	mg/m <sup>3</sup>	fibres/ml	Category	Origin
94-36-0	Dibenzoyl peroxide	-	5		TWA (8 h)	WEL
		-	-		STEL (15 min)	WEL
80-62-6	Methyl methacrylate	50	208		TWA (8 h)	WEL
		100	416		STEL (15 min)	WEL

## DNEL/DMEL values

CAS No	Substance	Exposure route	Effect	Value
94-36-0	dibenzoyl peroxide; benzoyl peroxide			
Consumer DNEL, long-term	inhalation	systemic	2,9 mg/m <sup>3</sup>	
Worker DNEL, long-term	dermal	systemic	6,6 mg/kg bw/day	
Consumer DNEL, long-term	oral	systemic	1,65 mg/kg bw/day	
Worker DNEL, long-term	inhalation	systemic	11,75 mg/m <sup>3</sup>	
Consumer DNEL, long-term	dermal	systemic	3,3 mg/kg bw/day	

## PNEC values

CAS No	Substance	Value
94-36-0	dibenzoyl peroxide; benzoyl peroxide	
Freshwater		0,000602 mg/l
Marine water		0,000602 mg/l
Freshwater sediment		0,338 mg/kg
Micro-organisms in sewage treatment plants (STP)		0,35 mg/l
Soil		0,0758 mg/kg

## Additional advice on limit values

Monitoring and observation processes:

## 8.2. Exposure controls



## Appropriate engineering controls

If handled uncovered, arrangements with local exhaust ventilation should be used if possible. Hand tool with integrated exhaust. Use an "effective exhaust ventilation system" according to 2001/59/EC (Annex 7A)..

Preferably use hand tools or low-speed tools equipped, if necessary, with an appropriate dust-extraction facility.

If high-speed tools are used, they should always be equipped with such a facility.

## Protective and hygiene measures

When using do not eat, drink, smoke, sniff. Do not breathe dust. General health and safety measures. Wash hands before breaks and after work. When using do not eat or drink.

## Eye/face protection

Dust protection eye glasses. Wear eye/face protection.

## Hand protection

Tested protective gloves must be worn. Wear cotton undermitten if possible. Take recovery periods for skin

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

## Tizian Blank PMMA Splint

Revision date: 29.08.2019

Product code: 214

Page 5 of 11

regeneration.

### Skin protection

Wear suitable protective clothing.

### Respiratory protection

If technical exhaust or ventilation measures are not possible or insufficient, respiratory protection must be worn. Respiratory protection necessary at: exceeding exposure limit values. Observe the wear time limits according GefStoffV in combination with the rules for using respiratory protection apparatus (BGR 190).

### Environmental exposure controls

No special environmental measures are necessary.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state:	solid
Colour:	transparent, colourless
Odour:	odourless

#### Test method

pH-Value:	not applicable
-----------	----------------

#### Changes in the physical state

Melting point:	not applicable
----------------	----------------

Initial boiling point and boiling range:	not applicable
--	----------------

Sublimation point:	not applicable
--------------------	----------------

Softening point:	121 °C	ISO 11357
------------------	--------	-----------

Pour point:	not applicable
-------------	----------------

not applicable:	
-----------------	--

Flash point:	> 250 °C	ASTM D 1929-68
--------------	----------	----------------

Sustaining combustion:	No data available
------------------------	-------------------

#### Flammability

Solid:	not determined
--------	----------------

Gas:	not applicable
------	----------------

#### Explosive properties

The product in the delivered form is not dust explosion capable; the enrichment of fine dust however leads to the danger of dust explosion.

Lower explosion limits:	not determined
-------------------------	----------------

Upper explosion limits:	not determined
-------------------------	----------------

Ignition temperature:	> 400 °C	ASTM D 1929-68
-----------------------	----------	----------------

#### Auto-ignition temperature

Solid:	not determined
--------	----------------

Gas:	not applicable
------	----------------

Decomposition temperature:	> 250 °C
----------------------------	----------

#### Oxidizing properties

No data available

Vapour pressure:	not applicable
------------------	----------------

Vapour pressure:	not applicable
------------------	----------------

Density:	1,189 – 1,198 g/cm³
----------	---------------------

Bulk density:	not applicable
---------------	----------------

Water solubility:	insoluble
-------------------	-----------

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**Tizian Blank PMMA Splint**

Revision date: 29.08.2019

Product code: 214

Page 6 of 11

**Solubility in other solvents**

practically insoluble

Partition coefficient:

not applicable

Viscosity / dynamic:

not applicable

Viscosity / kinematic:

not applicable

Flow time:

not applicable

Vapour density:

not applicable

Evaporation rate:

not applicable

Solvent separation test:

not applicable

Solvent content:

not applicable

**9.2. Other information**

Solid content:

100 %

none

**SECTION 10: Stability and reactivity****10.1. Reactivity**

No hazardous reaction when handled and stored according to provisions.

**10.2. Chemical stability**

The product is chemically stable under recommended conditions of storage, use and temperature.

Decomposition takes place from temperatures above: &gt; 250 °C

**10.3. Possibility of hazardous reactions**

No hazardous reaction when handled and stored according to provisions.

**10.4. Conditions to avoid**

Harmful dust is produced during dry-state pulverisation.

**10.5. Incompatible materials**

Oxidising agent, strong

**10.6. Hazardous decomposition products**

Thermal decomposition can lead to the escape of irritating gases and vapours. Hazardous decomposition products: SECTION 8: Exposure controls/personal protection (Control parameters)

**Further information**

none

**SECTION 11: Toxicological information****11.1. Information on toxicological effects****Toxicokinetics, metabolism and distribution**

No information available.

**Acute toxicity**

Based on available data, the classification criteria are not met.

No information available.

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

## Tizian Blank PMMA Splint

Revision date: 29.08.2019

Product code: 214

Page 7 of 11

CAS No	Chemical name				
	Exposure route	Dose	Species	Source	Method
80-62-6	methyl 2-methylprop-2-enoate; methyl 2-methylpropenoate; methyl methacrylate				
	oral	LD50 > 5000 mg/kg	Rat	OECD 401	
	dermal	LD50 > 5000 mg/kg	Rabbit		
	inhalative (4 h) vapour	LC50 29,8 mg/l	Rat		
94-36-0	dibenzoyl peroxide; benzoyl peroxide				
	oral	LD50 > 2000 mg/kg	Mouse	Nier, Korea 2001	OECD 401
	inhalative aerosol	LC50 > 24300 mg/l	Rat		

### Irritation and corrosivity

Based on available data, the classification criteria are not met.  
No information available.

### Sensitising effects

Based on available data, the classification criteria are not met.  
May cause sensitisation especially in sensitive humans. People who suffer from skin sensitization problems, asthma, allergies, chronic or recurring respiratory illnesses should not be deployed in any process using this preparation.

### Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.  
Based on existing data the substance does not fulfill the criteria of CMR-substances Cat. 1 and 2 according 67/548/EEC. No indications of human germ cell mutagenicity exist. No indications of human reproductive toxicity exist.

### STOT-single exposure

Based on available data, the classification criteria are not met.  
No information available.

### STOT-repeated exposure

Based on available data, the classification criteria are not met.  
No information available.

### Aspiration hazard

Based on available data, the classification criteria are not met.  
No information available.

### Specific effects in experiment on an animal

No information available.

### Additional information on tests

none

### Practical experience

### Observations relevant to classification

No information available.

### Other observations

No information available.

### Further information

Toxicological data are not available. May cause allergy or asthma symptoms or breathing difficulties if inhaled.  
Do not breathe dusts or mists.

## SECTION 12: Ecological information

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

## Tizian Blank PMMA Splint

Revision date: 29.08.2019

Product code: 214

Page 8 of 11

**12.1. Toxicity**

Due to the consistency along with the low water solubility of the product a bioavailability is unlikely. According to the criteria of the European classification and labelling system, the substance/the product has not to be labelled as "dangerous for the environment".

CAS No	Chemical name					
	Aquatic toxicity	Dose	[h]   [d]	Species	Source	Method
80-62-6	methyl 2-methylprop-2-enoate; methyl 2-methylpropenoate; methyl methacrylate					
	Acute fish toxicity	LC50 > 79 mg/l	96 h	Oncorhynchus mykiss (Rainbow trout)	OECD 203	
	Acute algae toxicity	ErC50 > 110 mg/l	72 h	Selenastrum capricornutum	OECD 201	
	Acute crustacea toxicity	EC50 69 mg/l	48 h	Daphnia magna (Big water flea)	OECD 202	
94-36-0	dibenzoyl peroxide; benzoyl peroxide					
	Acute fish toxicity	LC50 0,24 mg/l	96 h	Oryzias latipes (Ricefish)	Nier, Korea 2002c	OECD 203
	Acute algae toxicity	ErC50 0,44 mg/l	72 h	Selenastrum capricornutum	Nier, Korea 2002f	OECD 201
	Acute crustacea toxicity	EC50 0,07 mg/l	48 h	Daphnia pulex (water flea)	Nier, Korea 2002g	OECD 202
	Algae toxicity	NOEC 0,02 mg/l	3 d	Pseudokirchneriella subcapitata	REACH Dossier	EU Method C.3
	Crustacea toxicity	NOEC 0,001 mg/l	21 d	Daphnia pulex (water flea)	REACH Dossier	OECD 211
	Acute bacteria toxicity	0,30 g O <sub>2</sub> /g (0,35 mg/l)	0,5 h	activated sludge	REACH Dossier	OECD 209

**12.2. Persistence and degradability**

There are no data available on the preparation/mixture itself. According to experiences this product is inert and not degradable.

CAS No	Chemical name			
	Method	Value	d	Source
	Evaluation			
80-62-6	methyl 2-methylprop-2-enoate; methyl 2-methylpropenoate; methyl methacrylate			
	OECD 301C/ ISO 9408/ EEC 92/69/V, C.4-F	94 %	14	SDB Evonik
	Readily biodegradable (according to OECD criteria).			
94-36-0	dibenzoyl peroxide; benzoyl peroxide			
	OECD 301D/ EEC 92/69/V, C.4-E	71 %	28	REACH Dossier
	Readily biodegradable (according to OECD criteria).			

**12.3. Bioaccumulative potential**

Water solubility < 1 mg/L at 20°C. The product has not been tested. On the basis of existing data about the elimination/degradation and bioaccumulation potential longer term damage to the environment is unlikely.

**Partition coefficient n-octanol/water**

CAS No	Chemical name	Log Pow
80-62-6	methyl 2-methylprop-2-enoate; methyl 2-methylpropenoate; methyl methacrylate	1,38
94-36-0	dibenzoyl peroxide; benzoyl peroxide	3,2

**BCF**

CAS No	Chemical name	BCF	Species	Source
94-36-0	dibenzoyl peroxide; benzoyl peroxide	47,4	n/n	EpiSuite QSAR tool

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

## Tizian Blank PMMA Splint

Revision date: 29.08.2019

Product code: 214

Page 9 of 11

### 12.4. Mobility in soil

No data available.

### 12.5. Results of PBT and vPvB assessment

No information available.

### 12.6. Other adverse effects

No data available.

### Further information

none

## SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

#### Advice on disposal

Waste disposal according to EC directives 75/442/EEC and 91/689/EEC in the corresponding versions, covering waste and dangerous waste. Dispose according to legislation. Can be incinerated together with household waste in compliance with applicable technical regulations following consultation with approved waste disposal management companies and authorities in charge.

#### Waste disposal number of waste from residues/unused products

070213 WASTES FROM ORGANIC CHEMICAL PROCESSES; wastes from the MFSU of plastics, synthetic rubber and man-made fibres; waste plastic

#### Waste disposal number of used product

070213 WASTES FROM ORGANIC CHEMICAL PROCESSES; wastes from the MFSU of plastics, synthetic rubber and man-made fibres; waste plastic

#### Waste disposal number of contaminated packaging

150102 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); plastic packaging

#### Contaminated packaging

Non-contaminated packages may be recycled. Handle contaminated packages in the same way as the substance itself. Recycling possible without special treatment.. May be disposed of in household waste landfill. Recycle sales packaging via DSD (Duales System Deutschland).

## SECTION 14: Transport information

### Land transport (ADR/RID)

<u>14.1. UN number:</u>	Not restricted
<u>14.2. UN proper shipping name:</u>	Not restricted
<u>14.3. Transport hazard class(es):</u>	Not restricted
<u>14.4. Packing group:</u>	Not restricted

### Inland waterways transport (ADN)

<u>14.1. UN number:</u>	Not restricted
<u>14.2. UN proper shipping name:</u>	Not restricted
<u>14.3. Transport hazard class(es):</u>	Not restricted
<u>14.4. Packing group:</u>	Not restricted

### Marine transport (IMDG)

<u>14.1. UN number:</u>	Not restricted
<u>14.2. UN proper shipping name:</u>	Not restricted
<u>14.3. Transport hazard class(es):</u>	Not restricted
<u>14.4. Packing group:</u>	Not restricted

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

## Tizian Blank PMMA Splint

Revision date: 29.08.2019

Product code: 214

Page 10 of 11

### Air transport (ICAO-TI/IATA-DGR)

<b>14.1. UN number:</b>	Not restricted
<b>14.2. UN proper shipping name:</b>	Not restricted
<b>14.3. Transport hazard class(es):</b>	Not restricted
<b>14.4. Packing group:</b>	Not restricted

### 14.5. Environmental hazards

ENVIRONMENTALLY HAZARDOUS: no

### 14.6. Special precautions for user

Not restricted

### 14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not restricted

### Other applicable information

No dangerous good in sense of this transport regulation.

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### EU regulatory information

2010/75/EU (VOC): not applicable

2004/42/EC (VOC): not applicable

#### Additional information

none

#### National regulatory information

Employment restrictions: Observe restrictions to employment for juvenils according to the 'juvenile work protection guideline' (94/33/EC). Observe employment restrictions under the Maternity Protection Directive (92/85/EEC) for expectant or nursing mothers. Observe employment restrictions for women of child-bearing age.

Water contaminating class (D): 1 - slightly water contaminating

Skin resorption/Sensitization: Causes allergic hypersensitivity reactions.

### 15.2. Chemical safety assessment

Chemical safety assessments for substances in this mixture were not carried out.

## SECTION 16: Other information

### Changes

This data sheet contains changes from the previous version in section(s): 1,7,9.

### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service

LC50: Lethal concentration, 50%

LD50: Lethal dose, 50%

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**Tizian Blank PMMA Splint**

Revision date: 29.08.2019

Product code: 214

Page 11 of 11

**Relevant H and EUH statements (number and full text)**

H225	Highly flammable liquid and vapour.
H241	Heating may cause a fire or explosion.
H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
EUH208	Contains methyl 2-methylprop-2-enoate; methyl 2-methylpropenoate; methyl methacrylate, dibenzoyl peroxide; benzoyl peroxide. May produce an allergic reaction.
EUH210	Safety data sheet available on request.

**Further Information**

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.

Data sources:

MSDS - Manufacturer

ECHA - REACH dossier

---

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*